

Product datasheet for **SC312749**

TPRG1L (NM_182752) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TPRG1L (NM_182752) Human Untagged Clone
Tag:	Tag Free
Symbol:	TPRG1L
Synonyms:	FAM79A; h-mover; mover
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC312749 representing NM_182752. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
 GATCCGGTACCGAGGAGATCTGCCGCC**GCGATCGCC**
 ATGCTGCAACTGCGGGACTCGGTGGACTCGGCCGGTACGAGCCCCACGGCGGTGCTGGCGGCCGGCGAG
 GAGGTGGGGGACGGCGCGGCCCGGGCGGGGGCGGCCGGGGGCGGGGACGCCGCTGCGCCAGACACTC
 TGGCCTCTCAGCATCCACGACCCACGCGCCGCGCCCGCTCAAGGAGTACTTCGTGTTCCGGCCCGGC
 AGCATCGAGCAGGCAGTGGAGGAGATCCGCGTGGTGGTGGCGCCCGTGGAGGACGGCGAGATCCAGGGA
 GTGTGGCTGCTTACCGAGGTGGATCACTGGAACAATGAGAAGGAGCGGCTGGTGTGTCACGGAGCAG
 TCCCTGCTTATCTGTAATACGACTTCATCAGTCTCCAGTGCCAGCAGGTGGTGGGATAGCGCTCAAC
 GCAGTAGACACCATTTCTACGGAGAATCCAGTTTCCCCCTAAATCGCTCAACAAGCGAGAAGGTTTT
 GGGATTGCAATTCAGTGGGACAAGCAAAGTCGTCTTCTTCATAAACAGATGGAATCCCTGGTCTACC
 AACGTGCCCTATGCCACTTTACAGAACACCCGATGGCTGGCGCAGATGAGAAGACAGCATCTCTGTGT
 CAGTTGGAAGCTTCAAGGCTCTGTTAATCCAAGCTGTCAAAAAGCCAAAAAGAAAGCCCTTTGCCA
 GGACAGGCGAATGGCGTGCTGATCCTGGAGCGCCCCCTGCTCATCGAGACCTACGTGGGACTCATGTCC
 TTCATTAACAACGAGGCGAAACTGGGCTACTCCATGACCAGGGGCAAAATAGGCTTT**AG**
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
 TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC

Restriction Sites:	SgfI-MluI
ACCN:	NM_182752
Insert Size:	819 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_182752.3</u>
RefSeq Size:	2440 bp
RefSeq ORF:	819 bp
Locus ID:	127262
UniProt ID:	<u>Q5T0D9</u>
Cytogenetics:	1p36.32
MW:	30.2 kDa
Gene Summary:	Presynaptic protein involved in the synaptic transmission tuning. Regulates synaptic release probability by decreasing the calcium sensitivity of release.[UniProtKB/Swiss-Prot Function]