

Product datasheet for **SC312795**

Endonuclease V (ENDO V) (NM_173627) Human Untagged Clone


Product data:

Product Type:	Expression Plasmids
Product Name:	Endonuclease V (ENDO V) (NM_173627) Human Untagged Clone
Tag:	Tag Free
Symbol:	ENDO V
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC312795 representing NM_173627. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGCCCTGGAGGCGGGGAGGGCCGCCGAGGAAACGCTGTCCTGTGGAAACGGGAGCAAGCTCGG
CTGAAGGCCACGTCGTAGACCGGGACACCGAGGCGTGGCAGCGAGACCCCGCTTCTCGGGTCTGCAG
AGGGTCGGGGCGTTGACGTGTCCTTCGTGAAAGGGGACAGTGTCCGCGCTTGTGCTTCCCTGGTGGTG
CTCAGTTCCTGAGCTCGAGGTGGTGTATGAGGAGAGCCGCATGGTCAGCCTCACAGCCCCCTACGTG
TCGGGCTTCTGGCCTCCGAGAGGTGCCCTTCTTGTGGAGTGGTGCAGCAGCTGCGGAGAAGGAG
CCGGCCCTATGCCCAAGTCTTCTTGTGGATGGAAACGGGGTACTCCACCACCGAGGCTTTGGGGTG
GCCTGCCACCTTGGCGTCTTACAGACCTGCCGTGTGTTGGGGTGGCCAAGAACTTCTGCAGGTGGAT
GGGCTGGAGAACAACGCCCTGCACAAGGAGAAGATCCGACTCCTGCAGACTCGAGGAGACTCATTCCCT
CTGCTGGGAGACTCTGGGACTGTCTGGGAATGGCCCTGAGGAGCCACGACCGCAGCACCAGGCCCTC
TACATCTCCGTGGGCCACAGGATGAGCCTGGAGGCCGCTGTGCGCCTGACTTGCTGCTGCAGGTTTC
CGGATCCCAGAGCCCGTGCAGGCTGACATCTGCTCCCGAGAGCACATCCGCAAGTCGCTGGGACTC
CCCGGCCACCCACCGAGGAGCCGAAGGCGCAGAGGCCAGTGGCATGCCCAAGGAGACTCCGGA
GAGTCTCAGCACTTTGTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```

Restriction Sites: Sgfl-MluI



Plasmid Map:


ACCN: NM_173627

Insert Size: 849 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_173627.4](#)

RefSeq Size: 2840 bp

RefSeq ORF: 849 bp

Locus ID: 284131

UniProt ID: [Q8N8Q3](#)

Cytogenetics: 17q25.3

Protein Families: Druggable Genome

MW: 30.8 kDa

Gene Summary: Endoribonuclease that specifically cleaves inosine-containing RNAs: cleaves RNA at the second phosphodiester bond 3' to inosine. Has strong preference for single-stranded RNAs (ssRNAs) toward double-stranded RNAs (dsRNAs). Cleaves mRNAs and tRNAs containing inosine. Also able to cleave structure-specific dsRNA substrates containing the specific sites 5'-IIUI-3' and 5'-UIUU-3'. Inosine is present in a number of RNAs following editing; the function of inosine-specific endoribonuclease is still unclear: it could either play a regulatory role in edited RNAs, or be involved in antiviral response by removing the hyperedited long viral dsRNA genome that has undergone A-to-I editing. Binds branched DNA structures.[UniProtKB/Swiss-Prot Function]
Transcript Variant: This variant (1) represents the longest transcript and encodes the longest isoform (1).