

Product datasheet for **SC312878**

Ataxin 3 (ATXN3) (NM_030660) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	Ataxin 3 (ATXN3) (NM_030660) Human Untagged Clone
Tag:	Tag Free
Symbol:	ATXN3
Synonyms:	AT3; ATX3; JOS; MJD; MJD1; SCA3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC312878 representing NM_030660. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTT TAGTGAACCGTCAGAATTTTGT AATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCC CGGATCGCC
ATGGAGTCCATCTTCCACGAGAAACAGCCTTCTGGAAATATGGATGACAGTGGTTTTTCTCTATT
CGTTATAAGCAATGCCTTGAAAGTTGGGGTTTAGAACTAATCCTGTTCAACAGTCCAGAGTATCAGAGG
CTCAGGATCGATCCTATAAATGAAAGATCATTTATATGCAATTATAAGGAACACTGGTTTACAGTTAGA
AAATTAGGAAAACAGTGGTTAACTTGAATTCTCTCTTGACGGGTCAGAAATTAATATCAGATACATAT
CTTGCACCTTTCTGGCTCAATTACAACAGGAAGTTATTCTATATTTGTCGTTAAGGGTGATCTGCCA
GATTGCGAAGCTGACCAACTCCTGCAGATGATTAGGGTCCAACAGATGCATCGACCAAACTTATTGGA
GAAGAATTAGCACAACAAAAGAGCAAAGAGTCCATAAAACAGACCTGGAACGAGTGTTAGAAGCAAT
GATGGCTCAGGAATGTTAGACGAAGATGAGGAGGATTTGCAGAGGGCTCTGGCACTAAGTCGCCAAGAA
ATTGACATGGAAGATGAGGAAGCAGATCTCCGAGGGCTATTAGCTAAGTATGCAAGGTAGTTCAGAA
AACATATCTCAAGATATGACACAGACATCAGGTACAAATCTTACTTCAGAAGAGCTTCGGAAGAGACGA
GAAGCCTACTTTGAAAAACAGCAGCAAAAGCAGCAACAGCAGCAGCAGCAGCAGCAGGGGGACCTA
TCAGGACAGAGTTCACATCCATGTGAAAGGCCAGCCACCAGTTCAGGAGCACTTGGGAGTGATCTAGGT
GATGCTATGAGTGAAGAAGACATGCTTCAGGCAGCTGTGACCATGTCTTTAGAACTGTCAGAAATGAT
TTGAAAAACAGAAGGAAAAAAATAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI



[View online »](#)

Plasmid Map:


ACCN: NM_030660

Insert Size: 921 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_030660.4](#)

RefSeq Size: 6758 bp

RefSeq ORF: 921 bp

Locus ID:	4287
UniProt ID:	P54252
Cytogenetics:	14q32.12
Domains:	UIM, Josephin
Protein Families:	Druggable Genome, Transcription Factors
MW:	35 kDa

Gene Summary: Machado-Joseph disease, also known as spinocerebellar ataxia-3, is an autosomal dominant neurologic disorder. The protein encoded by this gene contains (CAG)_n repeats in the coding region, and the expansion of these repeats from the normal 12-44 to 52-86 is one cause of Machado-Joseph disease. There is a negative correlation between the age of onset and CAG repeat numbers. Alternatively spliced transcript variants encoding different isoforms have been described for this gene. [provided by RefSeq, Jul 2016]

Transcript Variant: This variant (h, also known as variant 2) is one of several transcript variants described in figure 2 of Bettencourt et al. (PMID: 19714377). This variant encodes isoform h (also known as isoform 2). Sequence Note: This RefSeq record was created from transcript and genomic sequence data because no single transcript was available for the full length of the gene. The extent of this transcript is supported by transcript alignments.