

Product datasheet for **SC313121**

HAPLN3 (AY037161) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	HAPLN3 (AY037161) Human Untagged Clone
Tag:	Tag Free
Symbol:	HAPLN3
Synonyms:	EXLD1; HsT19883
Vector:	<u>pCMV6 series</u>
Fully Sequenced ORF:	>NCBI ORF sequence for AY037161, the custom clone sequence may differ by one or more nucleotides

```

ATGGGCGTGTGCTCCTGGTCCCGTTGCTCCTGCTGCCCGGCTCCTACGGACTGCCCTTC
TACAACGGCTTCTACTACTCCAACAGCGCCAACGACCAGAACCTAGGCAACGGTCATGGC
AAAGACCTCCTTAATGGAGTGAAGCTGGTGGGAGACACCCGAGGAGACCCTGTTACCC
TACCAAGGGGCCAGTGTGATCCTGCCCTGCCGCTACCGCTACGAGCGGGCCCTGGTCTCC
CCGCGCGGTGTGCGTGTCAAATGGTGAAGCTGTGGGAGAACGGGGCCCCAGAGAAGGAC
GTGCTGGTGGCCATCGGGCTGAGGCACCGCTCCTTTGGGACTACCAAGCCGCGTGCAC
CTGCGGCAGGACAAAGAGCATGACGTCTCGCTGGAGATCCAGGATCTGCGGCTGGAGGAC
TATGGGCGTTACCGCTGTGAGGTATTGACGGGCTGGAGGATGAAAGCGGTCTGGTGGAG
CTGGAGCTGCGGGGTGTGGTCTTTCCTTACCAAGTCCCCAACGGGCGCTACCAAGTCAAC
TTCCACGAGGGCCAGCAGGTCTGTGCAGAGCAGGCTGCGGTGGTGGCCTCCTTTGAGCAG
CTCTTCCGGGCTGGGAGGAGGGCTGGACTGGTGAACGCGGGCTGGCTGCAGGATGCC
ACGGTGCAGTACCCCATCATGTTGCCCGGCAGCCCTGCGGTGGCCCGGGCTGGCACCT
GGCGTGCGAAGCTACGGCCCCGCCACCGCCGCTGCACCGCTATGATGTATTCTGCTTC
GCTACTGCCCTCAAGGGCGGGTGTACTACCTGGAGCACCTGAGAAGCTGACGCTGACA
GAGGCAAGGGAGGCCTGCCAGGAAGATGATGCCACGATCGCCAAGGTGGGACAGCTCTTT
GCCGCTGGAAGTTCCATGGCCTGGACCGCTGCGACGCTGGCTGGCTGGCAGATGGTAGC
GTCCGCTACCCCTGTGGTTCACCCGCATCCTAACTGTGGGCCCCAGAGCCTGGGGTCCGA
AGCTTTGGCTTCCCCGACCCGCAGAGCCGCTTGTACGGTGTCTTACTGCTACCGCCAGCAC

```

Restriction Sites:	Please inquire
ACCN:	AY037161
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	AY037161.1 , AAK67639.1
RefSeq Size:	1844 bp
RefSeq ORF:	1083 bp
Locus ID:	145864
Cytogenetics:	15q26.1
Domains:	Xlink, ig
Protein Families:	Secreted Protein
Gene Summary:	This gene belongs to the hyaluronan and proteoglycan binding link protein gene family. The protein encoded by this gene may function in hyaluronic acid binding and cell adhesion. [provided by RefSeq, Jul 2008]