

Product datasheet for SC313203

RSRC2 (NM_198261) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	RSRC2 (NM_198261) Human Untagged Clone
Tag:	Tag Free
Symbol:	RSRC2
Synonyms:	arginine/serine-rich coiled-coil 2; FLJ11021
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC313203 representing NM_198261. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTGTGTAACCGTCAGAATTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
 GATCCGGTACCGAGGAGATCTGCCGCC**CGCATCGCC**
 ATGATAAGAACCAACTTCTTCTTAAACAGGGAAGAAGACATGAATCCAAAGATAAATCCTCTAAGAAA
 CATAAGTCTGAGGAACATAATGACAAAGAACATTCTTCTGATAAAGGAAGAGAGCGACTAAATTCATCT
 GAAAATGGTGAGGACAGGCACAAACGAAAGAAAGTATCAAGAGGCAGAAGTCACTCAAGATCT
 AGGTCTCGTGAAAGACGCCATCGTAGTAGAAGCAGGGAGCGGAAGAAGTCTCGATCCAGGAGTAGGGAG
 CGGAAGAAATCGAGATCCAGAAGCAGAGAGAGGAAGAAATCGAGATCCAGAAGCAGGGAAGAAAACGG
 CGGATCAGGTCTCGTTCCCGCTCAAGATCAAGACACAGGCATAGGACTAGAAGCAGGAGTAGGACAAGG
 AGTAGGAGTCGAGATAGAAAGAAGAGAATTGAAAAGCCGAGAAGATTTAGCAGAAGTTTAAGCCGGACT
 CCAAGTCCACCTCCCTTCAGAGGCAGAAACACAGCAATGGATGCACAGGAAGCTTTAGCTAGAAGGTTG
 GAAAGGGCAAAGAAATTACAAGAACAGCGAGAAAAGGAAATGGTTGAAAAACAAAAACAAGAATA
 GCTGCAGCTGCAGCTACTGGAGTTCTGTTCTCAATGTTGCTGCCCTGTTGGCATCAGGAACACAAGTA
 ACACCTCAGATAGCCATGGCAGCTCAGATGGCAGCCCTGCAAGCTAAAGCTTTGGCAGAGACAGGAATA
 GCTGTTCTAGCTACTATAACCCAGCCGCTGTTAATCCAATGAAATTTGCTGAACAAGAGAAAAAAGG
 AAAATGCTTTGGCAGGGCAAGAAAGAGGGGACAAATCCCAATCTGCTGAAATATGGGAAAAATTGAAT
 TTTGGAACAAGGACCAAAATGTCAAATTTAGGAAATTTGATGGGTATTAAGAGTGAAGATGAAGCTGGA
 TGTAGCTCAGTTGATGAAGAAAGTTACAAGACTCTGAAGCAGCAGGAAGAAGTATTTCGAAATTTAGAT
 GCTCAGTATGAAATGGCAAGATCACAACCCACACAAAGAGGAATGGGTTTGGGTTTCACATCTTCA
 ATGCGAGGAATGGATGCAGTT**TGA**
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
 TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC

Restriction Sites: SgfI-MluI



[View online »](#)

ACCN:	NM_198261
Insert Size:	1128 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_198261.2</u>
RefSeq Size:	2352 bp
RefSeq ORF:	1128 bp
Locus ID:	65117
Cytogenetics:	12q24.31
MW:	43.8 kDa