

Product datasheet for SC313212

GM CSF Receptor alpha (CSF2RA) (NM_172246) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	GM CSF Receptor alpha (CSF2RA) (NM_172246) Human Untagged Clone
Tag:	Tag Free
Symbol:	CSF2RA
Synonyms:	alphaGMR; CD116; CDw116; CSF2R; CSF2RAX; CSF2RAY; CSF2RX; CSF2RY; GM-CSF-R-alpha; GMCSFR; GMCSFR-alpha; GMR; GMR-alpha; SMDP4
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC313212 representing NM_172246. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCTTCTCCTGGTGACAAGCCTTCTGCTCTGTGAGTTACCACACCCAGCATTCTCCTGATCCCAGAG
AAATCGGATCTGCGAACAGTGGCACCAGCCTCTAGTCTCAATGTGAGTTTGACTCCAGGACGATGAAT
TTAAGCTGGGACTGCCAAGAAAACAACCTTCAGCAAGTGTTCCTTAACTGACAAGAAGAACAGAGTC
GTGGAACCCAGGCTCAGTAACAACGAATGTTCTGTCACATTTCTGAAATTTGTCTGCATGAAGGAGTC
ACATTTGAGGTTACGTGAATACTAGTCAAAGAGGATTTCAACAGAACTGCTTTATCCAAATTCAGGA
AGGGAGGGTACCGCTGCTCAGAATTTCTCTGTTTCATCTACAATGCGGATTTAATGAACTGTACCTGG
GCGAGGGGTCCGACGGCCCCCGTGACGTCCAGTATTTTTGTACATACGAACTCAAAGAGAAGGAGG
GAGATCCGGTGTCTTATTACATACAAGACTCAGGAACCCATGTGGGATGTCACCTGGATAACCTGTCA
GGATTAACGTCTCGCAATTACTTTCTGGTTAACGGAACCGAGAAATGGCATCCAATTCCTTTGAT
TCACTTTTGACACAAAGAAAATAGAACGATTC AACCTCCCAGCAATGTACCGTACGTTGCAACACG
ACGCACTGCCTCGTACGGTGAAACAGCCAGGACCTATCAGAAGCTGTCGTACCTGGACTTTCAGTAC
CAGCTGGACGTCCACAGAAAGAAATACCCAGCCTGGCAGGAAAACCTACTGATTAATGTTTCTGGTAT
TTGAAAATAGATACAATTTCCAAGCTCTGAGCCAGAGCAAAACACAGTGTGAAGATCAGAGCTGCA
GACGTCCGCATCTTGAATGGAGCTCCTGGAGTGAAGCCATTGAATTTGATCATCTGGGAGGAATTCAC
CCCAGAGGAAGGAAAGGCTACCGCGAAGAGGCTTGACCGTGAAGGAAATTACCTGAGACCCAGAGGG
TGTAGGAATGGCATGGACATCTCCGCTCCGCGACACGGGGAACTGTTTTCTTGATGATGCTGTGAAC
CTTTATATCATTCTATGTTTTTATTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI



[View online >](#)

Plasmid Map:


ACCN: NM_172246

Insert Size: 1134 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_172246.2](#)

RefSeq Size: 1676 bp

RefSeq ORF: 1134 bp

Locus ID: 1438

UniProt ID: [P15509](#)

Cytogenetics: X;Y

Protein Families:	Druggable Genome, Secreted Protein, Transmembrane
Protein Pathways:	Cytokine-cytokine receptor interaction, Hematopoietic cell lineage, Jak-STAT signaling pathway, Pathways in cancer
MW:	43.5 kDa
Gene Summary:	<p>The protein encoded by this gene is the alpha subunit of the heterodimeric receptor for colony stimulating factor 2, a cytokine which controls the production, differentiation, and function of granulocytes and macrophages. The encoded protein is a member of the cytokine family of receptors. This gene is found in the pseudoautosomal region (PAR) of the X and Y chromosomes. Multiple transcript variants encoding different isoforms have been found for this gene, with some of the isoforms being membrane-bound and others being soluble. [provided by RefSeq, Jul 2008]</p> <p>Transcript Variant: This variant (3) lacks two coding exons in the 3' region, which causes a frameshift, compared to variant 1. The resulting isoform (b), which is likely soluble, is shorter and has a distinct C-terminus compared to isoform a.</p>