

Product datasheet for SC313219

DAZAP1 (NM_170711) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	DAZAP1 (NM_170711) Human Untagged Clone
Tag:	Tag Free
Symbol:	DAZAP1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC313219 representing NM_170711. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAACAACCTCGGGCGCCGACGAGATCGGGAAGCTCTTCGTGGGCGGTCTTGACTGGAGCAGACCCAA
GAGACTCTGCGCAGCTACTTTTCCCAATATGGAGAAGTCGTAGATTGTGTATCATGAAAGATAAAACC
ACCAACCACTCTCGAGGCTTTGGGTTTGTCAAATTTAAGACCCAACTGTGTGGGGACGGTGCTGGCC
AGCAGACCCGACACGCTAGATGGCCGAACATCGACCCCAAGCCATGCACACCCCGGGGGATGCAGCCG
GAGAGAACACGGCCGAAGGAAGGATGGCAGAAAGGACCCAGGAGCGATAACAGTAAATCAAATAAGATA
TTTGTCCGTGGAATTCACAATTTGGTGTGACAGAGCTCAGGGAATACTTCAAGAAGTTCGGAGTG
GTCACGGAGGTAGTCATGATCTATGACGCCGAGAAGCAGAGGCCCGAGGTTTTGGATTTATTACTTTC
GAGGACGAACAATCAGTGGACCAGGCTGTCAACATGCATTTTTCAGCAGATCATGGCAAAAAAGTGAA
GTTAAACGAGCTGAGCCTCGGGACAGCAAGAGCCAAGCGCCGGGACAGCCAGGTGCCAGCCAGTGGGGG
AGCCGGGTTGTGCCAACGCTGCCAATGGCTGGGACAGCCAGCCCGCCACGTGGCAGCAAGGATAT
GGCCCGCAAGGAATGTGGGTGCCGCGAGGACAGGCGATTGGTGGCTATGGACCGCCCTGCAGGAAGA
GGAGCCCCCGCCACCCCAACCGTTCACCTCTACATCGTGTCCACCCCTCTGGAGGCTTTCCCTT
CCCCAGGGCTTCCCTCAGGGCTACGGTGCCTCCGACAGTTTCAAGTTTGGCTACGGGCTCCACCTCA
CCGCCAGATCAGTTTGCCTCCGGGGTTCCTCTCCACAGCCACTCCGGGGCAGCACCTCTGGCT
TTCCACCGCCTCCGTCTCAGGCTGCCCGGACATGAGCAAGCCCGACAGCTCAGCCAGACTTCCC
TATGGTCAAGTATGGCTGGGTTCTATTCTCCAGCCCGCGGGTTCGGGCCACACTTTGTTTACAGT
CTTATGGTCAGGCTGAGCAGTGTATGGCTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: Sgfl-MluI



[View online »](#)

Plasmid Map:


ACCN: NM_170711

Insert Size: 1137 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_170711.2](#)

RefSeq Size: 2304 bp

RefSeq ORF: 1137 bp

Locus ID: 26528

UniProt ID: [Q96EP5](#)

Cytogenetics: 19p13.3

Protein Families: Stem cell - Pluripotency

MW: 40.5 kDa

Gene Summary: In mammals, the Y chromosome directs the development of the testes and plays an important role in spermatogenesis. A high percentage of infertile men have deletions that map to regions of the Y chromosome. The DAZ (deleted in azoospermia) gene cluster maps to the AZFc region of the Y chromosome and is deleted in many azoospermic and severely oligospermic men. It is thought that the DAZ gene cluster arose from the transposition, amplification, and pruning of the ancestral autosomal gene DAZL also involved in germ cell development and gametogenesis. This gene encodes a RNA-binding protein with two RNP motifs that was originally identified by its interaction with the infertility factors DAZ and DAZL. Two isoforms are encoded by transcript variants of this gene. [provided by RefSeq, Jul 2008]
Transcript Variant: This variant (1) represents the longer transcript and encodes the shorter isoform (a).