

Product datasheet for SC313260

OSCP1 (NM_145047) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	OSCP1 (NM_145047) Human Untagged Clone
Tag:	Tag Free
Symbol:	OSCP1
Synonyms:	C1orf102; NOR1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC313260 representing NM_145047. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTCCGGTGGGACGCTACCGCTGCTCTTCTTGAACCTGGGCGGGGAGATGCTTTACATCCTCGACCAA
CGGCTGCGGGCCAGAACATCCCGGGAGACAAGGCCCGCAAAGTTCTGAATGACATCATCTCCACCATG
TTCAATAGAAAAGTTTATGGAGGAATTATTCAAGCCTCAAGAGCTCTACTCCAAGAAGGCCCTGAGGACT
GTCTATGAGCGCCTGGCTCATGCCTCCATTATGAACTGAACCAGGCCAGCATGGATAAGCTCTATGAC
CTGATGACCATGGCTTTCAAATATCAAGTATTGCTGTGTCCCGACCCAAGGATGTGCTGCTGGTCACT
TTCAATCACTTGGATACCATCAAGGGATTATCCGAGACTCCCCAACCATCCTGCAGCAAGTGGACGAG
ACTTTGCGGCAGCTGACAGAAATATATGGTGGTCTCTCTGCAGGGGAGTTCCAGCTGATCCGGCAGACA
CTCCTCATCTTCTTCCAAGACCTGCACATCCGAGTATCCATGTTTCTAAAGGACAAAGTTTCAAATAAT
AACGGTCGCTTTGTGTTGCCGGTGTCCGGGCTGTTTCTTGGGAACTGAAGTTCAGGACTCATCAGA
ATGTTCAACAACAAGGTGAAGAAGTGAAGAGGATAGAATTCAAGCATGGTGGAACTATGTCCCTGCA
CCCAAAGAAGTTCTTTTGAACCTTATGGAGACCGAGTCTGAACTGGGAACTAACATGTACAGCGTG
AATCAGCCTGTGGAACTCATGTGTCTGGATCATCAAAGAACTAGCCTCATGGACCCAGGAAAGCATT
GCTCCAAACCCTCTTGTAAAGAAGAGCTGAATTTCTTGGCCAGGCTGATGGGAGGGATGGAGATTAAG
AAACCCAGTGGCCCTGAGCCCGGATTCCGGTTGAATCTTTACCACCGATGAAGAAGAGGAACAAGCA
GCGCTAACCAAGGCCAGAAGAGTTATCCTATGAAGTTATCAACATACAAGCCACCCAGGACCAGCAACGG
AGCGAGGAGCTGGCTCGAATCATGGGGAGTTTGAGATCACGGAGCAGCCAAGGCTGAGCACCAAGCAA
GGGACGATTTGCTCGCCATGATGGATGAGTTA TAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI



[View online »](#)

Plasmid Map:



ACCN: NM_145047

Insert Size: 1140 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_145047.4](#)

RefSeq Size: 1505 bp

RefSeq ORF: 1140 bp

Locus ID: 127700

UniProt ID: [Q8WVF1](#)

Cytogenetics: 1p34.3

MW: 43.3 kDa

Gene Summary: May be involved in drug clearance in the placenta.[UniProtKB/Swiss-Prot Function]
Transcript Variant: This variant (1) lacks an alternate in-frame exon compared to variant 3. The resulting isoform (1) has the same N- and C-termini but is shorter compared to isoform 3.