

Product datasheet for SC313281

MEIS2 (NM_170674) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	MEIS2 (NM_170674) Human Untagged Clone
Tag:	Tag Free
Symbol:	MEIS2
Synonyms:	CPCMR; HsT18361; MRG1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC313281 representing NM_170674. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTGTGAACCGTCAGAATTTTGTAAACGACTACTATAGGGCGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCGCAAAGGTACGATGAGCTGCCCCATTACGGCGGGATGGACGGAGTAGGGTTCCCGCTTCCATG
TACGGAGACCCTCACGCGCCGCGCCGATCCCCCGGTTACCCACCTGAACCACGGGCGCCGCTCCAC
GCCACACAGCACTACGGCGCGCACGCCCGCACCCCAATGTCATGCCGGCCAGTATGGGATCCGCTGTC
AACGACGCTTGAAGCGGGACAAGGACGCGATCTATGGGACCCGTTGTTTCTCTGTTAGCTCTGGTC
TTTGAGAAAGTGCAGCTGGCGACCTGCACTCCCCGGGAACCTGGAGTGGCTGGCGGAGACGTCTGCTCC
TCCGACTCCTTCAACGAGGACATCGCGGTCTTCGCCAAGCAGGTTCCGCGCGAAAAGCCACTTTTTTCC
TCAAATCCAGAGCTGGACAATTTGATGATACAAGCAATACAAGTACTAAGTTTCATCTTTTGGAGTTA
GAAAAGGTCCACGAAGTGTGCGATAACTTCTGCCACCGATACATTAGCTGTTTGAAGGGGAAAATGCC
ATCGACCTCGTCAATTGATGAAAGAGACGGCAGCTCCAAGTCAGATCATGAAGAACTTTCAGGCTCCTCC
ACAAATCTCGTGACCATAACCCCTTCTTGGCGAGACCACGATGATGCAACCTCAACCCACTCAGCA
GGCACCCAGGGCCCTCCAGTGGGGCCATGCTTCCAGAGCGGAGACAACAGCAGTGAAGGGGAT
GGTTTAGACAACAGTGTAGCTTACCTGGTACAGGTGACGATGATGATCCGGATAAAGGACAAAAACGC
CAGAAGAAAAGAGGCATTTTCCCAAAGTAGCAACAAATATCATGAGAGCATGGCTTCCAGCATCTC
ACACATCCGTACCCTTCCGAAGAGCAGAAGAAACAGTTAGCGCAAGACACAGGACTTACAATTTCCAA
GTAACAACCTGGTTTATTAATGCCAGAAGAAGAAATAGTACAGCCCATGATTGACCAGTCAAATCGAGCA
GTGAGCCAAGGAGCAGCATATAGTCCAGAGGGTACGCCATGGGGAGCTTTGTGTTGGATGGTCAGCAA
CACATGGGGATCCGGCTGCAGGACCTATGAGTGAATGGGCATGAATATGGGCATGGATGGGCAATGG
CACTACATGTAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: Sgfl-Mlul



[View online »](#)

Plasmid Map:


ACCN: NM_170674

Insert Size: 1185 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_170674.4](#)

RefSeq Size: 3352 bp

RefSeq ORF: 1185 bp

Locus ID: 4212

UniProt ID: [O14770](#)

Cytogenetics: 15q14

Domains:	homeobox
Protein Families:	Transcription Factors
MW:	43.1 kDa
Gene Summary:	<p>This gene encodes a homeobox protein belonging to the TALE ('three amino acid loop extension') family of homeodomain-containing proteins. TALE homeobox proteins are highly conserved transcription regulators, and several members have been shown to be essential contributors to developmental programs. Multiple transcript variants encoding distinct isoforms have been described for this gene. [provided by RefSeq, Jul 2008]</p> <p>Transcript Variant: This variant (b) uses an alternate in-frame splice junction at the 5' end of a coding exon compared to variant a. It encodes a shorter isoform (b) that is missing an internal segment compared to isoform a.</p>