

Product datasheet for **SC313317**

C2CD3 (AK092089) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	C2CD3 (AK092089) Human Untagged Clone
Tag:	Tag Free
Symbol:	C2CD3
Synonyms:	OFD14
Vector:	<u>pCMV6 series</u>
Fully Sequenced ORF:	>NCBI ORF sequence for AK092089, the custom clone sequence may differ by one or more nucleotides

```

ATGAAACAACGAAAAGGCCAAGGGTCTGGGGGCGAGCCGTGGGCGCAAAAAAGAGGTTTA
AGTGACATTTCTCCATCTACAAGCCTGCCACCTCTGGTTGAAGGCCAGCTACGCTGTTTT
CTAAAACCTTACTGTTAATAGAGTCATATGGAAGATTGCAAAGCCTCCCACTTGTGTACTT
GTCCGAGTGAGATGGTGGGGAGAAACATCAGATGGAACCTCTTTTGTCCAGGGATGCA
TTGCAGACTGAACCAAAAGCTGTGAGAACAACTACACGTTACGCTATTCTGTTGGTCCA
AAACAGTTTACCTCTTATCTAACAGATATGGCTGTGCTGGTGGTGAAGTAATCACCAAA
CTTGATGGTCTTCCAATTGGTAGAGTTCAGATCAATGGACTAGCTCAACTTTCTCCAACC
CATCAATCAATGGATTTTTTACCATTGTTTCATCAACGTCTAAGAACTTGGAGAACTC
CAGGTCTCACTTGCCCTGGAACCTCTGTGAGAACTTACGACAGCTACCATCCACTTCT
ACCACTGACATGACAGAAAATGTGCTTTTATCTAAGCAGGGATTGAGAGAGAATACTGAA
CCCAGCAGTACCCAGTTTCAGGTTCCATCAAGGCCTCGCGACATACATACCATCAAAAT
GATGGAAGAGAGTTAGCAGCCAAACAGCAGTAGATCAACCACTCCGAGGGGAAAAGACCAT
GTATGCTTTGAGAGAACCTGATACAATAAAGGATTCTTCTTTGGACTACAGCACAGT
CTTAATTCAGGACAGAGTTTAGAGTCTGTAACCTGAAAGGCAGAGCTCCACGGAAGCAG
ATGTCCCTTCTGAACAGTTCTGAATTCAGCCTCAAATTAGAACAGTTGCCAAGAGTCAC
AGTGACTCATGCAATTCTTCTTCAAACAACCTCCCTACCAAGGATCTTCTTTCAGCTCTG
TTAGAACAAGGCAATAAACTGCGTAATGCCATGGTGATTCTGCAATGAAATCAAGCCCA
GAGACCAGCATGTTGTTGGACCAAGTTCATCCTCTATTAATGAAGATTCTCTTAGAGCAT
CAACACAGATCAGAGCCTTTCTAGGAATCGGTTTAAAGACCACATTGAAGATCACCTCC
TCCCTTCAACTGAGAATACATTTTGAGACATGACACAAAAGCTGATACCAGAGCTATAC
AGCTGCTAT

```

Restriction Sites:	Please inquire
ACCN:	AK092089
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	AK092089.1 , BAC03802.1
RefSeq Size:	2586 bp
RefSeq ORF:	2586 bp
Locus ID:	26005
Cytogenetics:	11q13.4
Gene Summary:	This gene encodes a protein that functions as a regulator of centriole elongation. Studies of the orthologous mouse protein show that it promotes centriolar distal appendage assembly and is also required for the recruitment of other ciliogenic proteins, including intraflagellar transport proteins. Mutations in this gene cause orofaciodigital syndrome XIV (OFD14), a ciliopathy resulting in malformations of the oral cavity, face and digits. Alternative splicing of this gene results in multiple transcript variants. [provided by RefSeq, Nov 2014]