

Product datasheet for SC313326

GLEPP1 (PTPRO) (NM_030671) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	GLEPP1 (PTPRO) (NM_030671) Human Untagged Clone
Tag:	Tag Free
Symbol:	PTPRO
Synonyms:	GLEPP1; NPHS6; PTP-OC; PTP-U2; PTPROT; PTPU2; R-PTP-O
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC313326 representing NM_030671. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGTTACAGAGATGAATCCCAATGTGGTAGTGATCTCCGTGCTGGCCATCCTTAGCACACTTTTAATT
GGACTGTTGCTTGTACCTCATTATTCTTAGGAAAAAGCATCTGCAGATGGCTAGGGAGTGTGGAGCT
GGTACATTTGTCAATTTGCATCCTTAGAGAGGGATGAAAGCTCCATACAACTGGCGTAGGAGTATA
TTTGCTTTCTAACCTGCTACCCTCATGTCTTTGGACTGATTATCTTTGGCATTATATTAATCCT
TGGAGTAAAAATGGTTTAAAGAAGAGGAAACTGACAAACCCGGTCAACTGGATGACTTTGATGCCTAT
ATTAAGGATATGGCCAAAGACTCTGACTATAAATTTCTCTTCAGTTTGAGGAGTTGAAATTGATTGGA
CTGGATATCCCACTTTGCTGCAGATCTTCCACTGAATCGATGTAAAAACCGTTACACAAACATCCTA
CCATATGACTTCAGCCGTGTGAGATTAGTCTCCATGAATGAAGAGGAAGGTGCAGACTACATCAATGCC
AACTATATTCCTGGATACAACCTACCCAGGAGTATATTGCCACCCAGGGGCCACTGCCTGAAACCAGA
AATGACTTCTGGAAGATGGTCTGCAACAAAAGTCTCAGATTATTGTCATGCTCACTCAGTGAATGAG
AAAAGGAGGGTGAATGTGACCATTTACTGGCCATTACGGAAGAACCTATAGCCTATGGAGACATCACT
GTGGAGATGATTCAGAGGAAGAGCAGGACGACTGGGCCGTGAGACACTCCGGATCAACTATGCTGAC
GAGATGCAGGATGTGATGCATTTAACTACACTGCATGGCCTGATCATGGTGTGCCACAGCAAATGCT
GCAGAAAGTATCCTGCAGTTTGTACACATGGTCCGACAGCAAGCTACCAAGAGCAAAGGTCCCATGATC
ATTCAGTGCAGTGTGGCGTGGGACGGACAGGAACATTATTGCCCTGGACAGGCTCTTGACGACATT
CGGGATCATGAGTTTGTGACATCTTAGGGCTGGTGTGCAAAATGAGGTCATACCGGATGTCTATGGTA
CAGACAGAGGAGCAGTACATTTTATCCATCAGTGTGTGCAACTGATGTGGATGAAGAAGAAGCAGCAG
TTCTGCATCAGTGTGATATACGAGAATGTTAGCAAGTCTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: Sgfl-Mlul



[View online »](#)

Plasmid Map:


ACCN: NM_030671

Insert Size: 1218 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_030671.2](#)

RefSeq Size: 5002 bp

RefSeq ORF: 1218 bp

Locus ID: 5800

UniProt ID: [Q16827](#)

Cytogenetics: 12p13-p12

Protein Families: Phosphatase, Transmembrane

MW: 47.2 kDa

Gene Summary: This gene encodes a member of the R3 subtype family of receptor-type protein tyrosine phosphatases. These proteins are localized to the apical surface of polarized cells and may have tissue-specific functions through activation of Src family kinases. This gene contains two distinct promoters, and alternatively spliced transcript variants encoding multiple isoforms have been observed. The encoded proteins may have multiple isoform-specific and tissue-specific functions, including the regulation of osteoclast production and activity, inhibition of cell proliferation and facilitation of apoptosis. This gene is a candidate tumor suppressor, and decreased expression of this gene has been observed in several types of cancer. [provided by RefSeq, May 2011]

Transcript Variant: This variant (5) differs in the 5' UTR and 3' UTR, lacks a portion of the 5' coding region and initiates translation at a downstream, in-frame start codon, compared to variant 1. Variants 3 and 5 encode the same isoform (c), which has a shorter N-terminus compared to isoform a. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.