

Product datasheet for **SC313403**

JNK2 (MAPK9) (NM_139070) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	JNK2 (MAPK9) (NM_139070) Human Untagged Clone
Tag:	Tag Free
Symbol:	JNK2
Synonyms:	JNK-55; JNK2; JNK2A; JNK2ALPHA; JNK2B; JNK2BETA; p54a; p54aSAPK; PRKM9; SAPK; SAPK1a
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC (PS100020)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

Fully Sequenced ORF: >OriGene sequence for NM_139070 edited
 CCACGCGTCCGGAGGGATCTGAAACTTGCCACCCCTTCGGGATATTGCAGGACGCTGCAT
 CATGAGCGACAGTAAATGTGACAGTCAGTTTTATAGTGTGCAAGTGGCAGACTCAACCTT
 CACTGTCTAAAACGTTACCAGCAGCTGAAACCAATTGGCTCTGGGGCCCAAGGGATTGT
 TTGTGCTGCATTTGATACAGTCTTGGGATAAATGTTGCAGTCAAGAACTAAGCCGTCC
 TTTTCAGAACCAAACTCATGCAAAGAGAGCTTATCGTGAACCTTGCTCTTAAAAATGTG
 CAATCATAAAAATAAATTAGTTTGTAAATGTGTTTACACCACAAAAAACTCTAGAAGA
 ATTTCAAGATGTGATTTGGTTATGGAATTAATGGATGCTAACTTATGTCAGGTTATTCA
 CATGGAGCTGGATCATGAAAGAATGTCCTACCTCTTTACCAGATGCTTTGTGGTATTAA
 ACATCTGCATTACAGTGGTATAATTCATAGAGATTTGAAGCCTAGCAACATTGTTGTGAA
 ATCAGACTGCACCCTGAAGATCCTTGACTTTGGCCTGGCCCGGACAGCGTGCCTAACTT
 CATGATGACCCCTTACGTGGTGACACGGTACTACCGGGCGCCGAAGTCATCCTGGGTAT
 GGGCTACAAGAGAACGTTGATATCTGGTCAGTCGGGTGCATCATGGCAGAAATGGTCTC
 CCATAAAGTCTGTTCCCGGGAAGAGACTATATTGATCAGTGAATAAAGTTATTGAGCA
 GCTGGGAACACCATCAGCAGAGTTCATGAAGAACTTCAGCCAAGTGTGAGGAATTATGT
 CGAAAACAGACCAAAGTATCCTGGAATCAAATTTGAAGAAGTCTTCCAGATTGGATATT
 CCCATCAGAATCTGAGCGAGACAAAATAAAAAAAGTCAAGCCAGAGATCTGTTATCAAA
 AATGTTAGTGATTGATCCTGACAAGCGGATCTCTGTAGACGAAGCTCTGCGTCACCCATA
 CATCACTGTTTGGTATGACCCCGCCGAAGCAGAAGCCCCACCCTCAAATTTATGATGC
 CCAGTTGGAAGAAAAGAGAACATGCAATTGAAGAATGGAAGAGCTAATTTACAAAGAAGT
 CATGGATTGGGAAGAAAAGCAAGAATGGTGTGTAAGATGACGCTTCAGATGCAGC
 AGTAAGTAGCAACGCCACTCCTTCTCAGTCTTCATCGATCAATGACATTTATCCATGTC
 CACTGAGCAGACGCTGGCCTCAGACACAGACAGCAGTCTTGATGCCTCGACGGGACCCCT
 TGAAGGCTGTCGATGATAGGTTAGAAAATAGCAAACCTGTGAGCATTGAAGGAAGTCTCAC
 CTCGTTGGCCTGAAATGCTTGGGAGTTGATGGAACCAAATAGAAAAACTCCATGTTCTG
 CATGTAAGAAACACAATGCCTTGCCCTATTGACACCTGATAGGATTGCTGCTTAGATGA
 TAAAATGAGGCAGAAATGTCTGAAGAAAAAATTGCAAGCCACACTTCTAGAGATTTTG
 TTCAAGATCATTTCAGGTGAGCAGTTAGAGTAGGTGAATTTGTTTCAAATGTTACTAGTG
 ACAGTTTCTCATCATCTGTAAGTGTGAGATGTATGTGCATGTGACCACAAATGCTTGCT
 TGGACTTGCCCATCTAGCACTTTGGAAATCAGTATTTAAATGCCAAATAATCTTCCAGGT
 AGTGCTGCTTCTGAAGTTATCTTAACTCTTAAGTAAATTTGGTGTCTGTCCAGAAAA
 AGTCGATTTATGTGATTAATTGGCCATCATGATGTTATCATATCTTATCCCTTTTATG
 CTATGATTTATTCTATCTTTTGTATTTTCAAGACATATAATTAATCTATTTAATAAAT
 AAAAAATATATAGCTTTTCAAAAAAAAAAAAAAAAAA

Restriction Sites: Please inquire

ACCN: NM_139070

Insert Size: 1900 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: DNA sequence of MG2702_B09 matches to NM_139070.1.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_139070.1](#), [NP_620709.1](#)

RefSeq Size: 1942 bp

RefSeq ORF: 1275 bp

Locus ID: 5601

UniProt ID: [P45984](#)

Cytogenetics: 5q35.3

Domains: pkinase

Protein Families: Druggable Genome, ES Cell Differentiation/IPS, Protein Kinase

Protein Pathways: Adipocytokine signaling pathway, Colorectal cancer, Epithelial cell signaling in Helicobacter pylori infection, ErbB signaling pathway, Fc epsilon RI signaling pathway, Focal adhesion, GnRH signaling pathway, Insulin signaling pathway, MAPK signaling pathway, Neurotrophin signaling pathway, NOD-like receptor signaling pathway, Pancreatic cancer, Pathways in cancer, Progesterone-mediated oocyte maturation, RIG-I-like receptor signaling pathway, T cell receptor signaling pathway, Toll-like receptor signaling pathway, Type II diabetes mellitus, Wnt signaling pathway

Gene Summary:

The protein encoded by this gene is a member of the MAP kinase family. MAP kinases act as an integration point for multiple biochemical signals, and are involved in a wide variety of cellular processes such as proliferation, differentiation, transcription regulation and development. This kinase targets specific transcription factors, and thus mediates immediate-early gene expression in response to various cell stimuli. It is most closely related to MAPK8, both of which are involved in UV radiation induced apoptosis, thought to be related to the cytochrome c-mediated cell death pathway. This gene and MAPK8 are also known as c-Jun N-terminal kinases. This kinase blocks the ubiquitination of tumor suppressor p53, and thus it increases the stability of p53 in nonstressed cells. Studies of this gene's mouse counterpart suggest a key role in T-cell differentiation. Several alternatively spliced transcript variants encoding distinct isoforms have been reported. [provided by RefSeq, Sep 2008]

Transcript Variant: This variant (JNK2-b2) encodes the longer of the two JNK2 beta isoforms (JNK2 beta2). The JNK2-b2 variant differs from the JNK2-a2 variant in the use of an alternate internal coding exon of the same length. Thus, JNK2 beta2 isoform is the same length as JNK2 alpha2 isoform, with a few aa differences in an internal protein segment. Variants JNK-b2 and 9 both encode the same isoform (beta2). **Sequence Note:** This RefSeq record was created from transcript and genomic sequence data because no single transcript was available for the full length of the gene. The extent of this transcript is supported by transcript alignments.