

Product datasheet for **SC313423**

Nesprin 2 (SYNE2) (NM_182910) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Nesprin 2 (SYNE2) (NM_182910) Human Untagged Clone
Tag:	Tag Free
Symbol:	SYNE2
Synonyms:	EDMD5; KASH2; Nesp2; Nesprin-2; NUA; NUANCE; SYNE-2; TROPH
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC313423 representing NM_182910. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTGTAGAACCGTCAGAATTTTGTAAACGACTACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGAAGGTGACAGGAATGTTCCACCTGTTCCCCCTGCGTCCAGCACCCCTTATAAACACCCTATGGA
AAGCTACTATTACCTCCAGGCACGGATGGTGGCAAAGAAGGCCCGCGAGTCTGAATGGCAACCCACAG
CAGGAAGACGGGGACTGGCCGGTATCACAGAGCAGCAGTCCAGGTGCCTTCGACAGATGGGAGATGATT
CAAGCACAGGAGCTTACAATAAGCTCAAATAAAACAAAATTTGCAACAGCTGAACTCTGATATCAGC
GCCATCACTACTGGCTGAAAAAACTGAAGCAGAGCTGGAAATGTTAAAGATGGCAAAGCCTCCCTCT
GATATCCAGGAAATAGAAGTGAAGAGACTGCAGGAGATACTGAAAGCCTTTGACACTTACAAG
GCATTAGTGGTCTCTGTCAACGTGAGCAGCAAGGAATTTCTGAAACCGAGAGCCCCGAATCCACAGAG
CTCCAAAGTAGACTCCGCCAGCTGAGCCTGCTCTGGGAAGCAGCAGGGCGCAGTGGACAGCTGGAGA
GGGGGCTTACGACAGTCGCTCATGCAGTGCCAGGGCTCTAAGACACGGCCGCGCTCTGATGTGTTGTTT
TTAAGGACTTCCACCAGTTGAGTCAAAATCTGCTGCTGTGGTTAGCGAGTGCCAAGAACCGGAGGCAG
AAGGCTCATGTCACCGATCCAAAGGCAGACCCCGGGCTCTCCTAGAGTGTGGGAGGAACTAATGCAA
CTGAAAAGGAGCTGGTAGAACGTC AACCTCAAGTGGACATGTTACAGGAGATTTCAAACAGCCTTCTC
ATTAAGGGACATGGAGAAGACTGTATTGAAGCTGAAGAAAAGTGCATGTTATTGAGAAGAACTCAAA
CAGTTACGGGAGCAAGTGTCCCAAGATTTAATGGCCTTGCAAGGAAACCCAGAACCCAGCCTCACCCTG
CCAGCTTCGACGAGGTAGACTCGGGGGACAGCCTCCTGCAACATCCGTGCCAGCTCCCCGAGCAAAG
TTCAGAGCAGTGAGAACTACAGAAGCGAGGAGGACAGAGAGCAGGGTCCCCGGCAGCACACGGCCA
CAGCGCTCCTTCTCTCAAGGGTGGTCCGGGCAGCCCTACCCCTGCAGCTGCTCCTCCTGCTGCTGCTG
CTCCTGGCCTGCCTGCTGCCCTCCTCCGAAGAAGACTACAGCTGCACTCAGGCCAACAACTTTGCCCGG
TCCTTTTACCCCATGCTGAGGTACACCAATGGGCCACCCCCACATAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```



[View online >](#)

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_182910

Insert Size: 1290 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_182910.2](#)

RefSeq Size: 2586 bp

RefSeq ORF: 1290 bp

Locus ID: 23224

UniProt ID: [Q8WXH0](#)

Cytogenetics: 14q23.2

Protein Families: Druggable Genome, Transcription Factors

MW: 48.2 kDa

Gene Summary: The protein encoded by this gene is a nuclear outer membrane protein that binds cytoplasmic F-actin. This binding tethers the nucleus to the cytoskeleton and aids in the maintenance of the structural integrity of the nucleus. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Mar 2009]
Transcript Variant: This variant (2) differs in the 5' UTR and coding sequence compared to variant 5. The resulting isoform (2) is shorter at the N-terminus compared to isoform 5.