

Product datasheet for SC313436

NRXN3 (NM_138970) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	NRXN3 (NM_138970) Human Untagged Clone
Tag:	Tag Free
Symbol:	NRXN3
Synonyms:	C14orf60
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC313436 representing NM_138970. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTT TAGTGAACCGTCAGAATTTTGT AATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCC CGGATCGCC
ATGCACCTGAGAATCCACGCGAGACGGAGCCCTCTCGCCGGCCGGCCTGGACGCTTGGGATCTGGTTCC
CTGTTCTGGGGATGTATCGTCAGCTCTGTATGGAGTTCTTCTAATGTAGCTTCTCCTCCTCCACCTCT
TCCTCGCCGGGGTCTCACTCTCAGCAGGACCACTTTCCATGGCAGCAAGCATCACTCAGTGCCTATT
TCTATCTATCGTTCCTGTTCCCTTCGAGGAGGACACGCTGGCGCTACGTACATCTTTGGGAAAAGT
GGTGGGCTTATCCTCTACACCTGGCCAGCAATGACAGGCCAGCAGCGGTCTGACCGCTTGGCGTG
GGCTTCAGCACCCTGTGAAGGATGGCATCTTGGTCCGCATCGACAGTGTCCAGGACTTGGTGACTTC
CTCCAGCTTACATAGAACAGGGGAAAATTGGAGTTGTCTTCAACATTGGCAGATTGACATCTCCATC
AAAGAGGAGAGAACCCTGTAATGACGGCAAATACCATGTGGTACGCTTACCAGGAACGGCGGCAAC
GCCACCCTGCAGGTGGACAACCTGGCCAGTGAATGAACATTATCCTACAGGCCGGCAGTTAACCATCTTC
AACACTCAGGCGCAAATAGCCATTGGTGGAAAAGGACAAAGGACGCCTCTTCCAAGCCAACCTCTGGG
CTCTATTATGATGGTTTGGAAAGTACTGAACATGGCGGCTGAGAACAACCCCAATATTAATCAATGGA
AGTGTTCGGCTGGTTGGAGAAGTCCCATCAATTTGGGAACAACAGACGACCTCCATGCCACCAGAA
ATGTCTACTACTGTATGAAACCCTACTACAATGGCGACTACCACAACCCGTAAGAATCGCTCTACA
GCCAGATTGAGCAACATCAGATGATCTTGTTCATCTGCTGAATGTTCAAGTATGATGAAGACTTT
GTTGAATGTGAGCCGAGTACAGGTAGGTAGCAAAACCCACGGAGCCGGGAATCAGACGGGTTCCGGGG
GCCTCAGAGGTGATCCGGGAGTCGAGCAGCACAACAGGGATGGTCGTCGGCATTGTGGCTGCTGCCGCC
CTCTGCATCTTGATCCTCCTGTACGCCATGTACAAGTACAGGAACAGGGACGAGGGGTCCTATCAAGTG
GACGAGACGCGAACTACATCAGCAACTCCGCCAGAGCAACGGCAGCTCATGAAGGAGAAGCAGCAG
AGCTCGAAGAGCGGCCACAAGAAACAGAAAAACAAGGACAGGGAGTATTACGTGTAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```



[View online »](#)

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_138970

Insert Size: 1299 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_138970.4](#)

RefSeq Size: 8578 bp

RefSeq ORF: 1299 bp

Locus ID: 9369

UniProt ID:	Q9HDB5
Cytogenetics:	14q24.3-q31.1
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Cell adhesion molecules (CAMs)
MW:	47.3 kDa
Gene Summary:	<p>This gene encodes a member of a family of proteins that function in the nervous system as receptors and cell adhesion molecules. Extensive alternative splicing and the use of alternative promoters results in multiple transcript variants and protein isoforms for this gene, but the full-length nature of many of these variants has not been determined. Transcripts that initiate from an upstream promoter encode alpha isoforms, which contain epidermal growth factor-like (EGF-like) sequences and laminin G domains. Transcripts initiating from the downstream promoter encode beta isoforms, which lack EGF-like sequences. Genetic variation at this locus has been associated with a range of behavioral phenotypes, including alcohol dependence and autism spectrum disorder. [provided by RefSeq, Dec 2012]</p> <p>Transcript Variant: This variant (2) differs in the 5' UTR and contains multiple differences in the coding region, including the lack of multiple 5' exons, compared to variant 1. It initiates translation at an alternate start codon. The encoded isoform (2) is shorter and has a distinct N-terminus, compared to isoform 1. This variant encodes a beta isoform. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>