

Product datasheet for SC313471

RHBDL1 (NM_003961) Human Untagged Clone

Product data:

| | |
|---------------------------|---|
| Product Type: | Expression Plasmids |
| Product Name: | RHBDL1 (NM_003961) Human Untagged Clone |
| Tag: | Tag Free |
| Symbol: | RHBDL1 |
| Synonyms: | RHBDL; RRP |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pCMV6-Entry (PS100001) |
| E. coli Selection: | Kanamycin (25 ug/mL) |
| Fully Sequenced ORF: | >SC313471 representing NM_003961. Blue=Insert sequence Red=Cloning site Green=Tag(s) |


```
GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAATACGACTACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGGTAGGGTGAAGACGGGGAACTGAGGAGCTGGAGGACTGGGACCCAGGCACCAAGTGCCTG
CCAGCTCCTGGGATCAAGCAGGGTCCCAGGGAACAGACAGGCACGGGGCCCTGTCCAAAAGTGTGG
GAGCCTGAGCCTGATGCTCCAGCCAGCCTGGCCAGCCCTTTGGTCCAGGGGTGGGGCCGCACTCAG
GCCTTGGCTGGCGCTCCTCACTGCAGCAGCTGGACCCGAGAACACAGGCTTCATCGGTGCGGACACC
TTCCTGGCTGGTGCACAGCCATGAGCTGCCCTGGACCCGGAAGCTGGACATGCTGGTGGCCCTG
GCTCAGAGCAACGAGCAGGGCCAGTCTGCTACCAGGAGCTGGTGGACCTGATCAGCAGCAAGCGCTCC
AGCAGTTTCAAGCGGGCCATTGCTAACGGACAGCGGGCACTGCCCGGGACGGGCCGCTGGATGAGCCA
GGCCTAGGTGTCTACAAGCGTTTGTGCGTTACGTGGCCTACGAGATCCTGCCTTGTGAGGTGGACCGC
CGCTGGTACTTCTACCGTCAACCGAGCTGCCACCCCGTGTTCATGGCCTCGGTCACTTTGCCAG
ATCATCGTGTTCCTGTGTACGGGGCCCGCTCAACAAGTGGGTGCTGCAGACCTACCACCCGAGTAC
ATGAAGAGCCCCCTTGTGTACCACCCGGGACCGTGGCCGCGCTGGCGCTTCTCACCTACATGTTT
ATGCACGTTGGGCTGGAGCAGCTGGGGTCAACGCCCTCCTGCAGCTGATGATCGGGGTGCCCTGGAG
ATGGTGCACGGCTGCTCCGATCAGCCTGCTCTACCTGGCAGGCGTCTGGCAGGCTCCCTAACCGTC
TCCATCACCGACATGCGGGCCCGGTGGTGGGAGGCTCCGGCGGGGTCTACGCCCTGTGCTCGGCACAC
CTGGCCAACGTTGTGATGAAGTGGCTGGGATGAGATGTCCCTACAAGTTGCTGAGGATGGTGTGGCC
TTGGTGTGCATGAGCTCCGAGGTGGGCCGGGCGTGTGGCTGCGCTTCTCCCGCCGCTGCCCGCTCG
GGCCACAGCCAGCTTTCATGGCGACCTGGCAGGCGCGGTGGTGGGGTGGAGCATGGGCTGACCATC
CTGCGGAGCTACGAGGAGCGCTGCGGGACAGTGCAGGCTGGTGGGTGGTGTGCTGCTGACGACCC
TTCCTGCTTTCGCGCTTCTGGAACGCTTTCGCTACGACCTGCTGGGCGCCACATCCCCCACC
CCCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTAAACGGCCGCGC
```



[View online »](#)

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_003961

Insert Size: 1317 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_003961.1](#)

RefSeq Size: 1559 bp

RefSeq ORF: 1317 bp

Locus ID: 9028

Cytogenetics: 16p13.3

Protein Families: Protease, Transmembrane

MW: 48.3 kDa

Gene Summary: This gene encodes a protein similar to Rhomboid in *Drosophila* which is involved in signalling in the Spitz/epidermal growth factor receptor/mitogen-activated protein kinase pathway. The Rhomboid family of proteins consists of intramembrane serine proteases containing several transmembrane domains. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Jul 2013]