

Product datasheet for SC313490

CCDC91 (NM_018318) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CCDC91 (NM_018318) Human Untagged Clone
Tag:	Tag Free
Symbol:	CCDC91
Synonyms:	HSD8; p56
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC313490 representing NM_018318. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTAGTGAACCGTCAGAATTTGTAAACGACTACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGATGATGATGATTTTGGTGGTTTTGAGGCTGCGGAGACTTTTGATGGTGGAAGTGGTAAACCCAA
ACAACATCTCCTGCTATTCCTGGGCTGCCTTTCCTGCAGTATCTGGAGTCCATCTTTCACCATCTTCT
CCTGAGATTGACTGGACCGTGACCCTTCTTCCATTGGCTGCCTCTTCTGATGCCATTATTTCA
TCACCAGAGAATACACATGCAGCAAATAGCATTGTGAGTCAAATATTCCAAAAGCACAGATTACAGCA
TCAACACACACTCATCTGGATATCTCACTTTTTCCATTGGGTTAACTGATGAAAAAGTAAATGGAACA
ATTGCCCTTGTGGATGATTCTGAGGATCCTGGAGCCAATGTATCTAACATACAGTTCAGCAAAAAATT
TCAAGTCTGGAGATTAAGTCAAAGTATCTGAAGAAGAAAAACAGAGAATTAACAGGATGTGGAATCA
TTGATGAAAAAGCATAATGTCTTAGAAAAAGGCTTTCTAAAAGAAAAAGCAAGAGGCCATTTCTTTT
CAAGATAGATACAAAGAACTTCAGGAAAAACATAAAACAAGAATTGGAAGACATGAGGAAAGCTGGTAC
GAAGCCCTCAGCATTATTGTGGATGAATATAAGGCACTACTGCAGTCTTCAGTTAAGCAACAAGTAGAA
GCTATTGAAAAACAGTACATTTCTGCAATTGAGAAACAGGCACACAAGTGTGAGGAGTTGCTAAATGCT
CAGCATCAGAGGCTCCTTGAATGCTAGATACAGAGAAGGAAGTGTAAAAGAAAAATAAAGGAAGCT
TTGATTCAGCAATCTCAAGAACAGAAGGAAATATTGAAAAAGTGTGGAGGAAGAAAGGCAAGAAAT
AAAGAGGCATTAGTATCCGCTGCAAAGCTTGAAGAAAGCAAGTGAAGGATGCAGTTTTAAAGTTCGTA
GAAGAAGAAAGAAAAATTTAGAAAAAGCGCATGCTGAAGAAAGGGAATTAAGGAAGACAGAATCATGCA
AAAGATCAAGAAAAAGTATCTCAGGAAATTCAAAAAGCTATACAAGAACAAAGAAAAATAAGTCAGGAA
ACTGTTAAGGCAGCAATAATAGAAGAGCAGAAACGAAGTGAAGGCTGTGGAAGAGGCAGTGAAGAA
ACAAGAGATGAATTGATAGAGTATATAAAGAACAGAAAGGCTCGATCAAGTCAATCCGCCAAAGAGC
CTGTCCAGTTTGGAACTGTTCTCTCCTGTGCACAGAAACAGTTAAGTCTTTAATAGCTACGGAACCA
GTTGACATTGAATAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```



[View online »](#)

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_018318

Insert Size: 1326 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_018318.3](#)

RefSeq Size: 2352 bp

RefSeq ORF: 1326 bp

Locus ID: 55297

UniProt ID: [Q7Z6B0](#)

Cytogenetics: 12p11.22

MW: 50 kDa

Gene Summary: Involved in the regulation of membrane traffic through the trans-Golgi network (TGN). Functions in close cooperation with the GGAs in the sorting of hydrolases to lysosomes. [UniProtKB/Swiss-Prot Function]
Transcript Variant: This variant (1) encodes the longest isoform (1). Variants 1, 3, 4, 5, and 6 all encode the same isoform (1).