

## Product datasheet for SC313492

### EVA1C (NM\_058187) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	EVA1C (NM_058187) Human Untagged Clone
Tag:	Tag Free
Symbol:	EVA1C
Synonyms:	B18; B19; C21orf63; C21orf64; FAM176C; PRED34; SUE21
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC313492 representing NM_058187. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCTTCTGCCGGGACGCGCACGCCAACCGCCGACGCCAGCCCGTGCAGCATCCCGCCCTCCGCCG
CAGGTAGAGCCGCGGGGAGCTCTGCCCTCTTCTACTGACTGTCCTGGTCTGCTCAAAGAGATC
TCAGCGCTCACCAGCTTCTCTGGTTACCTAACCAACTCTGCAAAACCACACCCTATGCCTGTGAT
GGGGACTATTTGAATCTACAGTGCCTCGGCATTCTACGATAAGTGTCCAATCGGCATTTTATGGCAA
GATTACCAATGTGTAGTCCCAGAAGCCTGCCTCCAGAGGGAAGACAGCTTAACCTGTGTGCCAGCC
ACCACCTCCAGAAGGTGCTGGACGAATGCCAGAACCAGCGGGCCTGCCACCTCCTGGTCAATAGCCGT
GTTTTGGACCTGACCTTTGTCCAGGAAGCAGTAAATACCTCCTGGTCTCCTTTAAATGCCAACCTAAT
GAATTA AAAAACAAAACCGTGTGTGAAGACCAGGAGCTGAAACTGCACTGCCATGAATCCAAGTTCTC
AACATCTACTCTGCGACCTACGGCAGGAGGACCCAGGAAAGGGACATCTGCTCCTCCAAGGCAGAGCGG
CTCCCCCTTTTCGATTGCTTGTCTTACTCAGCTTTGCAAGTCTATCCCGAAGGTGCTATGGGAAGCAG
AGATGCAAAATCATCGTCAACAATCACCATTTGGGAAGCCCTGTTTGCCAGGCGTGAAAAATACCTC
ACTGTGACCTACGCATGTGTTCCCAAGAACATACTCACAGCGATTGATCCAGCCATTGCTAATCTAAAA
CCTTCTTTGAAGCAGAAAGATGGTGAATATGGTATAAACTTCGACCAAGCGGATCGAAGGTTCTGAGG
AAAGATGGAATTCTTGTAGCAACTCTCTGGCAGCCTTTGCTTACATTAGAGCCACCCGGAGAGAGCT
GCCCTGCTGTTTCGTGCCAGTGTCTGCATCGGCCCTGGCCCTCACACTGTGCGCCCTGGTCAATCAGAGAG
TCCTGTGCCAAGGACTTCCGCGACTTGCAGCTGGGGAGGAGCAGCTGGTGCCAGGAAGTACAAGGTC
GAGGAGGACAGCGAGGATGAAGAAGAGGAGGAGGACCCCTCTGAGTCTGATTTCCAGGGAACTGTGCG
GGGTTCTGTAGGACTTCATATCTATATACAGTTCATAGAAGCTGCAGAGCTCGCAGAAAGGATTGAG
CGCAGGGAGCAAATCATTAGGAAATATGGATGAACAGTGGTTTGGACACCTCGTCCCAAGAAACATG
GGCCAGTTCTACTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```



[View online »](#)

<b>Restriction Sites:</b>	Sgfl-Mlul
<b>ACCN:</b>	NM_058187
<b>Insert Size:</b>	1326 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<u><a href="#">NM_058187.4</a></u>
<b>RefSeq Size:</b>	1965 bp
<b>RefSeq ORF:</b>	1326 bp
<b>Locus ID:</b>	59271
<b>UniProt ID:</b>	<u><a href="#">P58658</a></u>
<b>Cytogenetics:</b>	21q22.11
<b>Domains:</b>	Gal_Lectin
<b>Protein Families:</b>	Transmembrane
<b>MW:</b>	49.5 kDa
<b>Gene Summary:</b>	<p>Binds heparin.[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (1) represents the predominant transcript and encodes the longer isoform (1).</p>