

Product datasheet for **SC313510**

ETEA (FAF2) (BC014001) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ETEA (FAF2) (BC014001) Human Untagged Clone
Tag:	Tag Free
Symbol:	FAF2
Synonyms:	ETEA; expressed in T-cells and eosinophils in atopic dermatitis; Fas associated factor family member 2; KIAA0887; UBXD8; UBX domain containing 8; UBX domain protein 3B; UBXN3B
Vector:	<u>pCMV6 series</u>
Fully Sequenced ORF:	>NCBI ORF sequence for BC014001, the custom clone sequence may differ by one or more nucleotides ATGGCGGCGCCTGAGGAGCGGGATCTAACCCAGGAGCAGACAGAGAAGCTGCTGCAGTTT CAGGATCTCACTGGCATCGAATCTATGGATCAGTGTGCCATACCTTGGAACAGCATAAC TGGAACATAGAGGCTGCTGTACAGGACAGATTGAATGAGCAAGAGGGCGTACCTAGTGTT TTCAACCCACCTCCATCACGACCCCTGCAGGTTAATACAGCTGACCACAGGATCTACAGC TATGTTGTCTCAAGACCTCAACCAAGGGGGCTGCTTGGATGGGGTTATTACTTGATAATG CTTCCATTCCGGTTTACCTATTACAGATACTTGATATATTTAGGTTTGCCTTCGTTTTT ATACGGCCTGACCCTCGCAGCCGGGTCCTGACCCCGTTGGGGACATTGTTTCATTTATG CACTCTTTTGAAGAGAAATATGGGAGGGCACACCCTGTCTTCTACCAGGGAACGTACAGC CAGGCACTTAACGATGCCAAAAGGGAGCTTCGCTTTCTTTTGGTTTATCTTCATGGAGAT GATCACCAGGACTCTGATGAGTTTTGTGCAACACACTCTGTGCACCTGAAGTTATTTCA CTAATAAACACTAGGATGCTCTTCTGGGCATGCTCTACAAACAAACCTGAGGGATACAGG GTCTCACAGGCTTTACGAGAGAACACCTATCCATTCTGGCCATGATTATGCTGAAGGAT CGAAGGATGACTGTGGTGGGACGGCTAGAAGGCCTCATTCAACCTGATGACCTCATTAAAC CAACTGACATTTATCATGGATGCTAACCCAGACTTACCTGGTGTGAGAACGCCTAGAAAGG GAAGAAAGAAACCAGACCCAAGTCTGAGACAACAGCAGGATGAGGCCTACCTGGCCTCT CTCAGAGCTGACCAGGAGAAAAGAAAAGAAAACGGGAGGAGCGGGAGCGTAAGCGGCGG AAGGAGGAGGAGGTGCAACAGCAAAAAGTTGGCAGAGGAGAGACGGCGGCAGAAATTTACAG GAGGAAAAGGAAAAGGAAGTTGGAATGCCTGCCCCCTGAACCTTCCCTGATGACCCTGAA AGTGTC AAGATCATCTTCAAATTACCTAATGATTCTCGAGTAGAGAGACGATTCCACTTT TCACAGTCTTAACAGTAATCCACGACTTCTTATTCTCCTTGAAGGAAAGCCAGAAAAG TTTCAGATTGAAGCCAATTTTCCAGGCGAGTGCTGCCCTGCATCCCTTCCAGAGGAGTG CCCAATCCCCCTACGCTACAGGAGGCCGGACTCAGCCACACAGAAGTTCTTTTTGTTTCAG GACCTAACTGACGAA
Restriction Sites:	Please inquire
ACCN:	BC014001



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>BC014001.2</u> , <u>AAH14001.2</u>
RefSeq Size:	2059 bp
RefSeq ORF:	1338 bp
Locus ID:	23197
Cytogenetics:	5q35.2
Domains:	UBX, UAS
Gene Summary:	The protein encoded by this gene is highly expressed in peripheral blood of patients with atopic dermatitis (AD), compared to normal individuals. It may play a role in regulating the resistance to apoptosis that is observed in T cells and eosinophils of AD patients. [provided by RefSeq, Jul 2008]