

## Product datasheet for SC313520

### PPP2R2C (NM\_020416) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	PPP2R2C (NM_020416) Human Untagged Clone
Tag:	Tag Free
Symbol:	PPP2R2C
Synonyms:	B55-GAMMA; B55gamma; IMYPNO; IMYPNO1; PR52; PR55G
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC313520 representing NM_020416. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTAGTGAACCGTCAGAATTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGGCGAGGACACGGACACGCGGAAAATTAACCACAGCTTCTGCGGGACCACAGCTATGTGACTGAA
GCTGACATCATCTCTACCGTTGAGTTCAACCACACGGGAGAGCTGCTGGCCACAGGTGACAAGGGCGGC
CGGGTCGTATCTCCAGCGGGAACAGAGAGTAAAAATGCGCCCCACAGCCAGGGCGAATACGACGTG
TACAGCACTTTCCAGAGCCACGAGCCGGAGTTTGACTATCTCAAGAGCCTGGAGATAGAGGAGAAGATC
AACAGATCAAGTGGCTCCACAGCAGAACGCCGCCACTCACTCCTGTCCACCAACGATAAAACTATC
AAATTATGGAAGATTACGAACGAGATAAAAGGCCGAAGGATACAACCTGAAGGATGAAGAGGGGAAA
CTTAAGGACCTGTCCACGGTACGTCAGGTGCCAGTGCTGAAGCCATGGATCTGATGGTGGAG
GTGAGCCCTCGGAGGATCTTTGCCAATGGCCACACCTACCACATCAACTCCATCTCCGTCAACAGTGAC
TGCGAGACCTACATGTCCGGGATGACCTGCGCATCAACCTCTGGCACCTGGCCATCACCACAGGAGC
TTCAACATCGTGGACATCAAGCCGGCAACATGGAGGACCTTACGGAGGTGATCAGCATCTGAGTTC
CATCCGCACCACTGCAACCTCTTCGTCTACAGCAGCAGCAAGGGCTCCCTGCGGCTCTGCGACATGCGG
GCAGTGCCCTGTGTGACAAGCATTCCAAGCTCTTTGAAGAGCCTGAGGACCCAGTAACCGCTCATTC
TTCTCGGAAATCATCTCCTCCGTGTCCGACGTGAAGTTCAGCCACAGCGGCCGCTACATGCTCACCAGG
GACTACCTTACAGTCAAGTCTGGGACCTGAACATGGAGGCAAGCCATAGAGACCTACCAGTCCAT
GACTACCTTCGGAGCAAGCTCTGTTCCCTGTACGAGAACGACTGCATTTTCGACAAGTTTGAATGTGCC
TGGAACGGGAGCGACAGCGTCATCATGACCGGGGCTACAACAACCTTCTCCGCATGTTTCGATCGGAAC
ACCAAGCGGGACGTGACCCTGGAGGCTCGAGGAAAAGCAGCAAGCCCCGGGCTGTGCTCAAGCCACGG
CGCGTGTGCGTGGGGGCAAGCGCCGGCGTGATGACATCAGTGTGGACAGCTTGGACTTCAACAAGAAG
ATCCTGCACACGGCTGGCACCCGGCTGAGAACATCATTGCCATCGCCGCCACCAACAACCTGTACATC
TCCAGGACAAGGTAACCTCTGACATGCACTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```



[View online >](#)

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM\_020416

Insert Size: 1344 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM\\_020416.3](#)

RefSeq Size: 4438 bp

RefSeq ORF: 1344 bp

Locus ID: 5522

UniProt ID: [Q9Y2T4](#)

**Cytogenetics:** 4p16.1

**Protein Families:** Druggable Genome, Phosphatase

**Protein Pathways:** Tight junction

**MW:** 51.5 kDa

**Gene Summary:** The product of this gene belongs to the phosphatase 2 regulatory subunit B family. Protein phosphatase 2 is one of the four major Ser/Thr phosphatases, and it is implicated in the negative control of cell growth and division. It consists of a common heteromeric core enzyme, which is composed of a catalytic subunit and a constant regulatory subunit, that associates with a variety of regulatory subunits. The B regulatory subunit might modulate substrate selectivity and catalytic activity. This gene encodes a gamma isoform of the regulatory subunit B55 subfamily. Alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Jul 2008]  
Transcript Variant: This variant (1) encodes isoform a. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.