

Product datasheet for **SC313529**

TEA domain family member 2 (TEAD2) (AK093236) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TEA domain family member 2 (TEAD2) (AK093236) Human Untagged Clone
Tag:	Tag Free
Symbol:	TEA domain family member 2
Synonyms:	ETF; TEAD-2; TEF-4; TEF4
Vector:	<u>pCMV6 series</u>
Fully Sequenced ORF:	>NCBI ORF sequence for AK093236, the custom clone sequence may differ by one or more nucleotides ATGGGGGAACCCCGGGCTGGGGCCGCCCTGGACGATGGCAGCGGCTGGACGGGCAGTGAG GAAGGCAGTGAGGAGGGTACCGCGGCAGTGAGGGGGCTGGGGGTGACGGGGGCCCGGAT GCAGAGGGGGTGTGGAGCCAGACATTGAGCAGAGCTTCCAGGAGGCCCTGGCCATCTAT CCACCCTGCGGCCCGGAAAATAATTTTGTCTGATGAAGCAAGATGTATGGTCGGAAT GAACTGATCGCCGCTACATCAAGCTGAGAACGGGAAGACCCGAACCGAAAACAGGTT TCTAGTCACATCCAGTTTTGGCCGAAGAAAATCAAGGGAAATCCAGTCCAAGTTGAAG GACCAGTTTTCCAAGACAAGGCTTTCCAGACAATGGCAACCATGTCCTCTGCCAGCTC ATCTCCGCGCTTCTCTGCAGGCCAAACTGGGTCCCACCTGGTCCCTCAGGTGGTCCAGGCC TCTGAGCTTTTCCAGTTTTGGTCTGGAGGATCTGGGCCCCCTGGAATGTTCCAGATGTG AAGCCATTCTCACAGACACCGTTACCTTGTCACTGACTCCCCATCTACTGACCTCCCA GGGTACGAGCCCCCAAGCCCTCTCACCCCTGCCCCACCTACCCCATCGCCCCAGCC TGGCAGGCTCGGGGCTGGGCACCGCCCGGTTGCAGCTGGTAGAGTTCTCAGCCTTCGTG GAACCGCCAGATGCAGTTGATTCTTACCAGAGGCACCTGTTCTGTGCACATCAGCCAGCAC TGCCCCAGCCCCGAGCGCCCGCTCGAGAGTGTGGACGTCCGGCAGATCTACGACAAA TTCCCTGAGAAAAGGGCGGCCCTCCGAGAGCTATATGATCGTGGCCCCCATGCCTTC TTCCTGGTCAAGTTCTGGGCGGACCTGAACTGGGGCCCAAGTGGTGAAGGAGGCAGGGGCC GGTGGCAGCATCAGCAGTGGTGGCTTCTACGGAGTGAGCAGCCAGTATGAGAGCCTGGAA CACATGACCTCACCTGTTCTCCAAGGTCTGCTCTTTTGGCAAGCAGGTGGTGGAGAAG GTGGAGACGGAACGGGCCAGCTGGAGGACGGCAGATTTGTGTACCGCTGCTGCGCTCG CCCATGTGCGAGTACCTGGTGAATTTCTTGACAAAGTTGCGGCAGCTGCCTGAGCGATAC ATGATGAACAGCGTCTGGAAAACCTTACCATCCTCCAGGTGGTGACAAACAGAGACACC CAGGAACTGCTGCTCTGCACCGCTATGTCTTCGAGGTCTCCACCAGCGAGCGTGGGGCC CAGCATCACATTTACCGCTGGTCAGGGAC
Restriction Sites:	Please inquire
ACCN:	AK093236



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>AK093236.1</u> , <u>BAC04104.1</u>
RefSeq Size:	2442 bp
RefSeq ORF:	2442 bp
Locus ID:	8463
Cytogenetics:	19q13.33
Domains:	TEA
Protein Families:	Transcription Factors
Gene Summary:	Transcription factor which plays a key role in the Hippo signaling pathway, a pathway involved in organ size control and tumor suppression by restricting proliferation and promoting apoptosis. The core of this pathway is composed of a kinase cascade wherein MST1/MST2, in complex with its regulatory protein SAV1, phosphorylates and activates LATS1/2 in complex with its regulatory protein MOB1, which in turn phosphorylates and inactivates YAP1 oncoprotein and WWTR1/TAZ. Acts by mediating gene expression of YAP1 and WWTR1/TAZ, thereby regulating cell proliferation, migration and epithelial mesenchymal transition (EMT) induction. Binds to the SPH and GT-IIC 'enhansons' (5'-GTGGAATGT-3'). May be involved in the gene regulation of neural development. Binds to the M-CAT motif. [UniProtKB/Swiss-Prot Function]