

## Product datasheet for **SC313607**

### ENAH (AK096246) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	ENAH (AK096246) Human Untagged Clone
Tag:	Tag Free
Symbol:	ENAH
Synonyms:	ENA; MENA; NDPP1
Vector:	<u>pCMV6 series</u>
Fully Sequenced ORF:	>NCBI ORF sequence for AK096246, the custom clone sequence may differ by one or more nucleotides ATGATGCATGCCTTAGAAGTGTTAAATTCACAGGAAACAGGGCCAACATTGCCTAGACAA AACTCACAACTACCTGCTCAAGTTCAAAATGGCCCATCCAAGAAGAATTGGAAATTCAA AGAAGACAACACAAGAACAGCAACGGCAAAAGGAGCTGGAGCGGAAAGGCTGGAGCGA GAAAGAATGGAAGAGAAAGGTTGGAGAGAGAGAGGTTAGAAAGGAAAGGCTGGAGAGG GAGCGACTGGAACAAGAACAGCTGGAGAGAGAGAGACAAGAACGGGAACGGCAGGAACGC CTGGAGCGGCAGGAACGCCTGGAGCGGCAGGAACGCCTGGAGCGGCAGGAACGCCTGGAT CGGGAGAGGCAAGAAAGACAAGAACGAGAGAGGCTGGAGAGACTGGAACGGGAGAGGCAA GAAAGGGAGCGACAAGAGCAGTTAGAAAGGGAACAGCTGGAATGGGAGAGAGAGCGCAGA ATATCAAGTGCTGCTGCCCTGCCTCTGTTGAGACTCCTCTAAACTCTGTGCTGGGAGAC TCTTCTGCTTCTGAGCCAGGCTTGCAGGCAGCCTCTCAGCCGGCCGAGACTCCATCCCAA CAGGGCATTGTCTTGGGACCACTTGACCTCCACCTCCTCCACCACTCCCACCAGGGCCT GCACAGGCTTCAGTAGCCCTCCCTCCTCCCCAGGGCCCCCTCCACCTCCTCCACTCCCA TCCACCGGGCCTCCACCGCCCCCTCCTCCCCCTCCTCCTCAATCAAGTACCCCTCCT CCTCCACCACTCCTGCCCCACCCTCCCTGCATCTGGATTCTTTTTGGCATCCATGTCA GAAGACAATCGCCCTTAACTGGACTTGACAGCTGCAATTGCCGGAGCAAACTTAGGAAA GTGTCACGGATGGAGGATACCTCTTCCCAAGTGGAGGGAATGCTATTGGTGTGAACTCC GCCTCATCTAAAACAGATACAGGCCGTGGAATGGACCCCTCCTTTAGGGGTAGTGGT TTAATGGAAGAAATGAGTGCCCTGCTGGCCAGGAGGAGAAGAATTGCTGAAAAGGGATCA ACAATAGAAACAGAACAAAAAGAGGACAAAGGTGAAGATTAGAGCCTGTAACCTTCTAAG GCCTCTTCAACAAGTACACCTGAACCAACAAGAAAACCTTGGGAAAGAACAATAACAATG AATGGCAGCAAGTACCTGTTATCTCCAGACCAAAATCCACACCCTATCACAGCCCAGT GCCAATGGAGTCCAGACGGAAGGACTTGACTATGACAGGCTGAAGCAGGACATTTTAGAT GAAATGAGAAAAGAAATTAACAAAGCTAAAAGAAGAGCTCATTGATGCAATCAGGCAGGAA CTGAGCAAGTCAATACTGCA
Restriction Sites:	Please inquire
ACCN:	AK096246



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<a href="#">AK096246.1</a> , <a href="#">BAC04736.1</a>
<b>RefSeq Size:</b>	3759 bp
<b>RefSeq ORF:</b>	3759 bp
<b>Locus ID:</b>	55740
<b>Cytogenetics:</b>	1q42.12
<b>Protein Pathways:</b>	Regulation of actin cytoskeleton
<b>Gene Summary:</b>	This gene encodes a member of the enabled/ vasodilator-stimulated phosphoprotein. Members of this gene family are involved in actin-based motility. This protein is involved in regulating the assembly of actin filaments and modulates cell adhesion and motility. Alternate splice variants of this gene have been correlated with tumor invasiveness in certain tissues and these variants may serve as prognostic markers. A pseudogene of this gene is found on chromosome 3. [provided by RefSeq, Sep 2016]