

## Product datasheet for **SC313635**

### **BSDC1 (NM\_018045) Human Untagged Clone**


#### Product data:

Product Type:	Expression Plasmids
Product Name:	BSDC1 (NM_018045) Human Untagged Clone
Tag:	Tag Free
Symbol:	BSDC1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC313635 representing NM_018045. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCGGAAGGGGAGGACGTGGGATGGTGGCGGAGCTGGCTGCAGCAGAGCTACCAAGCAGTCAAAGAG
AAGTCCTCTGAAGCCTTGGAGTTTATGAAGCGGGACCTGACGGAGTTTACCCAGGTGGTGCAGCATGAC
ACGGCCTGTACCATCGCAGCCACGGCCAGCGTGGTCAAGGAGAAGCTGGCTACGGAAGGCCTCTCAGGA
GCAACAGAGAAGATGAAGAAAGGGTTATCTGACTTCTAGGGGTGATCTCAGACACCTTTGCCCTTCG
CCAGACAAAACCATCGACTGCGATGTATCACCCTGATGGGCACACCGTCTGGCAGAGCTGAGCCCTAT
GATGGCACCAAGGCTCGCCTCTATAGCCTGCAGTCGGACCCAGCAACCTACTGTAATGAACAGATGGG
CCCCCGAATTGTTTGACGCCTGGCTTTCCAGTTCTGCTTGGAGGAGAAGAAGGGGAGATCTCAGAG
CTCCTTGTAGGCAGCCCTCCATCCGGGCCCTCTACACCAAGATGGTTCCAGCAGCTGTTTCCATTCA
GAATTCTGGCATCGGTATTTCTATAAAGTCCATCAGTTAGAGCAGGAGCAGGCCCGGAGGGACGCCCTG
AAGCAGCGGGCGGAACAGAGCATCTCTGAAGAGCCCGCTGGGAGGAGGAGGAAGAGGAGCTCATGGGC
ATTTACCCATATCTCCAAAAGAGGCAAAGGTTCTGTGGCCAAAATTTCTACATTCCCTGAAGGAGAA
CCTGGCCCCCAGAGCCCCTGTGAAGAGAATCTGGTGACTTCAGTTGAGCCCCCAGCAGAGGTGACTCCA
TCAGAGAGCAGTGAGAGCATCTCCCTCGTGACACAGATCGCCAACCCGGCCACTGCACCTGAGGCACGA
GTGCTACCAAGGACCTGTCCAAAAGCTGCTAGAGGCATCCTTGAGGGAACAGGGCCTGGCTGTGGAT
GTGGGTGAGACTGGACCCTCACCCCTATTCACCTCAAGCCCTAACGCCTGCTGGCCACACCCGGCGGC
CCAGAGCCAGGCCTCCAGCCAGAGTAGAGACTCTGAGGGAGGAGCGCCACAGACTTACGGGTGTTT
GAGCTGAACTCGGATAGTGGGAAGTCTACACCCTCCAACAATGGAAGAAAGGCTCAAGCACGGACATC
AGTGAGGACTGGGAGAAAGACTTTGACTTGGACATGACTGAAGAGGAGGTGCAGATGGCACTTTCCAAA
GTGGATGCCTCCGGGAGCTGGAAGATGTAGAGTGGGAGGACTGGGAGTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: Sgfl-Mlul



**Plasmid Map:**


**ACCN:** NM\_018045

**Insert Size:** 1293 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_018045.7](#)

**RefSeq Size:** 3329 bp

**RefSeq ORF:** 1293 bp

**Locus ID:** 55108

**UniProt ID:** [Q9NW68](#)

Cytogenetics: 1p35.1  
Domains: BSD  
MW: 47.2 kDa