

Product datasheet for **SC313662**

LBP (NM_004139) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	LBP (NM_004139) Human Untagged Clone
Tag:	Tag Free
Symbol:	LBP
Synonyms:	BPIFD2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


Fully Sequenced ORF: >SC313662 representing NM_004139.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGGGGCCTTGCCAGAGCCCTGCCGTCCATACTGCTGGCATTGCTGCTTACGCCACCCAGAGGCT
CTGGGTGCCAACCCCGCTTGGTCGCCAGGATCACCGACAAGGGACTGCAGTATGCCGCCAGGAGGGG
CTATTAGCTCTGCAGAGTGAGCTGCTCAGGATCAGCTGCCTGACTTCACCGGGACTTGAGGATCCCC
CACGTCGGCCGTGGGCGCTATGAGTTCCACAGCCTGAACATCCACAGCTGTGAGCTGCTTACTCTGCG
CTGAGGCCTGTCCCTGGCCAGGGCCTGAGTCTCAGCATCTCCGACTCCTCCATCCGGGTCCAGGGCAGG
TGAAGGTGCGCAAGTCATTCTTCAAACACAGGGCTCCTTTGATGTCAGTGTCAAGGGCATCAGCATT
TCGGTCAACCTCCTGTTGGGCAGCGAGTCTCCGGGAGGCCACAGTTACTGCCTCCAGCTGCAGCAGT
GACATCGCTGACGTGGAGGTGGACATGTCGGGAGACTTGGGTGGCTGTTGAACCTTTCCACAACCAG
ATTGAGTCCAAGTCCAGAAAGTACTGGAGAGCAGGATTTGCGAAATGATCCAGAAATCAGTGCCTCC
GATCTACAGCCTTATCTCAAACCTCGCCAGTTACAACAGAGATTGACAGTTTCGCCGACATTGATTAT
AGCTTAGTGGAAGCCCTCGGGCAACAGCCAGATGCTGGAGGTGATGTTAAGGGTGAAATCTTTCAT
CGTAACCACCGTTCTCCAGTTAACCTCCTTGCTGCAGTCATGAGCCTTCCTGAGGAACACAAACAAATG
GTCTACTTTGCCATCTCGGATTATGCTTCAACACGGCCAGCCTGGTTTATCATGAGGAAGGATATCTG
AACTTCTCCATCACAGATGACATGATACCGCCTGACTCTAATA TCCGACTGACCACCAAGTCTTCCGA
CCCTTCGTCCCACGGTTAGCCAGGCTCTACCCCAACATGAACCTGGAACCTCAGGGATCAGTGCCCTCT
GCTCCGCTCCTGAACTTCAGCCCTGGGAATCTGTCTGTGGACCCCTATATGGAGATAGATGCCTTTGTG
CTCCTGCCAGCTCCAGCAAGGAGCCTGTCTCCGGCTCAGTGTGGCCACTAATGTGTCCGCCACCTTG
ACCTTCAATACCAGCAAGATCACTGGGTTCTTGAAGCCAGGAAAGTAAAAGTGGAACTGAAAGAAATCC
AAAGTTGGACTATTCAATGCAGAGCTGTTGGAAGCGCTCCTCAACTATTACATCCTTAACACCTTCTAC
CCCAAGTTCAATGATAAGTTGGCCGAAGCTTCCCCCTTCTCTGCTGAAGCGTGTTCAGCTCTACGAC
CTTGGGTGCAGATCCATAAGGACTTCTGTTCTTGGGTGCCAATGTCCAATACATGAGAGTTTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_004139

Insert Size: 1446 bp

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_004139.4
RefSeq Size:	1894 bp
RefSeq ORF:	1446 bp
Locus ID:	3929
UniProt ID:	P18428
Cytogenetics:	20q11.23
Domains:	BPI1
Protein Families:	Druggable Genome, Secreted Protein
Protein Pathways:	Toll-like receptor signaling pathway
MW:	53.4 kDa
Gene Summary:	The protein encoded by this gene is involved in the acute-phase immunologic response to gram-negative bacterial infections. Gram-negative bacteria contain a glycolipid, lipopolysaccharide (LPS), on their outer cell wall. Together with bactericidal permeability-increasing protein (BPI), the encoded protein binds LPS and interacts with the CD14 receptor, probably playing a role in regulating LPS-dependent monocyte responses. Studies in mice suggest that the encoded protein is necessary for the rapid acute-phase response to LPS but not for the clearance of LPS from circulation. This protein is part of a family of structurally and functionally related proteins, including BPI, plasma cholesteryl ester transfer protein (CETP), and phospholipid transfer protein (PLTP). [provided by RefSeq, Apr 2012]