

Product datasheet for **SC313675**

ZNF434 (ZSCAN32) (NM_017810) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF434 (ZSCAN32) (NM_017810) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF434
Synonyms:	HCCS-5; ZNF434
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC313675 representing NM_017810.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCGGTGACAACAGAGCTGTATCCCTCTGTCAGCAAGAATGGATGTGCCAGGCCCTGCACAAAGGGCC
CTCTACAGGGGTGCCACCCAGAGGAAGGACAGTCACGTCTCGCTGGCAACAGGTGTGCCCTGGGGCTAT
GAAGAGACCAAGACGCTCCTGGCTATTCTTAGTAGTTCTCAATTTTATGGAAAACCTCAGACCTGTCAG
CAGAACAGCCAGATCTACAGGGCCATGGCGGAAGGACTCTGGGAGCAGGGTTTTCTGCGGACCCAGAA
CAGTGTGCGACCAAGTTCAAAAGCCTACAGTTGAGTTACCGCAAAGTGAGGAGAGGCCGTGTGCCTGAG
CCTTGTATCTTTTATGAGGAAATGAATGCTCTTTCAGGCTCCTGGGCCTCTGCACCTCTATGGCAAGC
GATGCTGTTCTGGCCAAGAAGGAAGTGATATTGAGGCTGGAGAGCTGAATCACCAGAATGGGGAACCC
ACGGAGGTAGAAGATGGCACTGTGGATGGTGCAGACAGGGATGAAAAGGACTTCAGGAATCCTGGCCAG
GAAGTCAGGAACTAGACCTGCCAGTGCTGTTCCCAAACAGACTTGGTTTTGAGTTCAAGAACGAGATT
AAAAAAGAAAATCTAAAATGGGATGATTCAGAGGAAGTAGAAATAACAAGGCTTTACAGAGAAAGTCC
AGAGGAGTTTATTGGCACTCTGAGCTACAAAAGGCTTGGAGAGTGAGCCAACATCAAGAAGGCAATGT
AGAAATTTCTCCAGGGGAGAGTGAGGAGAAAACCCCATCCAGGAGAAGATGAGTCACCAGAGTTTTGT
GCCAGGGACAAAGCCTGTACACATATCCTCTGTGGGAAAACCTGCTCTCAGAGTGTGCACTCTCCCCAC
AAGCCAGCGCTCAAACCTGGAAAAAGTATCTCAATGTCTGAATGTGGGAAAACCTTTAGCCGAAGTTCT
TATCTTGTTCGGCATCAAAGAATCCACACAGGCGAGAAGCCTCACAAGTGCAGTGAGTGCGGGAAGGGC
TTAGTGAGCGCTCAAACCTCACTGCCACCTACGAACCTCACACAGGGGAGAGGCCCTATCAGTGTGGG
CAATGTGGGAAAAGCTTCAACCAGAGTTCCAGCCTCATTGTCCACCAGAGGACCCATACCGGGGAAAAG
CCTTACCAGTGCATTGTCTGTGGAAAGAGATTCAACAACAGTTCCAGTTCAGTGCTCACCAGGCGCATC
CACATGGGGAGAGGCCCATACAAGTGTGAGTGTGGGAAAATCTTCAACAATAGCTCCCACTTCAGT
GCCACCGAAAAACCCACACTGGTGAAGGCTTACAGGTGTTCTCACTGTGAGAGAGGCTTCACTAAG
AACTCTGCCCTCACCCGTCATCAGACAGTACACATGAAAGCAGTACTCTCATCAGGAAGGAAGAGAT
GCGTTATGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_017810

Insert Size: 1458 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_017810.3</u>
RefSeq Size:	3076 bp
RefSeq ORF:	1458 bp
Locus ID:	54925
UniProt ID:	<u>Q9NX65</u>
Cytogenetics:	16p13.3
Protein Families:	Transcription Factors
MW:	54.9 kDa
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]