

Product datasheet for **SC313697**

RSL1D1 (NM_015659) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	RSL1D1 (NM_015659) Human Untagged Clone
Tag:	Tag Free
Symbol:	RSL1D1
Synonyms:	CSIG; L12; PBK1; UTP30
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC313697 representing NM_015659.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGAGGATTCGGCCTCGGCCTCGCTGTCTTCTGCAGCCGCTACTGGAACCTCCACCTCGACTCCAGCG
GCCCCGACAGCACGGAAGCAGCTGGATAAAGAACAGGTTAGAAAGCAGTGGACGCTCTCTTGACGCAT
TGCAAGTCCAGGAAAAACAATTATGGGTTGCTTTTGAATGAGAATGAAAGTTATTTTTAATGGTGGTA
TTATGGAAAAATCCAAGTAAAGAACTGAGGGTCAGATTGACCTTGCCCTCATAGTATTCGATCAGATTCA
GAAGATATCTGTTTATTTACGAAGGATGAACCCAATTCAACTCCTGAAAAGACAGAACAGTTTTATAGA
AAGCTTTTAAACAAGCATGGAATTTAAACCGTTTCTCAGATTATCTCCCTCCAACTCTAAAGAAGGAA
TATAAATCCTATGAAGCCAAGCTCCGCCTTCTGAGCAGTTTTGATTTCTCCTTACTGATGCCAGAATT
AGGCGGCTCTTACCCTCACTCATTGGGAGACATTTCTATCAAAGAAAGAAAGTTCCAGTATCTGTAAC
CTTCTGTCCAAGAATTTATCAAGAGAGATCAATGACTGTATAGGTGGAACAGTCTTAAACATTTCTAAA
AGTGGTCTTGCAAGTGTATACGTATTGGTCACGTTGGAATGCAAATTGAGCACATCATTGAAAACATT
GTTGCTGTCAACAAAGGACTTTCAGAAAAATTGCCAGAGAAGTGGGAGAGCGTGAAACTCCTGTTTGTG
AAAAGTGAAGAAATCGGCTGCACTTCCCATCTTTTCTCGTTTGTGAGCAATTGGGATGAAGCCACAAA
AGATCTTTGCTTAAATAAGAAGAAAAAGAGGCAAGGAGAAAACGAAGAGAAAAGAAATTTTAAAAACAA
AAGGAGAGGAAGAAGAAGAGGCAGCAGGCTAGGAAGACTGCATCAGTTCTTAGTAAAGATGATGTGGCA
CCTGAAAGTGGTGATACTACAGTGAAGAAACCTGAATCAAAGAAGGAACAGACCCAGAGCATGGGAAG
AAAAACGTGGCAGAGGAAAAAGCCCAAGTTAAAGCAACAATGAATCCGAAGACGAAATCCCACAGCTG
GTACCAATAGGAAAGAAGACTCCAGCTAATGAAAAAGTAGAGATTCAAAAACATGCCACAGGAAAGAAG
TCTCCAGCAAAGAGTCTAATCCAGCACACCTCGTGGGAAGAAAAGAAAGGCTTTGCCAGCATCTGAG
ACCCCAAAGCTGCAGAGTCTGAGACCCAGGGAAAAGCCAGAGAAGAAGCCAAAAATCAAAGAAGAG
GCAGTGAAGGAAAAAGTCTTCCGCTGGGAAAAAGATGCGAGACAGACTCCAAAAAGCCAGAGGCC
AAGTTTTTACCCTCCTAGTAAATCTGTGAGAAAAGCTTCCACACCCCAAAAAATGGCCAAAAAA
CCCAAAGTACCCAGTCGACCTAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_015659

Insert Size: 1473 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_015659.2](#)

RefSeq Size: 5137 bp

RefSeq ORF: 1473 bp

Locus ID: 26156

UniProt ID: [O76021](#)

Cytogenetics: 16p13.13

MW: 55 kDa

Gene Summary: Regulates cellular senescence through inhibition of PTEN translation. Acts as a pro-apoptotic regulator in response to DNA damage.[UniProtKB/Swiss-Prot Function]