

Product datasheet for **SC313722**

A1BG (NM_130786) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	A1BG (NM_130786) Human Untagged Clone
Tag:	Tag Free
Symbol:	A1BG
Synonyms:	A1B; ABG; GAB; HYST2477
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


Fully Sequenced ORF: >SC313722 representing NM_130786.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGTCATGCTCGTGGTCTTTCTCTTGTGGGTGTACCTGGGGCCAGTGACAGAAGCAGCCATA
TTTTATGAGACGCAGCCAGCCCTGTGGCAGAGTCCGAATCACTGCTGAAACCCTTGGCCAATGTGACG
CTGACGTGCCAGGCCACCTGGAGACTCCAGACTTCCAGCTGTTCAAGAATGGGGTGGCCAGGAGCCT
GTGCACCTTGACTCACCTGCCATCAAGCACAGTTCCTGCTGACGGGTGACACCCAGGGCCGCTACCGC
TGCCGCTCGGGCTTGCCACAGGATGGACCCAGCTGAGCAAGCTCCTGGAGCTGACAGGGCCAAAGTCC
TTGCTGCTCCCTGGCTCTCGATGGCGCCAGTGTCTGGATCACCCCGGCCTGAAAACAACAGCAGTG
TGCCGAGGTGTGCTGCGGGTGTGACTTTTCTGCTGAGGCGGGAGGGCGACCATGAGTTTCTGGAGGTG
CCTGAGGCCAGGAGGATGTGGAGGCCACCTTCCAGTCCATCAGCCTGGCAACTACAGCTGCAGCTAC
CGGACCGATGGGAAGGCGCCCTCTCTGAGCCAGCGCTACTGTGACCATTGAGGAGCTCGCTGCACCA
CCACCGCTGTGCTGATGCACCATGGAGAGTCTCCAGGTCTGCACCCTGGCAACAAGGTGACCCCTC
ACCTGCGTGGCTCCCCTGAGTGGAGTGGACTTCCAGCTACGGCGCGGGGAGAAAGAGCTGCTGGTACCC
AGGAGCAGCACCCAGCCAGATCGCATCTTCTTTACCTGAACGCGGTGGCCCTGGGGGATGGAGGTAC
TACACCTGCCGCTACCGCTGCATGACAACAAAACGGCTGGTCCGGGGACAGCGCCGGTTCGAGCTG
ATTCTGAGCGATGAGACGCTGCCCGCGCCGAGTTCTCCCGGAGCCGGAGTCCGGCAGGGCCTTGGCG
CTGCGGTGCCTGGCGCCCTGGAGGGCGCGCTTCCGCTGGTGCAGGAGGACAGGGCGGGCGCCGC
GTGCACCGTTCCAGAGCCCGCTGGGACCGAGGCGCTTTCGAGCTGCACAACATTTCCGTGGCTGAC
TCCGCCAACTACAGTGCCTCTACGTGGACCTGAAGCCGCTTTCGGGGCTCCGCGCCAGCGAGCGC
TTGGAGCTGCACGTGGACGGACCCCTCCAGGCTCAGCTCCGGGCGACGTGGAGTGGGGCGGTCCTG
GCGGGCGGAGATGCCGTCTGCGCTGCGAGGGACCCATCCCCGAGCTCACCTTCGAGCTGCTGCGGAG
GGCAGAGCAAGGCCGTGAAGACGGTCCGACCCCGGGCCGCGGCGAACCTCGAGCTGATCTTCGTG
GGGCCCCAGCACCGGCAACTACAGGTCCGCTACCGCTCCTGGGTGCCCCACACCTTCGAATCGGAG
CTCAGCGACCCTGTGGAGCTCCTGGTGGCAGAAAGTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_130786

Insert Size:	1488 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_130786.3
RefSeq Size:	1766 bp
RefSeq ORF:	1488 bp
Locus ID:	1
UniProt ID:	P04217
Cytogenetics:	19q13.43
Domains:	ig, IGc2, IG
Protein Families:	Druggable Genome, ES Cell Differentiation/IPS, Secreted Protein
MW:	54.3 kDa
Gene Summary:	The protein encoded by this gene is a plasma glycoprotein of unknown function. The protein shows sequence similarity to the variable regions of some immunoglobulin supergene family member proteins. [provided by RefSeq, Jul 2008]