

Product datasheet for **SC313879**

ZNF572 (NM_152412) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF572 (NM_152412) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF572
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC313879 representing NM_152412.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGAGCAAGAAAAAACTGTTGGTCTCAGATTCTAACAGCTTTATGGAGAGGGAGAGTTTGAAAAGC
CCTTTACAGGAGATACAAGTATGAATAATTTGGAACTGTTACCACAATAATTCTAAGGCAGATAAAA
CTTAAAGAGAAACCTTCAGAATGGTCTAAAAGACATAGACCACAACATTATAAGCATGAGGATGCAAAA
GAAATGCCACTGACATGGGTTCAAGATGAGATTTGGTGTGATGTTCTATGAGAGTGATGGCAAGTCA
GAGAATTGGGAAATTTATAGCTAAAGAGGAGGAAAAACCCAATCACCAGGAATGGGACTCAGGAGAA
CATACCAATGCCTGTGTCCAGCAGAATTCATCCTTTGTAGACAGACCCTATAAATGTTCCGAATGTTGG
AAAAGCTTCAGTAATAGTTCTCATTGCGTACTCACCAGAGGCCACTCAGGAGAAAAGCCTTATAAAA
TGCTCTGAGTGTGCAAAATGTTTTGTAAACAGTCTCACCTGATTCAGCATCTAAGAATGCACACAGGA
GAGAAGCCCTACCAGTGGTGAATGTGGGAAAAGCTTCAGCAATACCTCCCATCTTATTATCCATGAG
AGAACTCACACGGGAGAGAAACCTACAAATGTCCGAGTGTGGGAAGAGATTCAGCAGCAGCTCTCAC
CTTATTACAGCATCACAGATCACATACAGGTGAAAAACCATATGAATGTTCTGTCTGCGGAAAAGGCTTC
AGTCACAGCTATGTCCTAATAGAACATCAGAGGACTCACACTGGAGAAAAACCTTATAAGTGCCCTGAT
TGTGGGAAGAGTTTTAGTCAGAGTTCAGCCTCATTGCCACCAGCGGACACACAGGAGAGAAGCCC
TACAAATGTCTTGAGTGTGAAAAAGCTTTGGTTGTAATCTACTCTAATAAAACATCAGAGAATACAT
ACAGGAGAAAAAGCCTTATCAATGTCCAGAATGTGGGAAGAATTTAGTCGTAGTTCAAACCTTATTACA
CACCAGAAAATGCACACAGGAGAGAAAATCCTATGAAAATTCTGAATATGAGGAAAGTTTGGGTGAGAAC
TGCAATGTGATAGAAGAATGCAGAATCCAGTTAGGAGAGAAAACCATATAGATGTTGTGAATGTGGGAAG
AGTTTTGGCCTTAGCTCCCATCTCATTAGACATCAGAGAACACATACAGGAGAAAAACCTTACAGATGT
TCTGAGTGTGGAAAACCTTCAGTCAGAGTTCCACCCTGGTGATTCACCAAAGGACACATACAGGAGAG
AAACCTTATAAATGTCCTGATTGTGGTGAAGCTTCAGTCAGAGCTTTAACCTTATCAGGCACCGGAGG
ACCCACATAGGGGAAAAACCTTACAAATGTACCAGCTGTGAGAAATGCTTCAGCAGAAGTGCCTACCTC
AGTCAGCATCGGAAAATTCAGTAGAAAAGCCTTTGAGTCTCCCGACGTTGGGATTTTCTCATGAA
TGGACTTGGAAAACCTGTTCAAGGGAAATGCCCTTCATCTCTTCATTTCCGCTCAAATTCATCTTCC
TGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

- Restriction Sites:** Sgfl-Mlul
- ACCN:** NM_152412
- Insert Size:** 1590 bp
- OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
- OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
- Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_152412.2](#)

RefSeq Size: 3276 bp

RefSeq ORF: 1590 bp

Locus ID: 137209

UniProt ID: [Q7Z3I7](#)

Cytogenetics: 8q24.13

Domains: zf-C2H2

MW: 61.2 kDa

Gene Summary: May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]