

Product datasheet for **SC313911**

Syntrophin (SNTB1) (NM_021021) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Syntrophin (SNTB1) (NM_021021) Human Untagged Clone
Tag:	Tag Free
Symbol:	SNTB1
Synonyms:	59-DAP; A1B; BSYN2; DAPA1B; SNT2; SNT2B1; TIP-43
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC313911 representing NM_021021.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCGGTAGCGGCGCGCGCGCGCGCTGGCCCGCTGGCGCGGAGGCGGCCGGCGCAGCGGAGC
GGGCTGCTGGAAGTTTTGGTGGCGGATCGCTGGCACAAAGTTCTGGTGAACCTGAGCGAGGACGCCCTG
GTTCTGAGCAGCGAGGAGGGCGCTGCGGGTACAACGGCATCGGGACCGCCACCAATGGCTGTTCTGC
AGGGGCGCCGGGCTGGGCACCCGGGCGCGGGCGCGCAGCCCCGGACTCGCCCGCCGGGGTCCGC
ACCGCTTTCACCGACTGCCCGAGCAGGTGCCCGAGTCCATCTCGAACGAGAAGCGTGGCGTGAAGGTG
CTGAAGCAGGAGCTGGCGGGCTGGGATCAGCATCAAGGGGGCAAGGAGAACAAGATGCCCATCTC
ATCAGCAAGATCTTCAAGGGCTGGCGCGGACCAGACCAAGCCCTGTACGTGGCGACGCCATCTG
TCCGTGAACGGAGCCGACTGCGGGACGCCACCCAGCAGGCGGTGCAGGCGTTGAAGCGCGCGGGC
AAGGAAGTGTCTGGAAGTGAAGTACATGCGAGAAGCCAGCCCTATGTGAAGAAAGGATCCCCAGTA
TCCGAGATTGGGTGGGAAACACCTCCGCTGAATCCCCTCGTTAGGGGGCAGCACCTCAGACCCCCG
TCATCGCAGTCTTCTCCTTCCACAGAGACCGAAAAGCATCCCCCTCAAATGTGTACGTCACTCGG
AGTATGGCCTTGGCCGACCTGAGAACAGGCAGCTTAAAATCCACTCTCCAGATGCTAAGCACACGGTG
ATCCTAAGGAGCAAGGACTCAGCCACGGCCAGGCATGGTTCACTGCCATCCATTCCAACGTTAATGAC
CTGCTGACCCGAGTGATTGCTGAGGTGACAGAGCAGCTGGGAAAAACAGGCATTGCTGGGAGCCGAGAG
ATTAGGCATCTTGGCTGGCTTGCAGAAAAGGTGCCAGGGGAGAGCAAGAAAACAGTGGAACACAGCCCTG
GTTGTGCTGACTGAGAAAAGCCTTTTAACTATGACAGCATGCCACGGAGGAAGGAAGCCTGGTTCAGC
CCAGTTCACACATACCCCTCTTCCACCAGGCTGGTCCATTCAAGTCCAGGAAAGGGATCACCCAG
GCTGGTGTGGATCTGCTTTTGAACGCGAATGGTACCAGGCAAGGGATTGAAACACATCTTTCAGA
GCAGAGACCAGCAGGGACCTCTCCCACTGGACAAGGAGCATAGTACAGGGTTGCCACAATTCTGCTGAA
CTCATTGCTGAAATCAGCACTGCTTGCACCTACAAAAACAGGAGTGCCGTTTGACCATACATTATGAG
AATGGATTTTCTATTACCACTGAACCACAGGAGGTGCCTTTCCCAAGACCATCATAAGCTCCTTAT
GAAAAGCTCAAAATGTCTTCAGATGATGGAATCAGGATGCTGATTTAGATTTTGGAGGCAAGATGGA
GAGATTCAACTGGACCTTATTCTGCCCAAGCCAATTGTTTTCATCATTCTCTCTGTCAGCT
AAGATTACAAGACTGGGCTTGGTGGCC TGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

Restriction Sites:

Sgfl-MluI

Plasmid Map:



ACCN:	NM_021021
Insert Size:	1617 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_021021.3</u>
RefSeq Size:	4975 bp
RefSeq ORF:	1617 bp
Locus ID:	6641
UniProt ID:	<u>Q13884</u>
Cytogenetics:	8q24.12
Domains:	PDZ, PH
MW:	58.1 kDa
Gene Summary:	Dystrophin is a large, rod-like cytoskeletal protein found at the inner surface of muscle fibers. Dystrophin is missing in Duchenne Muscular Dystrophy patients and is present in reduced amounts in Becker Muscular Dystrophy patients. The protein encoded by this gene is a peripheral membrane protein found associated with dystrophin and dystrophin-related proteins. This gene is a member of the syntrophin gene family, which contains at least two other structurally-related genes. [provided by RefSeq, Jul 2008]