

Product datasheet for **SC313939**

ZNF498 (ZSCAN25) (NM_145115) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF498 (ZSCAN25) (NM_145115) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZSCAN25
Synonyms:	ZNF498
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC313939 representing NM_145115.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCTTAAAGAGCATCCAGAGATGGCGGAAGCTCCTCAGCAGCAGTTGGGTATTCCTGTGGTGAAACTG
GAGAAAGAGTTGCCATGGGGCAGAGGAAGGGAGGACCCTAGTCCAGAGACTTTTCGGCTGAGGTTTCGG
CAGTTCGGCTACCAGGAGGCAGCTGGACCCAGGAAGCTCTTAGGGAGCTCCAGGAGCTCTGTCTGTCGG
TGGCTGAGGCCCGAGTTGCACACCAAGGAGCAGATCCTGGAGCTGCTGGTGTGGAGCAGTTCCTCACT
ATCCTGCCCGCGAGTTCTACGCTGGATCCGGGAGCATGGCCAGAGAGTGGCAAGGCCCTGGCCGCC
ATGGTGGAGGACCTGACAGAAAGAGCACTGGAGGCCAAGGCGGTTCCATGCCACAGGCAGGGAGAGCAG
GAGGAAACAGCACTTTGCAGAGGCGCTTGGGAGCCAGGCATCCAGCTGGGGCCAGTGGAGGTCAAGCCT
GAATGGGGATGCCCCCTGGGAAGGAGTTCAGGTCCAGACCCAGGTACCGAGGAGCAGCTCAGTCAG
GACCCTGGAGATGAGACACGGGCCTTCCAGGAGCAAGCACTACCTGTTCTGCAGGCGGGTCTGGCCTC
CCCGCAGTGAATCCAGAGACCAAGAGATGGCAGCTGGGTCTTTACTGCTGGATCGCAGGGGTTGGGG
CCATTTAAAGATATGGCCCTGGCCCTCCCTGAGGAGGAGTGGAGGCATGTGACCCAGCCAGATAGAC
TGCTTTGGGGAGTATGTGAACCGCAGGACTGCAGGGTCTCTCCAGGCGGTGGGAGCAAGGAAAAGGAG
GCAAAACCCACAGGAAGACCTGAAAGGGGCGCTGGTGGCACTGACATCAGAGAGGTTTGGGGAAGCC
TCTCTCCAGGGCCCTGGGCTCGGAAGGGTCTGTGAGCAGGAGCCTGGTGGCCCTGCAGGCAGTGCGCCT
GGGCTTCTCTCCCCAGCAGGTGCCATCCCCCTGCCTGACGAAGTCAAACCCACAGCTCCTTCTGG
AAGCCTTTCAGTGCCCTGAGTGTGGGAAAGGATTCAGTCGGAGCTCCAATCTCGTCAGGCACCAGCGA
ACCCACGAAGAGAAGTCTTATGGCTGTGTGGAGTGTGGGAAGGGCTTTACCCTGAGAGAATACCTGATG
AAGCACCAGAGAACCCACCTGGGAAAGAGGCCCTACGTGTGCAGCGAGTGTGGAAAACCTCAGCCAG
AGACACCACCTGGAGGTGCACCAGCGCAGCCACACTGGGGAGAAGCCCTACAAGTGGGGGACTGCTGG
AAGAGCTTCAGCCGAGCAGCACCTGCAGGTGCACCGGAGGACGCACACCGGGGAGAAGCCCTACACC
TGCGAGTGTGGCAAGAGCTTCAGCAGGAATGCCAATCTGGCGGTGCACCGGCGTGGCCACACTGGCGAG
AAGCCATATGGGTGCCAGGTGTGCGGGAAGCGGTTCAGCAAAGGGGAGCGGCTGGTCCGACACCAGAGA
ATCCATACAGGGGAGAAGCCCTACCACTGTCTGCCTGCGGGCGAAGCTTCAACCAGAGGTCCATCTC
AACCGGCACCAGAAGACCCAGCACCGCCAGGAGCCGCTGGTGCAGTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

Restriction Sites:

Sgfl-MluI

Plasmid Map:



ACCN:	NM_145115
Insert Size:	1635 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_145115.2</u>
RefSeq Size:	4359 bp
RefSeq ORF:	1635 bp
Locus ID:	221785
UniProt ID:	<u>Q6NSZ9</u>
Cytogenetics:	7q22.1
Domains:	zf-C2H2
Protein Families:	Transcription Factors
MW:	61.5 kDa
Gene Summary:	<p>This gene encodes a protein that bears some similarity to zinc finger proteins, which are involved in DNA binding and protein-protein interactions. Multiple alternatively spliced transcript variants have been identified, but the full-length nature for most of them has not been determined. [provided by RefSeq, Jul 2008]</p> <p>Transcript Variant: This variant (1), as well as variants 2 and 3, encodes the longest isoform (a).</p>