

Product datasheet for **SC313942**

MOXD1 (NM_015529) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MOXD1 (NM_015529) Human Untagged Clone
Tag:	Tag Free
Symbol:	MOXD1
Synonyms:	dJ248E1.1; MOX; PRO5780
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC313942 representing NM_015529.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGTGCTGCTGGCCGCTGCTCCTGCTGTGGGGCTGCTCCCGGGACGGCGGGGGGGCTCGGGCCGA
ACCTATCCGCACCGGACCCTCTGGACTCGGAGGCAAGTACTGGCTGGGCTGGAGCCAGCGGGGCAGC
CAGATCGCCTCCGCCTCCAGGTGCGCACTGCAGGCTACGTGGGCTTCGGCTTCTCGCCACCGGGGCC
ATGGCGTCCGCCGACATCGTCGTGGGCGGGGTGGCCACGGGCGGCCCTACCTCCAGGATTATTTTACA
AATGCAAAATAGAGAGTTGAAAAAGATGCTCAGCAAGATTACCATCTAGAATATGCCATGGAAAATAGC
ACACACACAATAATTGAATTTACCAGAGAGCTGCATACATGTGACATAAATGACAAGAGTATAACGGAT
AGCACTGTGAGAGTGATCTGGGCTACCACCATGAAGATGCAGGAGAAGCTGGTCCCAAGTACCATGAC
TCCAATAGGGGCACCAAGAGTTTGGGTTATTGAATCCTGAGAAAAGTGTGTCTATCTACAGCCTTA
CCATACTTTGATCTGGTAAATCAGGACGTCCTCCATCCAAACAAGATACAACATATTGGTGCCAAATG
TTTAAGATTCTGTGTTCCAAGAAAAGCATCATGTAATAAAGTTGAGCCAGTGATACAGAGAGGCCAT
GAGAGTCTGGTGACCACATCCTGCTCTATCAGTGCAGCAACAACCTTAAACGACAGCGTTCTGGAGTCC
GGCCACGAGTGCTATCACCCCAACATGCCCGATGCATTCTCACCTGTGAACTGTGATTTTTGCCTGG
GCTATTGGTGGAGAGGGCTTTTCTTATCCACCTCATGTTGGATTATCCCTTGGCACTCCATTAGATCCG
CATTATGTGCTCCTAGAAGTCCATTATGATAATCCCACTTATGAGGAAGGCTTAATAGATAAATCTGGA
CTGAGGTTATTTACACAATGGATATAAGGAAATATGATGCTGGGGTATTGAGGCTGGCCTCTGGGTG
AGCCTCTTCCATACCATCCCTCCAGGGATGCCTGAGTTCAGTCTGAGGGTCACTGCACTTTGGAGTGC
CTGGAAGAGGCTCTGGAAGCCGAAAAGCCAAGTGAATTCATGTGTTTGTGTTCTTCTCCATGCTCAC
CTGGCTGGCAGAGGCATCAGGCTGCGTCATTTTCGAAAAGGGAAGGAAATGAAATTAATCTGCCTATGAT
GATGATTTTGACTTCAATTTCCAGGAGTTTCAGTATCTAAAGGAAGAACAACAATCTTACCAGGAGAT
AACCTAATTAAGTGTGCTACAACACGAAAGATAGAGCTGAGATGACTTGGGAGGACTAAGCACC
AGGAGTAAAATGTGCTCTCATACCTTCTTTATTACCAAGAATTAATCTTACTCGATGTGCAAGTATT
CCAGACATTATGGAACAACCTCAGTTCATTGGGGTTAAGGAGATCTACAGACCAGTCACGACCTGGCCT
TTCATTATCAAAGTCCCAAGCAATATAAAAACCTTTCTTTTCATGGATGCTATGAATAAGTTTAAATGG
ACTAAAAGGAAGGTCTCCTTCAACAAGCTGGTCTCAGCCTGCCAGTGAATGTGAGATGTTCCAAG
ACAGACAATGCTGAGTGGTGGTCAAGGAATGACAGCATTACCTCCAGATATAGAAAGACCTATAAAA
GCAGAACCTTTGGTGTGGCAGCTCTTCTTCTTCCCTGCACAGAGATTTCTCCATCAACTTGCTT
GTTTGCCTTCTGCTACTCAGCTGCACGCTGAGCACCAAGAGCTTGTGA
ACGCGTACGCGGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTAAACGGCCGGC
  
```

- Restriction Sites:** SgfI-MluI
- ACCN:** NM_015529
- Insert Size:** 1842 bp
- OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
- OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
- Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_015529.3](#)

RefSeq Size: 3063 bp

RefSeq ORF: 1842 bp

Locus ID: 26002

UniProt ID: [Q6UVY6](#)

Cytogenetics: 6q23.2

Domains: Cu2_monoox_C

Protein Families: Druggable Genome

MW: 69.7 kDa