

Product datasheet for **SC313947**

Alpha Taxilin (TXLNA) (NM_175852) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Alpha Taxilin (TXLNA) (NM_175852) Human Untagged Clone
Tag:	Tag Free
Symbol:	TXLNA
Synonyms:	IL14; TXLN
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC313947 representing NM_175852.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAAGAACCAAGACAAAAAGAACGGGGCTGCCAAACAATCCAATCCAAAAAGCAGCCCAGGACAACCG
GAAGCAGGACCCGAGGGAGCCAGGAGCGGCCAGCCAGGCGGCTCCTGCAGTAGAAGCAGAAGGTCCC
GGCAGCAGCCAGGCTCCTCGGAAGCCGGAGGGGGCTCAAGCCAGAACGGCTCAGTCTGGGGCCCTTCGT
GATGTCTCTGAGGAGCTGAGCCGCAACTGGAAGACATACTGAGCACATACTGTGTGACAAATAACCAG
GGGGGCCCGGCGAGGATGGGGCACAGGGTGAGCCGGCTGAACCCGAAGATGCAGAGAAGTCCCGGACC
TATGTGGCAAGGAATGGGGAGCCTGAACCAACTCCAGTAGTCAATGGAGAGAAGGAACCCCTCAAGGGG
GATCCAAACACAGAAGAGATCCGGCAGAGTGACGAGGTGCGGAGCCGAGACCATCGAAGGCCACAGGAG
AAGAAAAAGCCAAGGGTTTGGGAAGGAGATCACGTTGCTGATGCAGACATTGAATACTCTGAGTACC
CCAGAGGAGAAGCTGGCTGCTCTGTGCAAGAAGTATGCTGAACTGCTGGAGGAGCACCGGAATTCACAG
AAGCAGATGAAGCTCCTACAGAAAAAGCAGAGCCAGCTGGTGAAGAGAAGGACCACCTGCGCGGTGAG
CACAGCAAGGCCGCTCCTGGCCCGCAGCAAGCTTGAGAGCCTATGCCGTGAGCTGCAGCGGCAACAACGC
TCCTCAAGGAAGAAGGTGTGCAGCGGGCCCGGGAGGAGGAGAAGCGCAAGGAGGTGACCTCGCAC
TTCCAGGTGACTGAATGACATTAGCTGCAGATGGAACAGCACAATGAGCGCAACTCCAAGTGCGC
CAAGAGAACATGGAGCTGGCTGAGAGGCTCAAGAAGCTGATTGAGCAGTATGAGCTGCGCGAGGAGCAT
ATCGACAAAGTCTTCAAACACAAGGACCTACAACAGCAGCTGGTGGATGCCAAGCTCCAGCAGGCCCAG
GAGATGCTAAAGGAGGCAGAAGAGCGGCACCAGCGGGAGAAGGATTTTCTCCTGAAAGAGGCAGTAGAG
TCCCAGAGGATGTGTGAGCTGATGAAGCAGCAAGAGACCCACCTGAAGCAACAGCTTGCCCTATACACA
GAGAAGTTTGAGGAGTTCCAGAACACACTTTCCAAAAGCAGCGAGGTATTACCACATTCAAGCAGGAG
ATGGAAAAGATGACTAAGAAGATCAAGAAGCTGGAGAAAAGAAACCACCATGTACCGGTCCCGGTGGGAG
AGCAGCAACAAGCCCTGCTTGAGATGGCTGAGGAGAAAACAGTCCGGGATAAAGAACTGGAGGGCTG
CAGGTA AAAATCCAACGGCTGGAGAAGCTGTGCCGGGCACTGCAGACAGAGCGCAATGACCTGAACAAG
AGGGTACAGGACCTGAGTGTGGTGGCCAGGGCTCCCTACTGACAGTGGCCCTGAGAGGAGGCCAGAG
GGGCTGGGGCTCAAGCACCCAGCTCCCCAGGGTACAGAAAGCGCCTTGCTACCCAGGAGCACCCGAGC
ACAGAAGCATCAGGCCAGACTGGCCTCAAGAGCCACCTCCGCCAGGGCTTAG
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

Restriction Sites:

Sgfl-MluI

Plasmid Map:



ACCN:	NM_175852
Insert Size:	1641 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_175852.3</u>
RefSeq Size:	4843 bp
RefSeq ORF:	1641 bp
Locus ID:	200081
UniProt ID:	<u>P40222</u>
Cytogenetics:	1p35.2
Protein Families:	Druggable Genome
MW:	61.9 kDa
Gene Summary:	May be involved in intracellular vesicle traffic and potentially in calcium-dependent exocytosis in neuroendocrine cells.[UniProtKB/Swiss-Prot Function]