

Product datasheet for **SC313948**

ARHGAP36 (NM_144967) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ARHGAP36 (NM_144967) Human Untagged Clone
Tag:	Tag Free
Symbol:	ARHGAP36
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

Fully Sequenced ORF:	>OriGene ORF sequence for NM_144967 edited <pre> ATGGGTGGCTGCATTCCTTTTCTGAAGGCAGCAAGGGCACTGTGCCCCAGAATCATGCC CCTTTGCTGTTGTTGTCCGCCTTCATTTTTTTAGTGAGTGTCTTGGGAGGAGCCCCAGGA CACAACCCCGACCGCAGGACGAAGATGGTATCGATACACAGCCTCTCTGAGCTGGAGCGT CTGAAGCTGCAAGAGACTGCTTACCACGAACCTCGTGGCCAGACATTTCTCTCCGAATTC AAACCTGACAGAGCTCTGCCTATTGACCGTCCGAACACCTTGGATAAGTGGTTTCTGATT TTGAGAGGACAGCAGAGGGCTGTATCACACAAGACATTTGGCATTAGCCTGGAAGAGGTC CTGGTGAACGAGTTTACCCGCCGCAAGCATCTTGAAGTACAGCCACGATGCAGGTTGAA GAAGCCACCGGTCAAGGCTGCGGGCCGTCGTGCGGGAAACGTGGTGCAAGGGTGTGGC CGCATCCGGCGCTTTTTCAGTCGAGGCGGAATGAGCCACCTTGCCCCGGGAGTTCACT CGCCGTGGGCGTCGAGGTGCAGTGTCTGTGGATAGTCTGGCTGAGCTGGAAGACGGAGCC CTGCTGCTGCAGACCCTGCAGCTTTCAAAAATTCCTTTCCAATTGGCCAACGACTTCTG GGATCCAAAAGGAAGATGAGTCTCAATCCGATTGCGAAACAAATCCCCAGGTTGTTGAG GCTTGCTGCCAATTCATTGAAAAACATGGCTTAAGCGCAGTGGGATTTTTACCCTTGAA TACTCCGTGCAGCGAGTGCCTCAGCTCCGTGAAGAATTTGATCAAGGTCTGGATGTAGTG CTGGATGACAATCAGAATGTGCATGATGTGGCTGCACTCCTCAAGGAGTTTTCCGTGAC ATGAAGGATTCTCTGCTGCCAGATGATCTGTACATGTCTTCTCCTGACAGCAACTTTA AAGCCCCAGGATCAGCTTTCTGCCCTGCAGTTGCTGGTCTACCTGATGCCACCCTGCCAC AGTGATACCCTGGAGCGTCTGCTGAAGGCCCTGCATAAAATCACTGAGAACTGCGAGGAC TCAATTGGCATTGATGGACAGTTGGTCCAGGCAACCGTATGACTTCCACTAACTTGGCC TTGGTGTGTTGGATCTGCTCTCCTGAAAAAGGAAAGTTTGGCAAGAGAGAGTCCAGGAAA ACAAAGCTGGGGATTGATCACTATGTTGCTTCTGTCAATGTGGTCCGTGCCATGATTGAT AACTGGGATGTCCTCTTCCAGGTGCCTCCCATATTCAGAGGCAGGTTGCTAAGCGCGTG TGGAAGTCCAGCCGGAAGCACTTGATTTTATCAGACGCAGGAACCTTGAGGAAGATCCAG AGTGACGCATAAAGATGGAAGAGGATGCACTACTTTCTGATCCAGTGGAAACCTCTGCT GAAGCCCGGGCTGCTGTCTTGCTCAAAGCAAGCCTTCTGATGAAGGTTCTCTGAGGAG CCAGCTGTGCCTTCCGGCACTGCCCGTTCCCATGACGATGAGGAAGGAGCGGGTAACCT CCCATTCCGGAGCAAGACCGCCATTGCTCCGTGTGCCCCGGGAGAAGGAGGCCAAACT GGCGTCAGCTACTTCTTCCTTAG </pre>
Restriction Sites:	Please inquire
ACCN:	NM_144967
Insert Size:	3000 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_144967.2, NP_659404.2</u>
RefSeq Size:	3082 bp
RefSeq ORF:	1644 bp
Locus ID:	158763
UniProt ID:	<u>Q6ZRI8</u>
Cytogenetics:	Xq26.1
Domains:	RhoGAP
Protein Families:	Transmembrane
Gene Summary:	GTPase activator for the Rho-type GTPases by converting them to an inactive GDP-bound state.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) represents the longest transcript and encodes the longest isoform (1).