

Product datasheet for **SC314022**

SLC26A10 (NM_133489) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SLC26A10 (NM_133489) Human Untagged Clone
Tag:	Tag Free
Symbol:	SLC26A10
Synonyms:	FLJ32854
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC314022 representing NM_133489.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGAGGCTTGATTTAGCATCCTTGATGTCAGTCCTAAGAGTCTGGGAAGTGCATTTAAGTCCTGGAGG
TTGGACAAGGCCCTCCACAGCACACCTTTCCATCCACTTCTATCCCAGGCATGGCTTTTGCTCTC
CTGGCCTCCGTGCCCGGTTGGACTCTACACTTCTTTCTTCCCGTCTCATCTACAGCTTGCTA
GGTACTGGGAGACACTGTCCACAGGAACCTTCGCCATACTCAGCCTCATGACAGGCTCGGCCGTGAG
CGGCTGGTGCCGGAACCCCTCGTGGGAATCTGAGCGGAATCGAGAAGGAGCAGCTGGACGCTCAACGG
GTTGGGGTAGCCGCGCCGTGGCCTTCGGGAGCGGGCGTTGATGCTGGGGATGTTGCTGCTGCAGCTC
GGCGTCTGTCCACCTTTTGTCCGAGCCTGTGGTCAAGGCGCTGACCAGCGGGCCGCGCTGCACGTG
CTTTGTCCCAGCTGCCAGCCTTTGGGGTGTCCCTCCGCGCCAGATCGGCTGCTTCTCTCTCTTC
AAGACGCTGGCCTCCTTGCTGACTGCGCTGCCTCGGAGCAGTCCGCGGAACTGACCATCTCCGCGCTC
AGCCTGGCGCTGCTCGTGCCGGTCAAGGAATTGAACGTGAGATCCGAGACCGGCTACCCACGCCGATC
CCGGGGGAAGTCGTCTTGGTGCTTCTGGCTCCGTGCTGCTTACCTTCTGTGGACACAAGATAC
CAAGTCCAGATAGTGGGGCTGTTGCCTGGAGGATTTCCCAACCCCTCTCCCAACCTGGCTGAGCTG
CCCAGGATTCTGGCTGACTCGTGCCATTGCACTGGTTAGTTTTGCGGTGTCTGCCCTCCCTGGCCTCC
ATCCATGCAGACAAGTATAGCTACACTATTGACTCCAACCAGGAGTTCCTGGCACATGGTGCCTCCAAC
CTCATCTCCTCCCTCTTCTTGTCTTCCCAACTCGGCTACGCTGGCCACCACCAATCTACTGGTGGAT
GCTGGTGGGAAAACACAGCTGGCAGGCCCTTCTCCTGCACAGTGGTCTGTGCGGTGCTGCTGGCTG
GGGCCCTCTTTTACTATCTGCCAAGGCTGTCTGGCTTGATCAACATCTCCAGCATGCGCCAGGTG
TTCTGCCAGATGCAGAACTCCACAACATATGGCACATCAGCCGAGTGGACTTTCTCCTCAGTCCCG
GGGCTCTGCATCCTGAGCTATCCAACACCACTGTACTTTGGGACCGTGGGCAGTTTCGCTGCAACCTG
GAGTGGCACCTGGGGCTCGGAGAAGGAGAAAAGGAGACTTCAAAGCCAGATGGCCCAATGTTGAGTT
GCTGAGCCTGTCAGGGTGGTGGTCTAGACTTCAGTGGTGTACCTTTGCAGATGCTGCTGGGGCCAGA
GAAGTGGTGCAGGTGAGGGAGAGGCTGGCCAGCCGATGTCGAGATGCTAGGATCCGCTCCTCTGGCT
CAGTGTAAATGCCTTGGTGCAGGGACACTGACCCGGTAGGACTCCTGGACAGGGTACTCCAGATCAG
CTGTTTGTGAGTGTGAGGATGCAGCTGCTTATGCCCTGGGGAGCCTGTTAAGGGGCAGTAGCACCAGG
AGCGGGAGCCAGGAGGCACTGGGCTGCGCAAGTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```

Restriction Sites: Sgfl-Mlul

Plasmid Map:



ACCN: NM_133489

Insert Size: 1692 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_133489.2](#)

RefSeq Size: 2446 bp

RefSeq ORF: 1692 bp

Locus ID: 65012

UniProt ID: [Q8NG04](#)

Cytogenetics: 12q13.3

Protein Families: Transmembrane
MW: 60.1 kDa
Gene Summary: Chloride/bicarbonate exchanger.[UniProtKB/Swiss-Prot Function]