

Product datasheet for **SC314077**

ZNF791 (NM_153358) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF791 (NM_153358) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF791
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC314077 representing NM_153358.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGATCGCC
ATGGACTCAGTGGCTTTTGAGGATGTGTCTGTGAGCTTCAGCCAGGAGGAGTGGGCTCTGCTGGCTCCT
TCACAGAAGAACTCTACAGAGATGTGATGCAGGAAACATTCAAGAACCCTGGCATCTATAGGGGAAAAA
TGGGAAGACCCGAATGTTGAAGATCAACACAAAAACCAAGGACGAAATCTAAGAAGCCATACGGGAGAG
AGACTCTGTGAAGGTAAAGAAGGTAGTCAATGTGCAGAAAACTTCAGTCCCAATCTCAGTGTGACGAAG
AAGACTGCCGGAGTAAAACCATATGAGTGTACTATCTGTGAAAAAGCCTTCATGCGTCTCTCATCCCTT
ACTAGACACATGAGGTCTCACACTGGATACGAGCTATTTGAGAAGCCATATAAATGTAAGGAGTGTGAG
AAAGCCTTTAGTTATCTCAAATCCTTTCAAAGACATGAAAGGAGTCACACTGGAGAAAAACCTATAAA
TGTAACAATGTGAAAAACCTTCATATATCACCAGCCCTTTCAAAGACATGAGCGGACTCACATTGGA
GAAAAACCTATGAATGTAAGCAATGTGAAAAAGCTCTTAGTTGTTCCAGTTCGCTTCGAGTTCATGAA
AGGATTCACACTGGAGAAAAACCTATGAATGTAACAATGTGGGAAAGCCTTCAGTTGTTCCAGTTCT
ATTCGAGTACACGAAAGAACTCACACTGGAGAGAAACCTATGCATGTAAGGAATGTGGGAAAGCCTTC
ATTTCCACACAAGTGTCTAACACACATGATAACACACAACGGAGATAGACCTTATAAATGCAAGAA
TGTGGAAGGCATTCAATTTTCCAGTTTTTACGAGTACATGAAAGAATTCACACTGGAGAGAAACCC
TATAAATGTAAACAATGTGGTAAAGCCTTCAGATGTTCCACCTCCATTCAAATTCATGAAAGAATTCAT
ACTGGAGAGAAGCCCTATAAATGTAAAGAATGTGGGAAATCTTTCAGTGCACGCCAGCCTTTCGAGTA
CACGTGAGAGTGCATACTGGAGAGAAACCTATAAGTGTAAAGAATGTGGGAAAGCCTTTAGTAGAATC
AGTTACTTTTGAATACATGAAAGGACTCACACTGGAGAGAAACCTACGAATGTAAAAAATGTGGGAAA
ACTTTCAATTATCCTCTAGATTTGAAAAATCCACAAGAGAAATCACACTGGAGAAAAACCTATGAGTGT
AAGGAATGTGCAAAAACCTTCATTTCTCTTGAGAACTTTGGAAGACACATGATCACCCACACTGGAGAC
GGACCTTATAAATGTAGGACTGTGGGAAGGTGTTCAATTTTCTAGTGCCTTACGAACACATGAAAGA
ACTCACACTGGAGAGAAACCTATGAATGTAACAATGTGAAAAAGCCTTTAGTTGTTCTAGTTACATT
CGGATACATAAAAGAACTCACACTGGGAGAAACCTTATGAATGTAAGGAATGCGGGAAGGCCTTTATT
TATCCACAAGCTTTCAAGGACACATGAGAATGCATACTGGAGAGAAACCTATAAATGTAAGAATGT
GGGAAGGCCTTTAGTCTTCACAGTTCCTTTCAAAGACATACAAGAATTCACAATTATGAGAAACCTCTT
GAATGTAAGCAATGTGAAAAAGCCTTCAGTGTGCCACATCCTTAAAAAACATATGAGAATGCACAAT
CGATAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_153358

Insert Size: 1731 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_153358.2</u>
RefSeq Size:	2495 bp
RefSeq ORF:	1731 bp
Locus ID:	163049
UniProt ID:	<u>Q3KP31</u>
Cytogenetics:	19p13.13
Domains:	KRAB, zf-C2H2
Protein Families:	Transcription Factors
MW:	66.9 kDa
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]