

Product datasheet for **SC314124**

NF2 (NM_181825) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	NF2 (NM_181825) Human Untagged Clone
Tag:	Tag Free
Symbol:	NF2
Synonyms:	ACN; BANF; merlin-1; SCH
Vector:	<u>pCMV6 series</u>



[View online »](#)

Fully Sequenced ORF: >NCBI ORF sequence for NM_181825, the custom clone sequence may differ by one or more nucleotides

```

ATGGCCGGGGCCATCGCTTCCCGCATGAGCTTCAGCTCTCTCAAGAGGAAGCAACCCAAG
ACGTTACCGTGAGGATCGTCACCATGGACGCCGAGATGGAGTTCAATTGCGAGATGAAG
TGGAAAGGAAGGACCTTTTGTATTTGGTGTGCCGACTCTGGGGCTCCGAGAAACCTGG
TTCTTTGGACTGCAGTACACAATCAAGGACACAGTGGCCTGGCTCAAAATGGACAAGAAG
GTACTGGATCATGATGTTTTCAAAGGAAGAACCAGTCACCTTTCACTTCTTGCCCAATTT
TATCCTGAGAATGCTGAAGAGGAGCTGGTTCAGGAGATCACACAACATTTATTCTTCTTA
CAGGTAAGAAGCAGATTTTAGATGAAAAGATCTACTGCCCTCCTGAGGCTTCTGTGCTC
CTGGCTTCTTACGCCGTCCAGGCCAAGTATGGTGACTACGACCCAGTGTTCACAAGCGG
GGATTTTTGGCCCAAGAGGAATTGCTTCCAAAAAGGGTAATAAATCTGTATCAGATGACT
CCGAAATGTGGGAGGAGAGAATTACTGCTTGGTACGCAGAGCACCGAGGCCGAGCCAGG
GATGAAGTGAATGGAATATCTGAAGATAGCTCAGGACCTGGAGATGTACGGTGTGAAC
TACTTTGCAATCCGGAATAAAAAGGGCACAGAGTCTGCTTGGAGTGGATGCCCTGGGG
CTTCACATTTATGACCCTGAGAACAGACTGACCCCAAGATCTCCTTCCCGTGAATGAA
ATCCGAAACATCTCGTACAGTGACAAGGAGTTTACTATTAACCCTGGATAAGAAAATT
GATGTCTTCAAGTTAACTCCTCAAAGCTTCGTGTTAATAAGCTGATTCTCCAGCTATGT
ATCGGGAACCATGATCTATTTATGAGGAGAAGGAAAGCCGATTCTTTGGAAGTTCAGCAG
ATGAAAGCCCAGGCCAGGGAGGAGAAGGCTAGAAAGCAGATGGAGCGGCAGCCCTCGCT
CGAGAGAAGCAGATGAGGGAGGAGGCTGAACGCACGAGGGATGAGTTGGAGAGGAGGCTG
CTGCAGATGAAAGAAGAAGCAACAATGGCCAACGAAGCACTGATGCGGTCTGAGGAGACA
GCTGACCTGTTGGCTGAAAAGGCCAGATCACCGAGGAGGAGGCAAACTTCTGGCCAG
AAGGCCGAGAGGCTGAGCAGGAAATGCAGCGCATCAAGGCCACAGCGATTTCGACGGAG
GAGGAGAAGCGCCTGATGGAGCAGAAGGTGCTGGAAGCCGAGGTGCTGGCACTGAAGATG
GCTGAGGAGTCAGAGAGGAGGGCCAAAGAGGCAGATCAGCTGAAGCAGGACCTGCAGGAA
GCACGCGAGGCGGAGCGAAGAGCCAAGCAGAAGCTCCTGGAGATTGCCACCAAGCCACG
TACCCGCCATGAACCAATTCCAGCACCGTTGCCTCCTGACATACCAAGCTTCAACCTC
ATTGGTGACAGCCTGTCTTTGACTTCAAAGATACTGACATGAAGCGGCTTTCCATGGAG
ATAGAGAAAGAAAAGTGAATACATGGAAAAGAGCAAGCATCTGCAGGAGCAGCTCAAT
GAACTCAAGACAGAAATCGAGGCCTTGAAGTGAAGAGAGGGAGACAGCTCTGGATATT
CTGCACAATGAGAACTCCGACAGGGGTGGCAGCAGCAAGCACAATACCATTAAAAAGCCT
CAAGCCCAAGGCAGAAGACCTATCTGCATT

```

- Restriction Sites:** Please inquire
- ACCN:** NM_181825
- OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
- OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
- Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_181825.1](#), [NP_861546.1](#)

RefSeq Size: 3106 bp

RefSeq ORF: 1773 bp

Locus ID: 4771

UniProt ID: [P35240](#)

Cytogenetics: 22q12.2

Protein Families: Druggable Genome

Gene Summary: This gene encodes a protein that is similar to some members of the ERM (ezrin, radixin, moesin) family of proteins that are thought to link cytoskeletal components with proteins in the cell membrane. This gene product has been shown to interact with cell-surface proteins, proteins involved in cytoskeletal dynamics and proteins involved in regulating ion transport. This gene is expressed at high levels during embryonic development; in adults, significant expression is found in Schwann cells, meningeal cells, lens and nerve. Mutations in this gene are associated with neurofibromatosis type II which is characterized by nervous system and skin tumors and ocular abnormalities. Two predominant isoforms and a number of minor isoforms are produced by alternatively spliced transcripts. [provided by RefSeq, Jul 2008] Transcript Variant: This variant (12) includes an alternate exon and differs in the 3' UTR, compared to variant 1. The resulting protein (isoform 2), one of the two predominant isoforms, has a distinct C-terminus, compared to isoform 1.