

Product datasheet for **SC314158**

Rotatin (RTTN) (AK096404) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Rotatin (RTTN) (AK096404) Human Untagged Clone
Tag:	Tag Free
Symbol:	Rotatin
Synonyms:	MSSP
Vector:	<u>pCMV6 series</u>



[View online »](#)

Fully Sequenced ORF:	>NCBI ORF sequence for AK096404, the custom clone sequence may differ by one or more nucleotides <pre> ATGGTCTGGCAGGGCTCATCAGGAACTCGGTTCATCAGCTGGCCGAGATCAGGGAGCGC GCTCTCAAGAGTATTCTCTGCAAGATTGAGCACAACTTAATCTGCTACGCTGATCTCATT CAGGAGAGGCAACTTTTTCTTCATTGCTGGAATGGTTCAATTTCCCGTCCGTTCCGATG AAGGAAGAGGTTCTGAACCTGTTAAGCAGATTGGTTAAGTATCCCCAGCAGTCCAACAT TTGGTTGACGTTGGTGCAGTAGAGTTCTTATCTAAGCTTCGGTCTAATGTGGAGCCAAAT CTGCAGGCTGAAATTGATGGCATTCTGGATGGACTTTTTCTTCTCCTTCGGAAGTTCCT GCACTATCTTCTGCCGATACCAAACCAATCAAACTGAATTGTCAAAAAACCCTGAAATC TTAACAGGATATTTCCCAAGACAAAAGTAATTTCCAGCAGATGGAAGTGCCGCCACGA CCAGTGGTAAATCAGACTGTGAAGTGCTGAAGTTTTCTACATTTCTTGGCTACCCCTG ACCACCACAGACAGACATGTCTCTCTCTAATGAAAGCTCTTAAAGAAGTAGTAACCAC ACTTTAATCTGGAACACCTGTGAAGTATTGAAGGATGTTATCATGCAAGATTTTCTGCT GAGATTTTCTTCAAAGGCCAAAAATTGTTGAGAGCCTTTTATCTCTGTTGAAACTGGCC TTTGGAGGAGATGAAAGCATCGCTGGCATTACAGTCGGTGTCTGCTGCAGCAGCTG TGCATGTATTTAAGAAACAGACTTAACCTTCACCGAGATCCAGGTTTTTCTCCAATAAA CACGACACAGTTTCCCAAATTCTTCTTGTCTTATTGTGATGAAGCAAGAGGTACTCAT CATTCCCAAGATCCTTCCCCAGGAAGCAGCAGCCCTCGGCCTTCTGTGGTTGGGCGCACA GGCCAGCGACCCAGAGGAGATGGCCAGGACTGGGATGCAGCGTCTCTAGTGGAAGTAGT AGTCATGCTCATGTAACTCCAGGATATCCGTTTCACTTACCTTTGGATATGGGACACATA GATCTGCCAGAGCTGGAACTGAAGACACATTGGAACACAATTCAGCAGCTCAGTCTT CCCCAGTTTTGTGCTCCATTCTGGAATCAGCTGTTCTCTCTTAAAGAACAGGTAGCAGA CAAGTGATAATAAGAGTTCTGGAATTGCTTACAGAGGATATGACACTTATTGGTGAAGCA ATTTCAACAGATATCTGGGATGACAGCAGCCTTTTTGGTATAGATATGAAGGAGAACTA CTCCTGGTGTGGGAGCCCTTGGAGAAACCATGTGCTACCATAAAAAGTAGTATCAGTTTG GAGCAGCCAGAGGTGATGCTTGTGCACCACAGAATGGCCTTTATCAGCATTTCCCTGTTT GCAGTTGACTGCTACAAACGCTTCTCCCTGTTGAAAAGGCAAGCGAGTTTTTATCAGAG CCTATGTCAACAGCATTATTTCTCTTTCTTGGACATGCCTATTTCTTGGAAATATCCA AATATTCATGAAGCTGTTGTGGCTATTTGGAACAGCTGAATTCTGAAACTACAGTATT TATAAACGAACTGCAGAGGCCGTTTATTCAATTGAATGCACCTGTAACCTCCTGTCTGAT ATTGGGAAAGAGGAGAGAAGAATCTCTTAGAATTAAGTGGAGCTGGCAGACCAAGCCTTG CGTAGCTTTTCTATCATCAGCATTTCCCGCTAATAAAGGAAATCATCAGCATTTGTTCA AAGATG </pre>
Restriction Sites:	Please inquire
ACCN:	AK096404
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>AK096404.1</u>
RefSeq Size:	2361 bp
RefSeq ORF:	2361 bp
Locus ID:	25914
Cytogenetics:	18q22.2
Gene Summary:	This gene encodes a large protein whose specific function is unknown. Absence of the orthologous protein in mouse results in embryonic lethality with deficient axial rotation, abnormal differentiation of the neural tube, and randomized looping of the heart tube during development. In human, mutations in this gene are associated with polymicrogyria with seizures. In human fibroblasts this protein localizes at the ciliary basal bodies. Given the intracellular localization of this protein and the phenotypic effects of mutations, this gene is suspected of playing a role in the maintenance of normal ciliary structure which in turn effects the developmental process of left-right organ specification, axial rotation, and perhaps notochord development. [provided by RefSeq, Jan 2013]