

Product datasheet for **SC314168**

AIFM3 (NM_144704) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	AIFM3 (NM_144704) Human Untagged Clone
Tag:	Tag Free
Symbol:	AIFM3
Synonyms:	AIFL
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC314168 representing NM_144704.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGGCGGCTGCTTCTCCAACCCAAACCAGTGGAGCTCAAGATCGAGGTGGTGTCTGCCTGAGAAGGAG
CGAGGCAAGGAGGAGCTGTCGGCCAGTGGGAAGGGCAGCCCCGGGCCTACCAGGCAATGGCACGGCC
CGCCACTTCCACACGGAGGAGCGCTGTCCACCCCTCACCCCTACCCAGCCCTCAGATTGCGTGGAG
GCTGCTGTCTGCCACGTCAAGGACCTCGAGAATGGCCAGATGCGGGAAGTGGAGCTGGGCTGGGGGAAG
GTGTTGCTGGTGAAGACAATGGGGAGTTCACACGCCCTGGGCCATAAGTGTCCGCACTACGGCGCACCC
CTGGTGAAGGCGTTCTGTCCCCTGGTGGTGGTGGTGGTGGTGGTGGTGGTGGTGGTGGTGGTGGTGGT
ACTGGGGACCTGGAGGACTTCCCTGGCCTGGACAGTCTACACAAGTTCAGGTGAAGATTGAGAAGGAG
AAGGTGTACGTCCGGGCCAGCAAGCAGGCCCTACAGCTGCAGCGAAGGACCAAGGTGATGGCCAAGTGT
ATCTCTCCAAGTCTGGGTACAGCAGTAGCACCAATGTGCTCATTGTGGTGCAGGTGCAGCTGGCCTG
GTGTGTGCAGAGACTGCGGCAGGAGGGCTTCTCCGACCGGATCGTCCGTGTGCAGCTAGACCGGCAC
TTCCCTACGACCGTCCCAAGCTCAGCAAGTCCCTGGACACACAGCCTGAGCAGCTGGCCCTGAGGCC
AAGGAGTTTTTCCGAGCCTATGGCATCGAGGTGCTCACCGAGGCTCAGGTGGTACAGTGGACGTGAGA
ACTAAGAAGTCTGTGTTCAAGGATGGCTTCAAGCTGGAGTACAGCAAGCTGCTGCTGGCACCAGGGAGC
AGCCCCAAGACTCTGAGCTGCAAAGGCAAAGAAGTGGAGAACGTGTTCACTATCCGGACGCCAGAGGAT
GCCAATCGCGTGGTGGGCTGGCCCGAGGCCGCAACGTGGTGGTGGTGGTGGTGGTGGTGGTGGTGGTGG
GAGGTGGCCGCTTACCTGACGGAGAAGGCCACTCTGTGTCTGTGGTGGAGCTGGAGGAGACGCCCTTC
AGGAGTTCTCTGGGGAGCGCGTGGTGGTGGTGGTGGTGGTGGTGGTGGTGGTGGTGGTGGTGGTGGTGG
TACATGCAGACGGAGGTGTCTGAGCTGCGGGGCCAGGAGGAAAGCTGAAGGAGTTGTGCTGAAGAGC
AGCAAGGTGCTGCGGGCTGACGTCTGCGTGGTGGCATTGGTGCAGTGCCCGCCACAGGCTTCTGAGG
CAAAGCGGCATCGGTTTGGATTCCCGAGGCTTCATCCCTGTCAACAAGATGATGCAGACCAATGTCCCA
GGCGTGTGTCAGCTGGCGATGCTGTACCTTCCCCCTTGCCTGGAGGAACAACCGCAAAGTGAACATT
CCACATTGGCAGATGGCTCATGCTCAGGGGCGCGTGGCAGCCAGAACATGTTGGCGCAGGAGGCGGAG
ATGAGCACTGTGCCCTACCTCTGGACCGCATGTTTGGCAAGAGCCTGCGCTACGCGGGCTACGGAGAA
GGCTTCGACGACGTCATCATCCAGGGGATCTGGAGGAGCTGAAGTTTGTGGCTTTTACACTAAAGGC
GACGAGGTGATCGCCGTGGCCAGCATGAACTACGATCCCATTGTGTCCAAGGTGCTGAGGTGCTGGCC
TCAGGCCGTGCCATCCGGAAGCGGGAGGTGGAGCTGTTTGTGCTGCACAGCAAGACTGGCGACATGTCC
TGGCTTACGGGAAAGGATCCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_144704

Insert Size: 1818 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_144704.2](#)

RefSeq Size: 2375 bp

RefSeq ORF: 1818 bp

Locus ID: 150209

UniProt ID: [Q96NN9](#)

Cytogenetics: 22q11.21

Domains:	pyr_redox, Rieske
Protein Families:	Druggable Genome
MW:	66.8 kDa
Gene Summary:	Induces apoptosis through a caspase dependent pathway. Reduces mitochondrial membrane potential.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) encodes the longest isoform (1).