

Product datasheet for **SC314175**

LINGO2 (NM_152570) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	LINGO2 (NM_152570) Human Untagged Clone
Tag:	Tag Free
Symbol:	LINGO2
Synonyms:	LERN3; LRRN6C
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC314175 representing NM_152570.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGCTTACACGGCCATATCATGCTGGCAGCCATTCTGGGTCTGGCTGTGGTGTAAATCTTCATGGGA
TCCACCATTTGGCTGCCCGCTCGCTGTGAGTGCTGCCAGAACAAATCTGTTAGCTGCACAGAAGG
CGATTGATCGCCATCCCAGAGGGCATTCCCATCGAAACCAAATCTTGGACCTCAGTAAAAACAGGCTA
AAAAGCGTCAACCCTGAAGAATTCATATCATATCCTCTGCTGGAAGAGATAGACTTGAGTGACAACATC
ATTGCCAATGTGGAACAGGAGCATTCAACAATCTCTTAACTGCGTTCCTCCGCCTAAAAGGCAAT
CGTCTAAAGCTGGTCCCTTTGGGAGTATTCACGGGGCTGTCCAATCTCACTAAGCTTGACATTAGTGAG
AATAAGATTGTCATTTACTAGACTACATGTTCCAAGATCTACATAACCTGAAGTCTCTAGAAGTGGGG
GACAATGATTTGGTTTATATACACACAGGGCATTGAGTGGGCTTCTAGCTTGGAGCAGCTCACCTG
GAGAAATGCAACTAACAGCAGTACCAACAGAAGCCCTCTCCACCTCCGAGCCTCATCAGCCTGCAT
CTGAAGCATCTCAATATCAACAATATGCCTGTGTATGCCTTTAAAAGATTGTTCCACCTGAAACACCTA
GAGATTGACTATTGGCCTTTACTGGATATGATGCCTGCCAATAGCCTCTACGGTCTCAACCTCACATCC
CTTTCAGTCACCAACCAATCTGTCTACTGTACCCTTCCTTGCCTTTAAACACCTGGTATACCTGACT
CACCTTAACTCTCCTACAATCCCATCAGCACTATTGAAGCAGGCATGTTCTCTGACCTGATCCGCCTT
CAGGAGCTTCATATAGTGGGGGCCAGCTTCGCACCATTGAGCCTCACTCCTTCCAAGGGCTCCGCCTT
CTACGCGTGCTCAATGTGTCTCAGAACCTGCTGGAACTTTGGAAGAGAATGTCTTCTCCTCCCTAGG
GCTCTGGAGGTCTTGAGCATTACAACAACCCTCTGGCCTGTGACTGCCGCTTCTCTGGATCTTGACAG
CGACAGCCACCCTGCAGTTTGGTGGCCAGCAACCTATGTGTGCTGGCCAGACACCATCCGTGAGAGG
TCTTTCGAAGATTCCATAGCACTGCCCTTCTTTTACTTTACCTGCAAAAAACCCAAAATCCGTGAA
AAGAAGTTGCAGCATCTGCTAGTAGATGAAGGGCAGACAGTCCAGCTAGAATGCAGTGCAGATGGAGAC
CCGACGCTGTGATTTCTGGGTGACACCCGAAGGCGTTTCATCACCACCAAGTCCAATGGAAGAGCC
ACCGTGTGGGTGATGGCACCTTGAAATCCGCTTTGCCAGGATCAAGACAGCGGGATGTATGTTTGC
ATCGCTAGCAATGCTGCTGGGAATGATACCTTCACAGCCTCCTTAACTGTGAAAGGATTGCTTCAGAT
CGTTTTCTTTATGCGAACAGGACCCCTATGTACATGACCGACTCCAATGACACCATTTCCAATGGCACC
AATGCCAATACTTTTCCCTGGACCTTAAACAATACTGGTGTCTACAGCTATGGGCTGCTTCACATTC
CTGGGAGTGGTTTTATTTGTTTTCTTCTCTTTTGTGTGGAGCCGAGGAAAGGCAAGCACAAAAAC
AGCATTGACCTTGAGTATGTGCCAGAAAAACAATGGTGTGTTGTGGAAGGGGAGGTAGCTGGACCC
AGGAGTTCAACATGAAAATGATTTGA
AGCGGACCGACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGGGATCTGGCAGCAAATGAT
ATCCTGGATTACAAGGATGACGACGATAAGGTTTAA
  
```

Restriction Sites: SgfI-RsrII

ACCN: NM_152570

Insert Size: 1821 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_152570.2](#)

RefSeq Size: 3236 bp

RefSeq ORF: 1821 bp

Locus ID: 158038

UniProt ID: [Q7L985](#)

Cytogenetics: 9p21.2-p21.1

Domains: LRRNT, LRRCT, LRR, ig

Protein Families: Druggable Genome, Transmembrane

MW: 68.1 kDa