

## Product datasheet for **SC314266**

### **SLC2A13 (NM\_052885) Human Untagged Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	SLC2A13 (NM_052885) Human Untagged Clone
Tag:	Tag Free
Symbol:	SLC2A13
Synonyms:	HMIT
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**Fully Sequenced ORF:** >SC314266 representing NM\_052885.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGTCCCGCAAGCAAGCGAGAATGTGGAGTACACGCTGCGGAGCCTGAGCAGCCTGATGGCGAGCGG
CGCAGGAAGCAGCCGGAGCCGGACGCGGGGAGCGCGGCCGGGGAGTGCAGCCTCCTGGCTGCCCGCGAA
TCGAGCACCGCCTGCAGAGCGCGGGCGCGGGCGCGGGCGTGGGGACCTGGAGCGCGCGCGCGG
CGGCAGTTCAGCAGGACGAGACCCCGCCTTCGTGTACGTGGTGGCCGTCTTCTCCGCGCTGGCGGGC
TTCCTGTTGGCTATGACACCGGGTGGTGTGAGGGCCATGCTGCTGCTCAAGCGGCAGCTCAGTCTG
GACGCGCTGTGGCAGGAGCTGCTGGTGTCCAGCACGGTGGGGCGGCTGCCGTCTCGGCGCTGGCCGGA
GGCGCCCTCAACGGCGTCTTCGGCCGCGCGCTGCCATCCTCTGGCCAGTGCCTCTTACCAGCCGGC
TCCGCGGTGCTGGCTGCGGCCAACAAACAGGAGACTGCTCGCCGGCCGCTGGTGTGGGACTCGGC
ATCGGCATTGCTTCTATGACAGTCCAGTGTACATTGCGGAGGTCTCACCACCAATTTAAGAGGCCGA
TTAGTACCATTAAATACCCTCTTCATCACAGGAGGGCAGTCTTTGCAAGTGTGTTGATGGAGCCTTC
AGTTATCTCCAGAAGGATGGATGGAGGTACATGTTGGGACTTGCAGCAGTTCGGCGGTTATACAGTTT
TTTGGCTTTCTTTTTGCCTGAAAGCCCTCGATGGCTTATTCAGAAAGGACAGACTCAGAAGGCCCGT
AGAATTTTATCTCAGATGCGTGGTAACCAGACCATTGATGAGGAATATGATAGCATCAAAAAACATT
GAAGAGGAGGAAAAAGAGGTTGGCTCAGCTGGACCTGTGATCTGCAGAATGCTGAGTTATCCCCAACT
CGCCGAGCTTTAATTGTGGGTTGTGGCCTACAAATGTTCCAGCAGCTCTCAGGCATTAACACCATCATG
TACTACAGTGAACCATCTGCAGATGTCTGGTGTGAAGATGATAGACTTGAATATGGCTGGCTTCA
GTTACAGCCTTCACAAATTTCAATTTTACACTTGTGGGAGTCTGGCTTGTGAGAAGTGGGCCGAGA
AAGCTTACCTTTGGTAGTTTAGCAGGTACCACCGTAGCACTATTATTCTTGCCTTGGGATTTGTGCTA
TCAGCCCAAGTTTCCCACGCATCACTTTTAAGCCAATAGCTCCGTCAGGTCAGAACGCCACTTGCACA
AGATACAGTACTGTAATGAATGTATGTTGGATCCAGACTGCGGTTTCTGCTACAAGATGAACAAATCA
ACTGTCATTGACTCCTCCTGTGTTCCAGTTAATAAAGCATCTACAAATGAGGCAGCCTGGGGCAGGTGT
GAAAATGAAACCAAGTTCAAAACAGAAGATATATTTGGGCTTACAATTTCTGCCCTACTCCATACTCC
TGGACTGCACTTCTGGCCTATTTTATATCTTGTCTTTTGCACCTGGAATGGGACCAATGCCTTGG
ACTGTGAATTCTGAAATATATCCCCTTTGGCAAGAAGTACAGGAAATGCATGTTTCTGGAATAAAC
TGGATTTTCAATGCTCCTGGTTTCACTAACATTTTACACACAGCAGAGTATCTTACATACTATGGAGCT
TTCTTCTCTATGCTGGATTTGCTGCTGTGGGACTCCTTTTTCATCTATGGCTGTCTTCTGAGACAAA
GGCAAAAATAGAGGAAATGAATCACTCTTTGACAACAGGCTATGTACATGTGGCACTTCAGATTCT
GATGAAGGGAGATATATTGAATATATTCGGGTAAGGGAAGTAACTATCATCTTTCTGACATGATGCT
TCTGATGTGGAATAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:** SgfI-MluI

**ACCN:** NM\_052885

**Insert Size:** 1947 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_052885.3</a>
<b>RefSeq Size:</b>	7018 bp
<b>RefSeq ORF:</b>	1947 bp
<b>Locus ID:</b>	114134
<b>UniProt ID:</b>	<a href="#">Q96QE2</a>
<b>Cytogenetics:</b>	12q12
<b>Domains:</b>	sugar_tr
<b>Protein Families:</b>	Transmembrane
<b>MW:</b>	70.4 kDa
<b>Gene Summary:</b>	H(+)-myo-inositol cotransporter (PubMed:11500374). Can also transport related stereoisomers (PubMed:11500374).[UniProtKB/Swiss-Prot Function]