

Product datasheet for **SC314289**

C5orf34 (AK056925) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	C5orf34 (AK056925) Human Untagged Clone
Tag:	Tag Free
Symbol:	C5orf34
Vector:	<u>pCMV6 series</u>



[View online »](#)

Fully Sequenced ORF:	>NCBI ORF sequence for AK056925, the custom clone sequence may differ by one or more nucleotides <pre> ATGGCAGCTGAACTGCGAATGATACTTTATGAAGATGATTCAGTACAAGTACAATATGTT GATGGTTCCACATTGCAACTTTCTCCCTGTGGCTCTGAATTTTATTTGAAAAGTCACCA CCTGTTTCAGCACATCCTTTAGAACAACCAGAAAGAAATTCGTCAAAGGACACATTTTGTC ATTAGCACTTACAGAGAGCAACTACAGCGAGCCCTAGATTTTCGAAACTCTTCAGCTACT TGCCCTTTTTTATCTGAAACCATCATACCTTCTGAAAGAAAAAGCATATCTTCATTGAC ATAACAGAAAGTGAGATGGCCAGTCTTGATACAGATGGTACCATGATATATATGGAGAGT GGCATTGTGAAGATAACATCTTTAGATGGTCATGCATACCTCTGCCTGCCAGATCTCAG CATGAATTTACAGTACATTTTTTGTGTAAGTTAGCCAGAAGTCAGACTCATCTGCAGTG TTGTCCAGAAACAAATAATAAGCCCCAAAAGATAAACTAGTTGAAAAAAGTGGCAAAATC TGTATACGTGGAAATTTACCAGGACAGAGACTGAAGAATAAGAAAAATGAGTTTCATTGC CAGATCATGAAATCCAAAGAACTTTAAGAAGATGAGTTGTGTAATGGAAGTGAAGGG AGGGAAGAGCTGCCTTCGCTGGTACAAAGCACACATGTGTATACACATGGGTCAAGCAG TGCTGGTCTGTGGCTGCCTGTCCAGAGGAATGGAATATCCTTTGTCTTTAGCACTTCAT TTTCATAACAAAATCAGCAATATGTCTAAAATTGATGCACATATAACTCAGAGTAGATTT TTAACTTCTGATATTTACAGAGGAAAGAGGAAAAGTGGTTTCTGTTCTTCCAGGGCCCTG TCACTCAGCTGTCCAGTCCCACACCTACACAGGTGGAATTTTTGTGATTCACTTTTACAG AGACAATCTGATGAATATTCCTATCCTGAACTAGTGAAATGGTTTGGTACAAAGGTGTT ACATATAGACTTACCCATCAAAATATGAAGTCAATAGAGATTTATTCTGGGGATGGATCT GTTTTCAAATCAGAAGGGGCTTATTTGGGAAGTATTTTACTTATTATTCTATTCAAGAA GGATCAGGAAAGAGAGAAGAGAAAACTTACTCAGTAAATAACCTTCCCTCCAGATAGACCG GGAAGTCCATTCACTGTCTGGTTCTTAATTAACAGGCAACAAGAAATTTCAACATTGT GTGAAGATGAGACTTTCTTTAAGCCATAACTATCGTATATGCTGCTGGAAGAAATGGTACCT GGGATAAATGATAGCAATATACTGCCTTTGGTTTTGAAAGAGTCACTCATACCCAGTGTG GGAAGATTTCTTGCTACTCTGATGACAAAGTACATGCTATCTTTTAGATGGCATTACT CTAACCCTAAATTTGGAATTTAGCTCTCCTATTGAGAAGAGACAGGTAAATCAAGGTCTT AACTTAGGTTGGTGAAGTTAACTTTTCTGATGGACAAGAGCAGTTAATTCAGATTGAA CACCTGAACCATATGAAAGATATGTGACAACAGTAACATCATGGTGCAGGAGACTGACC CAGACTAGTCCCAGAGAGATGCCACTCATTTCATCTTCTGTTCTCCAAGAAAATTGG TCTGTTGCTTCTGAAGTGAAGAGATACAGAAGTTAACTTGCTACTAGAAAATAGTGGT ATCCTAAACCAGATTTCTAACAAGAAAAATGAACAACAGTCTTTTGATCATTATAAACCA GGATCTTCAGAAACCTTGCTAGGAGAAAGTTAATGAAAATAGAGTATCAATAGCATTGAAA AAACCTCTGAAATCCTTCACGATATTGACTGTCTTCTATCAAACTCTAAAAAG </pre>
Restriction Sites:	Please inquire
ACCN:	AK056925
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [AK056925.1](#), [BAB71311.1](#)

RefSeq Size: 2441 bp

RefSeq ORF: 1917 bp

Locus ID: 375444

Cytogenetics: 5p12