

Product datasheet for **SC314301**

Uromuroid (UMOD) (NM_003361) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Uromuroid (UMOD) (NM_003361) Human Untagged Clone
Tag:	Tag Free
Symbol:	Uromuroid
Synonyms:	ADMCKD2; ADTKD1; FJHN; HNFJ; HNFJ1; MCKD2; THGP; THP
Vector:	<u>pCMV6 series</u>



[View online »](#)

Fully Sequenced ORF: >NCBI ORF sequence for NM_003361, the custom clone sequence may differ by one or more nucleotides

```

ATGGGGCAGCCATCTCTGACTTGGATGCTGATGGTGGTGGCCTCTTGGTTCATCACAACGCGAGCC
CTGACACCTCAGAAGCAAGATGGTGTCTGAATGTCACAGCAATGCCACCTGCACGGAGGATGAGGCCGT
TAGCAGCTGCACCTGTCAGGAGGGCTTACCGGCGATGGCCTGACCTGCGTGGACCTGGATGAGTGCGCC
ATTCTTGAGCTCACAACGCTCCGCCAACAGCAGCTGCGTAAACACGCCAGGCTCCTTCTCTGGTCT
GCCCCGAAGGCTTCCGCCTGTCGCCCGGTCTCGGCTGCACAGACGTGGATGAGTGCCTGAGCCTGGGT
TAGCCACTGCCACGCCCTGGCCACATGTGTCAATGTGGTGGCAGCTACTTGTGCGTATGCCCGCGGGC
TACCGGGGGATGGATGGCACTGTGAGTGTCTCCCGGGCTCTGCGGGCCGGGTTGGACTGCGTCCCCG
AGGGGACGCGCTCGTGTGCGGGATCCGTGCCAGGCGCACCGCACCTGGACGAGTACTGGCGCAGCAC
CGAGTACGGGGAGGGCTACGCTGCGACACGGACCTGCGCGGTGGTACCGCTTCGTGGCCAGGGCGGT
GCGCGCATGGCCGAGACCTGCGTCCAGTCTGCGCTGCAACACGGCCGCCCATGTGGCTCAATGGCA
CGCATCCGTCCAGCGACGAGGGCATCGTGAAGCCGCAAGGCTGCGCGCACTGGAGCGGCCACTGCTGCT
GTGGGATGCGTCCGTCCAGGTGAAGGCCTGTGCCGGCGGCTACTACGTCTACAACCTGACAGCGCCCCC
GAGTGTACCTGGCGTACTGCACAGACCCAGCTCCGTGGAGGGGACGTGTGAGGAGTGCAGTATAGACG
AGGACTGCAAATCGAATAATGGCAGATGGCACTGCCAGTGCAAACAGGACTTCAACATCACTGATATCTC
CCTCTGGAGCACAGGCTGGAATGTGGGGCCAATGACATGAAGGTGTCGCTGGGCAAGTGCCAGCTGAAG
AGTCTGGGCTTCGACAAGGTCTTGTGTACCTGAGTGACAGCCGGTGTCTGGGCTTCAATGACAGAGACA
ACCGGGACTGGTGTCTGTAGTGACCCAGCCGGGATGGCCCTGTGGGACAGTGTGACGAGGAATGA
AACCCATGCCACTTACAGCAACACCCTTACCTGGCAGATGAGATCATATCCGTGACCTCAACATCAAA
ATCAACTTTGCATGCTCCTACCCCTGGACATGAAAGTCAGCCTGAAGACCGCCCTACAGCCAATGGTCA
GTGCTCTAAACATCAGAGTGGGCGGACCGCATGTTACCGTGCAGGATGGCGCTTCCAGACCCCTTC
CTACACGACGCCCTACCAAGGCTCCTCCGTGACACTGTCCACTGAGGCTTTTCTTACGTGGGACCATG
TTGGATGGGGGCGACCTGTCCCATTGCACTGCTCATGACCAACTGCTATGCCACACCCAGTAGCAATG
CCACGGACCCCTGAAGTACTTCACTATCCAGGACAGATGCCACACACTAGAGACTCAACTATCCAAGT
GGTGGAGAATGGGAGTCTCCAGGGCGGATTTCCGTCCAGATGTTCCGGTTTGTGGAACTATGAC
CTAGTCTACCTGCACTGTGAAGTCTATCTGTGACACCATGAATGAAAAGTGAAGCCTACCTGCTCTG
GGACCAGATCCGAAGTGGGAGTGTATAGATCAATCCCGTGTCTGAACTGGGTCCCATCACACGGAA
AGGTGTCCAGGCCACAGTCTCAAGGGCTTTTAGCAGCTTGGGGCTCTGAAAGTCTGGCTGCCTCTGCT
CTCTCGGCCACCTTGACCCTGACTTTTCAGTGA
    
```

- Restriction Sites:** Please inquire
- ACCN:** NM_003361
- OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
- OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
- Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_003361.2 , NP_003352.2
RefSeq Size:	2327 bp
RefSeq ORF:	1923 bp
Locus ID:	7369
UniProt ID:	P07911
Cytogenetics:	16p12.3
Domains:	zona_pellucida, EGF_CA, EGF, EGF
Gene Summary:	<p>The protein encoded by this gene is the most abundant protein in mammalian urine under physiological conditions. Its excretion in urine follows proteolytic cleavage of the ectodomain of its glycosyl phosphatidylinositol-anchored counterpart that is situated on the luminal cell surface of the loop of Henle. This protein may act as a constitutive inhibitor of calcium crystallization in renal fluids. Excretion of this protein in urine may provide defense against urinary tract infections caused by uropathogenic bacteria. Defects in this gene are associated with the renal disorders medullary cystic kidney disease-2 (MCKD2), glomerulocystic kidney disease with hyperuricemia and isosthenuria (GCKDHI), and familial juvenile hyperuricemic nephropathy (FJHN). Alternative splicing of this gene results in multiple transcript variants. [provided by RefSeq, Jul 2013]</p> <p>Transcript Variant: This variant (1) encodes the shorter isoform (a). Both variants 1 and 2 encode the same isoform.</p>