

Product datasheet for **SC314317**

SF4 (SUGP1) (NM_172231) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SF4 (SUGP1) (NM_172231) Human Untagged Clone
Tag:	Tag Free
Symbol:	SF4
Synonyms:	F23858; RBP; SF4
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC314317 representing NM_172231.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGAGTCTCAAGATGGACAACCGGGATGTTGCAGGAAAGGCTAACCGGTGGTTTGGGGTTGCTCCCCCT
AAATCTGGAAAAATGAACATGAACATCCTTACCAGGAAGAGCTCATCGCTCAGAAGAAACGGAAATT
GAAGCCAAAAATGGAACAGAAAGCCAAGCAGAATCAGGTGGCCAGCCCTCAGCCCCACATCCTGGCGAA
ATCACAAATGCACACAACCTTCTCTGCATTTCCAACAAGTTTGCCAACGATGGTAGCTTCTTGCAGCAG
TTTCTGAAGTTGCAGAAGGCACAGACCAGCACAGACGCCCGACCAGTGCGCCACGCGCCCTCCAGC
ACACCCACCCCGAGCGTGGGAAGAGGTCCCTGCTCATCAGCAGGCGGACAGGCTGGGGCTGGCCAGC
CTGCCGGGCCCTGTGAAGAGCTACTCCCACGCCAAGCAGCTGCCCGTGGCGCACCCCGAGTGTCTTC
CAGTCCCCTGACGAGGACGAGGAGGAGGACTATGAGCAGTGGCTGGAGATCAAAGTTTACCCCCAGAG
GGAGCCGAGACTCGAAAGTGATAGAGAAATTGGCCCGCTTGTGGCAGAAGGAGGCCCGAGTTAGAA
AAAGTAGCTATGGAGGACTACAAGGATAACCCAGCATTTGCATTTTGCACGATAAGAATAGCAGGGAA
TTCTCTACTACAGGAAGAAGGTGGCTGAGATAAGAAAGGAAGCACAGAAGTCCGAGGCAGCCTCTCAG
AAAGTTTACCCCGAGAGGACGAAGAGGTCAAGAACCTTGCGAAAAAGTTGGCCAGGTTTATAGCGGAC
GGGGTCCCGAGGTGAAACATTGCCCTCCAGAAACCGTGAGAACCAGGCATTAGCTTTCTGTAT
GAGCCCAATAGCCAAGGGTACAAGTACTACCGACAGAAGCTGGAGGAGTCCGGAAAGCCAAGGCCAGC
TCCACAGGCAGCTTACAGCACCTGATCCCGCCCTGAAGCGCAAGTCCCCTCCTGAGGCCCTGTCAAGG
TCCTTACCCCGAGCCACCTGCCCGCCCTCGTCCACGCCTGCGCCACTATCATCCCTGCTCCAGCT
GCCCCCGGAAGCCAGCCTCCGCAGCCACCGTGAAGAGGAAGCGGAAGAGCCGGTGGGGCCTGAAGAG
GATAAGGTAGAGCTCCCACCTGCTGAACCTGGTGCAGAGGGACGTGGATGCCTCTCCCTCGCCTCTGTCA
GTTCAGGACCTCAAGGGGCTCGGCTATGAGAAGGGGAAGCCTGTGGGTCTAGTGGGCGTCACAGAGCTT
TCAGACGCCCAGAAGAAGCAGCTGAAGGAGCAGCAGGAGATGCAGCAGATGTACGACATGATCATGCAG
CACAAGCGGGCCATGCAGGACATGCAGCTGCTGTGGGAGAAGGCAGTCCAACAGCACCAGCACGGCTAT
GACAGTGTAGGAGGTGGACAGCGAGCTGGGCACCTGGGAGCACCAGCTGCGGCGCATGGAGATGGAT
AAGACCAGGGAATGGGCCGAGCAGCTGACAAAGATGGGCCGGGGCAAGCACTTATCGGAGACTTCTG
CCTCCAGACGAGCTGGAAAAGTTTATGGAGACCTTCAAGGCCCTGAAGGAGGGCCGTGAGCCTGACTAC
TCAGAGTACAAGGAGTTCAAGCTGACTGTGGAGAACATCGGCTACCAGATGCTGATGAAGATGGGCTGG
AAGGAGGGCGAGGGGCTGGGCTCAGAGGGCCAGGGCATCAAGAACCAGTGAACAAGGGCACCACCACA
GTGGACGGCGCTGGCTTCGGCATTGACCGGCCGCGGAGCTCTCCAAGGAGGACGACGAGTATGAGGCG
TTCCGCAAGAGGATGATGCTGGCCTACCGCTTCCGGCCCAACCCCTGAACAATCCAGACGGCCTTAC
TACTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_172231

Insert Size: 1938 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_172231.3</u>
RefSeq Size:	2116 bp
RefSeq ORF:	1938 bp
Locus ID:	57794
UniProt ID:	<u>Q8IWZ8</u>
Cytogenetics:	19p13.11
MW:	72.5 kDa
Gene Summary:	SF4 is a member of the SURP family of splicing factors.[supplied by OMIM, Sep 2003]