

Product datasheet for **SC314350**

JPH1 (NM_020647) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	JPH1 (NM_020647) Human Untagged Clone
Tag:	Tag Free
Symbol:	JPH1
Synonyms:	CMT2K; JP-1; JP1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC314350 representing NM_020647.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGACGGGCGGAAGGTTTCGACTTCGACGATGGCGGCACCTACTGCGGCGGCTGGGAGGAGGCAAGGCG
CACGGGCATGGCATCTGCACGGGGCCCAAGGGCCAGGGCGAGTACTCGGGCTCCTGGTCGCACGGCTTC
GAGGTGGTTCGGAGGCTACACCTGGCCCAGCGGCAACACCTACCAGGGCTACTGGGCGCAGGCAAGCGG
CACGGGCTGGGGGTGGAGACGAAGGGCAAGTGGATGTACCGGGGGAGTGGTCACATGGTTTCAAGGGG
CGCTACGGGGTCCGGCAGAGCCTGTGCACCCCCGCTCGCTACGAGGGTACCTGGAGTAACGGGCTGCAA
GACGGGTACGGCGTGGAGACCTACGGGGACGGAGGTACCTACCAGGGCCAGTGGGCGGAGGCATGCGG
CATGGCTACGGCGTGGCCAGAGCGTGCCTACGGCATGGCCACGGTGATCCGCTACCGCTGCGTACC
TCGCTGGCCTCGTGGCAGCGAGCAGCAATGGCAGCGTGCTCCACGACGCCGAGCCGCCGCCGAC
AGCCCGGCCCGCACCCCGCGCGTTTCGTGCTCAACTTCCACGCAGACGCTGAGCTAGCGGGCAAGAAG
AAGGGCGGCTCTTCCGAGGGGCTCCCTTCTTGAAGCATGAACTTCGCAAGTCCGAATCCAAGTCT
TCCATCTCGAGCAAGCGAGCTCTGTCCGACGCGACGCGGCCATGAGCAGAATTAGTTCCAGCGATGCC
AACTCCACGATCAGCTTTGGCGATGTAGATTGTGATTTTGTCCCGGTGGAAGACCACGTGGACGCCACC
ACCACGGAACCTACATGGGCGAGTGAAGAACGACAAGCGCAACGGCTTCGGCGTTAGCGAGCGCTCC
AATGGCATGAAGTATGAAGGGGAGTGGGCAATAACAAGAGGCATGGATATGGTGTACCGTGTTTCTT
GACGGCTCCAAAGAAGAGGGAATAACAAAAATAATTTCTGGTCCGTGGGATAAGGAAGCAGCTTATA
CCAATAAGACATACAAAACTAGGGAGAAGGTGGACAGAGCAATTGAAGGCGCCAAAGGGCAGCTGCC
ATGGCCAGAACCAAAGTGGAAATAGCAAATCAAGGACTGCACATGCCAGAGCGAAGGCCGATGCCGCC
GACCAGGCCGCGCTGGCCGCTCGCCAGGAGTGCACATCGCGAGAGCTGTGGCCAGGGAGCTGTACCT
GATTTCTACCAACCAGGCCCTGATTACGTCAAACAGAGATTTTCAGGAAGGTGTAGATGCTAAAGAAAT
CCAGAAGAAAAGGTACCAGAAAAGCCACCTACACCAAAGGAGTCTCCTCATTTTTATCGAAAGGCACG
ACACCCCCAAGATCTCCTGAGGCAAGTCCAAACACAGCCACTCTCCTGTTCTCCCCAAAGCCCTG
AAGAAGCAAAACCCAGCTCAGGGGCGAGACTCAACCAAGACAAAAGGAGTGTGGCTGATGAGCAGGTG
ACGGCCATTGTCAATAAGCCCTTGATGTCAAAGGCTCCACGAAGGAGGAGGAGCGTTGTGCCCCAG
TCCAAGTACTCTGGCCGCCACCACATCCCCAACCCAGTAACGGGGAGCTGCATTCTCAGTATCACGGC
TACTACGTGAAGCTGAACGCCCCCAGCACCTCCAGTAGACGTGGAGGACGGCGATGGATCCAGCCAG
TCTTCTCAGCACTGGTGCACAAGCCATCCGCTAACAAGTGGAGTCCCTCAAATCTGTGACAAAACCA
GTTGCCAAAGAAAGCAAGCTGAGCCAAAGCTAAGAAGTCTGAACTTGTATACCAAGAATCCAGCA
AGCAACGATTATGCCCTGCTTTGGAAAAAGAAGCCAATTCAGGCCCTAATTCAATCATGATTGTCTT
GTCATGCTGTTGAATATCGGGTTGGCCATTCTTTTGTCACTTTCTAACTTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_020647

Insert Size: 1986 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_020647.3</u>
RefSeq Size:	4427 bp
RefSeq ORF:	1986 bp
Locus ID:	56704
UniProt ID:	<u>Q9HDC5</u>
Cytogenetics:	8q21.11
Domains:	MORN
Protein Families:	Transmembrane
MW:	71.7 kDa
Gene Summary:	Junctional complexes between the plasma membrane and endoplasmic/sarcoplasmic reticulum are a common feature of all excitable cell types and mediate cross talk between cell surface and intracellular ion channels. The protein encoded by this gene is a component of junctional complexes and is composed of a C-terminal hydrophobic segment spanning the endoplasmic/sarcoplasmic reticulum membrane and a remaining cytoplasmic domain that shows specific affinity for the plasma membrane. This gene is a member of the junctophilin gene family. [provided by RefSeq, Jul 2008] Transcript Variant: This variant (1) represents the longer transcript. Both variants 1 and 2 encode the same protein. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.