

Product datasheet for **SC314356**

ZNF800 (NM_176814) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF800 (NM_176814) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF800
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC314356 representing NM_176814.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCCTTTAAGGGACAAATACTGTGAGACTGACCACCATCATCACGGATGCTGTGAACCAGTTTATATC
CTGGAACCTGGAGATCCTCCTTTGTTACAGCAACCACTACAGACATCCAAATCTGGTATTCAACAAATA
ATTGAGTGCTTTTCGATCAGGAACCTAAACAACCTTAAGCATATTTTATTAAGATGTGGACACTATTTTT
GAATGTAAGTTATGCCGAGTCTCTTCAGAGGATTACCAATTTAATTACCCATAAAAAATTCTACTGC
CCACCAAGTCTCCAGATGGATGACAACCTTCTGATGTAATGATAAAACAAAGCCAAAGCCATAAATGAT
CTCCTAGAAGCCATATATCCAAGTGTGGACAAACGAGAATATATTATTAAGCTAGAACCATAGAAACT
AATCAAAATGCAGTATTTCAATATATTTTCGAGGACTGATAATCCTATTGAAGTCACAGAGTCAAGCAGT
ACTCCTGAACAAACCGAAGTTCAGATACAGGAACTAGCACTGAACAGTCAAAAACAGTACCGGTTACA
GATACAGAGGTGGAACCTGTAGAGCCCCCTCCTGTTGAGATTGTTACAGATGAAGTTCACCTACATCT
GATGAACAACTCAGGAGTCGAGGCTGACTTGGAACTTCTGACAATTCTGATTTTGGTACCAGTTG
ATATGTTGCTTTGTAGAAAAGAATTCAATTCTAGACGAGGTGTTCCGCTCACATTCGAAAAGTACAC
AAGAAAAAGATGGAAGAACTAAAAAGTACATTGAAACACGAAAGAAATCCAAACCAATCCTCTAAAGGA
CGCAGTAAGAATGTTCTAGTTCATTAAGTAGGAGTTGTCAGATGTTGTAATCATTTGCTACAAAA
GCGAATGTAAGGAGGCATTTTGTGAAGTTCATAGAGGACTAAGGAGGGATTCAATTACTCCTGATATA
GCAACAAAGCCTGGGCAACCTTTGTTCCCTGGATTCTATTTCTCTAAAAAATCTTTAAGACTCGAAAA
CAAAAGTCTTCTCAAAGGCTGAATACAATTTAACTGCATGCAATGCCTCCTTTGCAAGAGGAAATAT
AGTTCACAAATAATGCTTAAAAGACATATGCAAAATGTCCACAAGATAACTCTTTCTGGAACAACTCT
AAAAGAGAAAAAGGCCCTAATAACTGCCAACAGTTCAGAAATAAAAGTTAAAGTTGAACCAGCAGAT
TCTGTAGAATCTTCACCCCTTCCATTACCCATTCTCCACAGAATGAATTAAGGGAACAAATCATTCA
AATGAAAAAAGAACACACCGGCAGCACAGAAAAATAAAGTTAAACAAGACTCTGAAAGCCCTAAATCA
ACTAGTCCGTCGGCTGCAGGTGGCCAGCAAAAAACAGAAAAACAAAACCTTTCAGCTGGCTTTGACTTT
AAGCAACTTTACTGTAACTTTGTAAACGTCAGTTTACTTCCAAACAGAAGTACTGACTAAACACATCGAG
TTGCACACAGATGGAAATAACATTTATGTTAAATTTCTACAAGTGCCTCTTTGCACTTATGAACTCGT
CGGAAACGTGATGTGATACGACATATAACTGTGGTTCATAAAAAGTCATCTCGTTATCTTGGGAAAAATA
ACAGCCAGTTTAGAGATCAGAGCTATAAAAAAGCCTATTGATTTTGTCTAAATAAAGTGGCAAAAAAGA
GGCCCTTCGAGGGATGAAGCAAAACATAGTGATTCAAAACATGATGGCACTTCTAACTCTCCTAGTAAA
AAGTATGAAGTAGCTGACGTCGGTATTGAAGTAAAGTCACAAAAAATTTCTCTTACAGATGCAAT
AAATGTGGAAGGCATTTGCCAAAAAGACTTACCTTGAACATCATAAGAAAACTCATAAGGCAAAATGCT
TCCAATTCACCTGAAGGAAACAAAACCAAGGCCGAAGTACAAGATCTAAGGCTCTTGTCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_176814

Insert Size: 1995 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_176814.4](#)

RefSeq Size: 4489 bp

RefSeq ORF: 1995 bp

Locus ID: 168850

UniProt ID: [Q2TB10](#)

Cytogenetics: 7q31.33

MW: 75.2 kDa

Gene Summary: May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]