

## Product datasheet for **SC314361**

### ZNF333 (NM\_032433) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF333 (NM_032433) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF333
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**Fully Sequenced ORF:** >SC314361 representing NM\_032433.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGGAATCCGTCACCTTTGAGGATGTGGCCGTGGAGTTTCATCCAGGAGTGGGCATTGCTGGACAGCGCA
CGGAGGAGCCTGTGCAAATACAGGATGCTTGACCAGTGCAGGACCCTGGCCTCCAGGGAACTCCACCA
TGCAAACCCAGTTGTGTCTCCAGCTGGGGCAAAGAGCAGAGCCAAAGGCAACAGAACGAGGGATTCTC
CGTGCCACAGGTGTTGCCTGGGAATCTCAACTAAACCCGAAGAGTTGCCTTCTATGCAGGATCTTTTG
GAAGAAGCATCTCCAGGGACATGCAAATGGGGCCGGGGCTGTTCTGAGGATGCAGCTGGTGCCCTCC
ATAGAAGAGAGGGAGACACCATTGACTCGAGAGGACCGGCCAGCTCTCCAGGAGCCGCTTGGTCTCTG
GGATGCACGGGACTGAAGCCGCTATGCAGATTACAGGGTGGTGATACCAGTGCCTACTCTGGGCCAC
CGCAACCCATGGGTGGCCAGGGATTCTGCTGTGCCTGCACGTGACCCTGCCTGGCTTCAGGAGGACAAA
GTGGAGGAAGAAGCTATGGCTCCTGGGCTGCCAACCGCCTGTTACAGGAACCAGTACCTTTGCAGAT
GTGGCTGTGGTGTTCACCCAGAAGAATGGGTGTTTCTGGACTCTACTCAGAGGAGCCTGTATAGAGAT
GTGATGCTGGAGAACTACAGAACTGGCCTCTGTGGCTGATCAACTGTGCAAACCAATGCGTTGTCT
TATTTGGAAGAAAAGAGGAGAGCAGTGGACCACTGACAGGGCGTCTCTCAGACACCTGTGCAGAACCT
CAGTGTCAACCCCAAGAGGCAATTCCTAGCCAAGATACTTTTACAGAGATCCTGTCCATTGATGTGAAA
GGGAGAACCTCAGCCTGGAGAAAACTCTATAAATAAATGAACCTTGAGAAACCTTTTAAACAGCATT
GAACCACTTTTCCAGTACCAGAGAATTCATGCTGGAGAGGCATCCTGTGAATGTCAAGAGATTAGAAAT
TCCTTCTTCCAGAGTGCCACCTAATTGTGCCGAGAAAAATCCGTAGTGGGGATAAATCCTATGCATGT
AACAAATGTGAAAAATCCTTCAGATACAGCTCTGACCTTATCAGGCATGAGAAGACTCATACTGCAGAG
AAGTGCTTTGACTGTCAAGAATGTGGCAAGCCTTCAAATATTCCTCGAATCTCCGGCGACACATGAGA
ACCCATACCGGAGAGAAGCCATTTGAATGTAGTCAGTGTGGGAAAACCTTACAGGAACCTTTAACCTG
ATTTTGCACCAGAGAAACACACAGGAGAGAAGCCCTACGAGTGTAAAGATTGTGGGAAAGCCTTCAAT
CAGCCATCATCCCTCAGGAGCCACGTGAGAACTCACACTGGAGAGAAGCCCTTTGAATGCAGCCAGTGT
GGGAAAGCCTTACGGAACTCTTCACTGAAGACACATCTGCGAACCCATACCAGAGAGAAACCATAT
GAATGCAACCAGTGTGGCAAGCCTTCCGGACGAGCACTCATCTGAACGTGCACAAGAGGATACACACA
GGGAGAAACTGTATGAGTGCAGACTTGGGTCAGGTCTTGAGTCGCTTTCAACCTGAAGAGTCAC
ATGCGAACTCACACTGGAGAGAAGCCCTATGTGTGCCAGGAATGTGGCGAGCCTTCAGTGAACCTCA
TCCCTCAGGAAACATGCAAGGACTCACAGTGGCAAGAAGCCCTATGCATGCCAGGAATGCGGGCAGCC
TTTGGTCACTTTCACATCTTATTGTACATGTGAGAACACACAGTCCCGGAGACCCATCAATGTAAAT
CAGTGTGAGAAAGCCTTACGGACAGCTCCTCACTCACTGTACACAAAAGAACCCATGTGGGAAGAGAG
ACCATTAGGAATGGCAGCCTGCCTTTATCCATGTCTCATCCATACTGTGGGCCCTTGTCTAATTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:** SgfI-MluI

**ACCN:** NM\_032433

**Insert Size:** 1998 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_032433.2</a>
<b>RefSeq Size:</b>	3767 bp
<b>RefSeq ORF:</b>	1998 bp
<b>Locus ID:</b>	84449
<b>UniProt ID:</b>	<a href="#">Q96JL9</a>
<b>Cytogenetics:</b>	19p13.12
<b>Domains:</b>	KRAB, zf-C2H2
<b>Protein Families:</b>	Transcription Factors
<b>MW:</b>	75.5 kDa
<b>Gene Summary:</b>	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) encodes the longest isoform (1).