

Product datasheet for **SC314363**

GIMAP8 (NM_175571) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	GIMAP8 (NM_175571) Human Untagged Clone
Tag:	Tag Free
Symbol:	GIMAP8
Synonyms:	IAN-9; IAN6; IAN9; IANT
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC314363 representing NM_175571.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTCAGAGCAGAGCTGCCAGATGTCCGAATCGCGCTCCTCCTCTGGGAAAAATGCCGCTCGGAAAA
AGTCCACAGGAAATGCCATTCTGGGCAACATGTGTTCAAGTCCAAGTTCAGTGATCAGACAGTGATC
AAAATGTCCAGAGAGAGATTGGGTCTGAGAGAAAGGAAGTTGTGGTAATTGACACCCTGACCTT
TTCTCCTCAATAGCTTGTGCTGAAGACAAGCAACGCAACATCCAACACTGCTTGGAGCTCTCTGCTCCC
AGCCTCCATGCTCTGCTTTGGTAATTGCCATCGGCCATTTACAAAGGGAGGATGAGGAAACAGCCAAG
GGCATCCAACAAGTGTGGAGCTGAAGCCAGGAGGCACATCATTATTGTCTTCACTCGGAAGGATGAT
TTGGGGGATGACTTGTGCAAGATTTTATTGAAAAAACAACCTCTCAAGCAGTTGGTTCAAGACTAT
GAGGGCCGATACTGCATTTTCAACAACAAGACCAATAGTAAGGATGAGCAGATCACCCAGGTGTTGGAG
CTCCTTCGCAAGTTGAGTCTTTGGTGAATACGAACGGAGGCCCTATCATGTGAACTTCAAACCTGAA
GGCAGCAGTTTTCAAGATTGTGTGAATGAAGCTGCATCTCAAGAGGGAGACAAGCCACAGGGCCCAAGG
GAAAGGCAGCTGCAGTCCACAGGACCCGAGCAGAATCCGGGGACATCAGAACTGACAGTCCCTCTGTG
GGGAAACGCGGTGCTGGAAAAAGTGCAGCAGGAAACAGCATTCTGGGGAGGCAGGCCCTTTCAGACCGGA
TTTAGTGAGCAGTCAGTAACCCAGAGCTTCTTGTCTGAGAGCAGAAGCTGGAGAAAAAGAAAGTTTCG
ATCATTGATGCTCCGGACATCTCATCTTTAAAGAACATTGACTCAGAAGTTAGAAAAACACATCTGTACA
GGCCCCATGCCTTCTGCTGGTACACCACTGGGCTTTTACACTAAGAATGATGAGGCAGTGTGAGC
ACCATCCAAAACAATTTGGAGAAAAATCTTTGAGTACATGATCATACTTCTTACCAGGAAAGAAGAT
TTAGGGGATCAGGATCTAGATACGTTCTTAAGAAACAGCAATAAAGCTCTCTATGGTCTCATCCAGAAG
TGTA AAAACAGATATAGTGCCTTCAACTACCGGCAACAGGAGAAGAAGAGCAAAGGCAGGCGGACGAG
CTCCTGGAAAAAATTGAGAGCATGGTGCATCAGAATGGGAACAAGCATTGTGTTTTTCAGAGAAAAAGAA
ACCCTGAACATTGCTTGTGGGAGAAGCGGGACTGGGAAGAGTGCGACCGGGAACCTATCTCTGGGG
AGCCTCGTCTTACCTCTCGGCTCCGGGCCAGCCAGTCACCAAGACCAGCCAGAGTGGCAGGAGGACA
TGGGACGGACAGGAGGTGGTGGTGTGGACACTCCTTCTTCAACCAGATGCTGGATGTCGAAAAGGAC
CCATCCCGTTAGAAGAGGAGGTCAAGCGCTGTTGTCTGCTGTGAAAAAGGGGACACATTTTTTGTG
CTGGTGTCCAGCTGGGACGATTCACTGAAGAGGACAAAACAGCTGTGGCGAACTGGAGGCCATCTTT
GGAGCAGACTTTACGAAATACGCGATTATGCTGTTCAACCAGGAAAGACCTAGGGCGGGGAATTTG
GAAGACTTCATGAAGAACTCAGATAACAAGCCCTTCGGCGCATTTTTAAAAAGTGTGGCGGCGAGTT
TGTGCTTTTAAACAACAAAGAAACAGGCCAGGCCAGGAAACCCAGGTGAAAGCTCTTTTAAACAAGGTC
AATGATCTGAGAAAAGAAAGTGGGTGGTCCGGGTATCCCATACACAGGAGAACGTCAGCAAACCTAAT
AAAAATGTCCAGGAAATGTCCCAAGCCGAAAACTCCTTAAAAATTTAATAGGTATTTTACAATAG
ACCGGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_175571

Insert Size: 1998 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_175571.3](#)

RefSeq Size: 4196 bp

RefSeq ORF: 1998 bp

Locus ID: 155038

UniProt ID: [Q8ND71](#)

Cytogenetics: 7q36.1

MW: 74.9 kDa

Gene Summary: This gene encodes a protein belonging to the GTP-binding superfamily and to the immuno-associated nucleotide (IAN) subfamily of nucleotide-binding proteins. In humans, the IAN subfamily genes are located in a cluster at 7q36.1. [provided by RefSeq, Jul 2008]