

Product datasheet for **SC314405**

ZNF665 (NM_024733) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Tag Free
Symbol:	ZNF665
Synonyms:	ZFP160L
Vector:	<u>pCMV6 series</u>



Fully Sequenced ORF: >NCBI ORF sequence for NM_024733, the custom clone sequence may differ by one or more nucleotides

```

ATGGCTCTTCTCAGGGACAGTTGACATTCAAGGATGTGGCCATAGAATTCTCTCAGGAG
GAGTGGACATGCCTGGACCCTGCTCAGAAAGACTTTGTACAGGGACGTCATGTTGGAGAAT
TATAGGAACCTGGTCTCCTGGATATCTCTTGTAAATGTGTAACACGGATTTGCCACCA
AAGGGGAAGAACATATGGGAGAAGCGTTCTACACGGTGAAGTTGGAGAGACTTGAAAGC
TGTGACACTGCAGGCTTGTCTTCCAGGAAGTTCAGAAAAATACATACGACTTTGAGTGT
CAGTGGAAAGATGATGAAGGAAATTATAAACAGTACTTATGTTGCAAAAAGAAAAATCTC
CCTGGTAGAAGAGCTCAACGTGATAGAAGGGCTGCAGGAAACAGGCATATTGAAAAATCAG
CTTGGAGTAAGCTTTTCACTCACATCTCCCTGAAGTGCAGCAATTTCAACGTGAAGGGAAA
ATTTATGAATACAATCAAGTTGAGAAGTCTCTAATAATCGAGGAAAACATTATAAATGT
GATGAATGTGGCAAGGCTTTCAGTCAAACTCACGGCTAACAAAGTCATAAGAGAATTCAT
ACTGGAGAGAAGCCTTACCAGTGAATAAGTGTGGCAAGCCTTTACTGTTCTGTTCAAAC
CTAACATCCATCAGGTCATCCATACTGGAGAAAAACCTTACAAATGTAATGAATGTGGA
AAGGTCTTCAGTCAACCTTCAAACCTTGCAGGTCATCAGAGAATTCATACTGGAGAGAAA
CCTTACAAGTGTAATGAGTGTGGCAAGCCTTTAGAGCACATTCAAACTAACTACACAT
CAGGTCATCCATACTGGAGAAAAACCTTACAAATGTAAGGAATGTGGCAAGTGCTTCACT
CAAAATTCACACCTTGAAGTCATCGAAGAATTCATACTGGGAGAAAGCCTTACAAGTGT
AATGAGTGTGGCAAGCCTTTAGTGTTCGCTCAAGCCTGACTACCCATCAGACAATCCAC
ACTGGAGAAAAACCTTACAAATGTAATGAATGTGGCAAGGCTTTCAGGCACAATTCATAC
CTTGCAAAGCATCGGCAATTCATACTGGTGAAGAACCTTACAAGTGAATGAGTGTGGG
AAAGCCTTCAGTATGCATTCAAACCTAACTAAGCATCAGATCATCCATACCGGAGAAAAG
CCTTTCAAATGTAATGAATGCGTCAAGGTTTCACTCAGTATTCACACTAGCAAATCAT
CGAAGAATTCATACTGGAGAGAAACCTTACAGGTGTGATGAGTGTGGCAAGCCTTTAGT
GTGCGATCAAGCCTAACTACCCATCAGGCAATTCATACTGGAGAAAAGCCTTACAAATGT
AATGACTGCGGTAAAGTCTTCACACAAAATTCACACCTTGAAGTCATCGGGGAATTCAT
TCTGGAGAGAAACCTTACAAGTGTGATGAATGTGGTAAAGCCTTCAGTCAAACATCACAA
CTTGCAAGGCATTGGAGAGTTTCACTGGAGAAAAACCTTACAAGTGAATGAGTGTGGC
AAAGCCTTTAGTGTTCATTCAAGCCTAACTATACATCAGACAATACATACTGGACAAAAA
CCATACAAATGTAATGATTGCGGCAAGGCTTTCAGACACAATTCGTACCTTGCAATTCAT
CAGAGAATTCACACTGGTGAAGAACCTTATAAGTGAATGAGTGTGGCAAGCATTAGT
GTACATTCAAACCTAGCTACCCATCAGGTCATCCATACTGGAGAAAAACCTTACAAATGT
AACGAATGTGGCAAGGCTTCACTCAAACTCACATCTTGCAATCATAGAAGAATTCAC
ACTGGAGAAAAACCTTATAGGTGTAATGAGTGTGGAAAGCCTTCAGTGTTCGTTCAACC
CTAACTACCCATATGGCAATCCATACTGGAGACAAACCTTACAAATGTAAACCAATGTGGC
AAGGTCTTTACTCAAAATTCAAACCTTGCAAAACATCGAAGAATTCATAGTGGA

```

Restriction Sites: Please inquire

ACCN: NM_024733

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_024733.2, NP_079009.2</u>
RefSeq Size:	3493 bp
RefSeq ORF:	2037 bp
Locus ID:	79788
UniProt ID:	<u>Q9H7R5</u>
Cytogenetics:	19q13.42
Domains:	zf-C2H2
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]