

Product datasheet for **SC314410**

MS4A14 (NM_032597) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MS4A14 (NM_032597) Human Untagged Clone
Tag:	Tag Free
Symbol:	MS4A14
Synonyms:	MS4A16; NYD-SP21
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC314410 representing NM_032597.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGAGTCAACATCCCAGGACAGAAGGGCAACTCAGTCATCACTATAAAACCAAACGAACTGTATTG
ACTGCATTTCCCTACAGACCTCATAGCTCTCTGCTGGATTTTCTGAAGGGAGAGCCAAGAGTCTTGGGG
GCTACCCAGATCCTGCTTGCTCTAATCATTGTGGGCTTTGGAAGTATATTTGCACTTAATTACATCGGT
TTCTCCCAAAGACTTCCCTTGTGTCTCACAGGATATCCATTCTGGGAGCACTTATTTTATTCTT
ACAGGATACCTCACAGTAACCGATAAGAAATCAAACTTCTGGGCAAGGTGTACGGGCATGAATGTT
ATCAGCTCCTTGTTGCGATAACTGGGATTACTTTCACCATTCAGCTACAGACATCAAGACAAGTAC
TGCCAGATGCCATCCTTTGAAGAAATATGTGTTTTAGTAGAACTCTTTTATTGTTCTGTTTTCTTG
CCTTCGGATGTTACTCAAATAGTGAACAACCTGCCCCAGAAGAAATGATCAATTACAATTTGTGCTT
CAAGAAGAGTTTTCCAGTGATGATTCAACAACAAATGCACAATCTGTTATCTTTGGAGGCTATGCTTTC
TTCAAGTTAACTCTCTAGGAGTCTTTAGTCTCCCAACAGGTAATAAAGGTAGAGAATTTGTGCCA
GATGAACAAAAGCAAAGTATCCTTCCATCTCCCAAATTTTCAGAGGAAGAAATTGAACCTTTGCCTCCC
ACACTAGAGAAAAAGCCCTCAGAAAATATGTCCATTACAGCTAGACTCTACATTTAAACAAATGAAAGAT
GAAGATCTACAATCTGCTATTGTACAACCTTCTCAATGCAAACCAAGCTTCTGCAGGACCAAGCTGCG
TCACTCCAAGTTTTTCCATCCCATTCTGCACTAAACTCGAAGATATATCACCTGAAGACTTGCCATCC
CAAGCTCTACAGTAGAAGGCCTGTGAGAACAAACCATGCCATCTAAGTCTACATCATCCCATGTCAA
CAGTCTTCTAATCTGACAGCTAATGACCTGCCCCCTCAAGGCATACTATCCCAAGACACATCATCTCA
GATATGCTGTTTCATGACATGACATCCCAAGATATGCAATCCCTAGATATGCTATCTCAAGACACACCA
TCCCAGCCATGCCACCTCAAGACATACCTTCCCAAGATATGCTATCCCAAGCTCTATCAGCGCATGCC
ATATTACCTGAAGCCTCAACATCCCATATTGTGCAGTTCCCTGAAATACAACACCTACTTCAGCAGCCC
CCAGATCTTCAACCAGAAAACTGAACCTCAAAACCAGCAAATTTTACAATGTCATATCAAGATATT
AGATCAGAAGTTATGGAAGAGACCAAGAATGGAATCTGAGGAGGAAGTCCATAGAAGAAAACTCTCA
AGACGGCATTCTTAAACCAGCAAACCAAGCCTTGCAATACTTAAGGAGACATTCTTTAGACGTGCAA
GCCAAAGGCCAGAAATCCTCAAAGAGGCATTCTTAGATCAGCAAAGCAAAGGCTGGCAATCTCCAAAG
CAGAAATCCTTAGACCAGCAAATCAAAGACTGGCTATCCCCAAAGAGGCACTCCGTAGATAAGCAAGCT
CAACTTAATCAAATAAGAGCAACTCCAGATCAGCAAGCTGAAGATCAGCAAGCCAAAGGGGAACAA
TACCCAGAAGGACAATCTAAAGATGGACAAGTTAAAGACCAGCAGACTGATAAGGAGCAAACTCAAAG
AAGCAAACCCAGGATCAGCAAACCTGAAGACCAGCGGCCCAAGAGAAGAAATCCCGAAAGGACAATTC
CAAAATGTTCAAGCCGAAGGACAGCAAGCTCAGGTGGAGAAAGTGCCAAACTGTTATGCCAAGATTCA
GAATCCCAATACAGCAATACCAATTCTGGCAATTCACAAAGGCAATCTCCAGGCTGGACAACCCAGG
ACTGTCAATCTTTGGCCAAGAATCCCCTGACTGGATAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCG
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_032597

Insert Size: 2040 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_032597.4</u>
RefSeq Size:	3353 bp
RefSeq ORF:	2040 bp
Locus ID:	84689
UniProt ID:	<u>Q96JA4</u>
Cytogenetics:	11q12.2
Domains:	CD20
Protein Families:	Transmembrane
MW:	76.6 kDa
Gene Summary:	May be involved in signal transduction as a component of a multimeric receptor complex. [UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) has multiple coding region differences, compared to variant 4, which results in a shorter protein (isoform 1), compared to isoform 4. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.