

Product datasheet for **SC314417**

MPHOSPH10 (NM_005791) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MPHOSPH10 (NM_005791) Human Untagged Clone
Tag:	Tag Free
Symbol:	MPHOSPH10
Synonyms:	CT90; MPP10; MPP10P; PPP1R106
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC314417 representing NM_005791.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGGCGCCGACAGTCTGGCGTCGACGGACCCTGGAGCGGTGTCTGACGGAAGTCGGCAAAGCCACGGGT
CGGCCCGAGTGCTTCTCAGATTCAAGAGGGATTGGCATCAAAGTTCATTCTTTAACAAAAGTGCTT
TATGACTTTAATAAAATATTAGAGAATGGTAGGATCCATGGAAGCCCTTGCAAAAACCTTGATAGAA
AATTTTGTATGATGAGCAGATTTGCAACAACCTGGAATTGCAAAATGAACCAATTTTACAATACTTTCAG
AATGCAGTTAGTGAAACAATTAATGATGAAGATATCAGTCTTCTCCAGAGAGTGAAGAACAGGAACGT
GAAGAGGATGGTTCAGAGATAGAGGCTGATGACAAGGAGGACCTAGAAGATTTAGAGGAGGAGGAAGTG
TCCGACATGGGTAAATGATGATCCTGAAATGGGTGAGAGAGCTGAAAACCTCAAGCAATCTGATCTGAGG
AAAAGCCCGTTTTCAGTGATGAGGATTCTGACCTTGACTTTGATATCAGCAATTTGGAACAGCAGAGC
AAGGTGCAAAACAAAGGACAGGAAAACCAAGAGAAAAGTCCATAGTAGATGATAAATCTTCAAATCTC
TCTGAAATGGAGCCTATTTAGAAAACATAGAAAAAGAAGAGGAACGAAAAGATGATAATGATGAGGAG
GAGGAAGATATTGATTTTTTTGAAGATATTGATTCTGATGAAGATGAAGGGGACTGTTTGAAGTAAA
AACTTAAGTCAGGTAAAAGTTCCAGAAATCTGAAATACAAAGATTTTTTTGATCCAGTTGAAAAGTGAT
GAAGACATAACAAATGTTTCATGATGATGAGCTGGATTCAAACAAAGAAGATGATGAAATTGCTGAAGAA
GAAGCAGAAGAACTAAGTATTTTCGGAACCGGATGAAGATGATGACCTTCAAGAAAATGAAGACAATAAA
CAACATAAAGAAAGCTTGAAAAGAGTGACCTTTGCTTTACCAGATGATGCGGAACTGAAGATACAGGT
GTTTTAAATGTAAAGAAAAATTTCTGATGAAGTTAAATCCTCCTTTGAAAAAGACAGGAAAAGATGAAT
GAAAAAATTCATCTTTAGAAAAAGAGTTGTTAGAAAAAAGCCGTGGCAGCTTCAGGGGGAAGTGACA
GCACAGAAGAGGCCAGAGAACAGCCTCCTGGAGGAGACCCTACACTTTGACCATGCTGTCCGGATGGCA
CCTGTGATTACAGAGGAACACCCTTCAACTGGAAGATATCATTAAACAGAGGATAAGAGATCAGGCT
TGGGATGATGTAGTACGTAAAGAAAAACCTAAAGAGGATGCATATGAATATAAAAGCGTTTAACCTTA
GACCATGAGAAGAGTAAATTGAGCCTTGCTGAAATTTATGAACAGGAGTACATCAAACCTCAACAGCAA
AAACAGCAGAAGAAGAAATCCAGAACATGTAGAAATTCAGAAGATGATGGATTCCCTCTTCTTAAAA
TTGGATGCCCTCTCAAATCTCACTTTATCCCTAAACCGCCTGTACCAGAGATTAAGTTGTGTCAAAT
CTGCCAGCCATAACCATGGAGGAAGTAGCCCCAGTGAGTGTTAGTGATGCAGCTCTCCTGCCCCAGAG
GAGATCAAGGAGAAAAATAAAGCTGGAGATATAAAACAGCTGCTGAAAAACAGCTACAGACAAGAAA
CGAGAGCGAAGGAAAAAGAAATATCAAAAGCGTATGAAAATAAAGAGAAGGAGAAGCGGAGAAAACCTG
CTTGAAAAGAGCAGTGTAGATCAAGCAGGAAATACAGCAAAACAGTAGCTTCGGAGAAGTTAAACAG
CTGACCAAACTGGCAAAGCTTCCTTCATAAAGGATGAAGGTAAAGACAAGGCCTTAAAGTCCTCTCAA
GCATTCTTTTCTAAATTACAAGATCAAGTAAAAATGCAATCAATGATGCAAGAAAAACAGAAAAGAAA
AAGAAGAAAAGACAGGATATTTCTGTTCAATTAAGCTGTA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_005791

Insert Size: 2046 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_005791.2</u>
RefSeq Size:	2488 bp
RefSeq ORF:	2046 bp
Locus ID:	10199
UniProt ID:	<u>O00566</u>
Cytogenetics:	2p13.3
MW:	78.9 kDa
Gene Summary:	This gene encodes a protein that is phosphorylated during mitosis. The protein localizes to the nucleolus during interphase and to the chromosomes during M phase. The protein associates with the U3 small nucleolar ribonucleoprotein 60-80S complexes and may be involved in pre-rRNA processing. [provided by RefSeq, Dec 2010]