

Product datasheet for **SC314423**

TXLNB (NM_153235) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TXLNB (NM_153235) Human Untagged Clone
Tag:	Tag Free
Symbol:	TXLNB
Synonyms:	C6orf198; dj522B19.2; LST001; MDP77
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC314423 representing NM_153235.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGGAGGCTAATCACTCTGAACAGCTCTCAGCGGAACGACAGTCAACACCTCCAGGTGACAGTTCATCA
TTACCCAGTCACAATGGCCTGGAGAAGGAAGATGGCCAGGATTCTCCAACCCAGTCCAACCCAGAGAG
AAAGAGGCAAGTGTGCACCCCGATATCTCTGAAGAGCTGAATCGACAGCTGGAAGACATCATTAACACT
TATGGGCTGCTGCCAGCACAGCAGGGAAGAGGGCTCTGCCAGGGCCAGTGAGCAGCCTGAGAATGCA
GAATCACCTGACAACGAGGATGGGGACTGTGAGGAAACAACGAAGAGGCTGGAAGAGAACCCTGTTGCT
TCTGGAGAGCCACCCACTGTCAAAGAGCCCGTCAGCAATAAGGAGCAAAAATTGGAAGAAAAATCTTA
AAAGGATTAGGCAAGAAGCCAACTGCTAATGCAAAATCTGAACAAGTTGCAACACCGGAAGAAAAAG
TTTGATTTTTTATTCAAGAAGTATGCTGAATTGCTGGATGAACATCGTACTGAGCAAAAAGATTAAAG
CTCCTCCAAAAGAAACAGGTACAAATTCAAAAAGAAAAGGACCAGTTACAAGGTGAACACAGCAGAGCT
ATCCTCGCTCGAAGCAAAATGGAGAGTCTGTGCCGGGAGCTGCAGAGACACAACAAGACTCTGAAGGAA
GAGGCGCTTCAGCGGGCAGCTGAGGAAGAAGAGAAAAGGAAGGAATCACAAGCCATTTCCAGAGTACC
CTCACGGACATCCAGGGCCAGATCGAGCAGCAGAGTGAGCGAAATATGAAGCTCTGTGAGGAGAACACA
GAGCTTGCAAAAAGCTGAAAAGCATCATCGATCAGTATGAGCTCAGAGAGGAGCATCTGGACAAAAATA
TTTAAACACAGAGAACTGCAGCAGAAGCTGGTGGATGCAAAGCTTGAGCAGGCCCAAGAAATGATGAAG
GAAGCGGAGGAGCGACACAAACGAGAAAAGGAATATTTGCTGAACCAGGCAGCAGAGTGGAAGTTTCA
GCGAAAGTGCTGAAGGAGCAAGAGACAGTCTGCAGGCTCAGCTCACTCTCTACTCAGGAAGTTTGAA
GAATTCAGAGCACACTAACTAAAAGCAACGAGGTGTTTGCCACGTTCAAACAGGAAATGGACAAAAACA
ACTAAGAAAAATGAAGAAGCTGGAAGAGGACACAGCCACATGGAAGGCCGATTTGAGAAGTGTAAACAA
GCTCTGTTGGACATGATTGAAGAGAAAGCACTGAGAGCTAAAGAATATGAGTGCTTTGTGATGAAAATC
GGGAGGCTAGAGAACCTTGCCGTGCTTTACAAGAAGAGAGAAACGAATCCACAAAAAATCAGAGAC
GCAGAAATATCTGAAAAGGATGACCAAGTCAGCACAACCTCCGATGAAGAGCCAGAGTCAACAGTCTCT
GTGGATCAAGAGATTGACGCAGAGGAGTTAATAGTGTCAAACCGCCGTGAAAAATCTGCCACAGCC
TTCATGATAATTCATCATCCAGAGTCAACCCCGCACCAGTCCAAAGAAACCAACCCGAAATAGGCAGT
TCTCAGGAGAGTGCTGACGCCGCTCTCAAGGAGCCAGAGCAACCCCTCTGATCCCTTACGCGGATTCA
GAGAGTCCCTGCCTCCCTAACTCCTCAGGCTGAAGCCGAAGGAGGCAGTGATGCTGAACCTCCCTCC
AAGGCCAGTAATTCTCTGCCGGTTGGGAGCAGAAACCAATGCGAGGGTCTCCCTGTTGGAGCACAG
GCTGATCAGGCGCTCTGGAAGCCAGAGGCAGAAGCTTCCGGTCAGGCCCCACAGGCTCCACCCAGAGCC
TCCCTACAGAAGATGGAGGCAGATGTGCCTGCTCCAGCATGCGCAGCAGAGAAGAGCACGTTGCAGCCATG
GTGCCTGCATGCGAGCCAGTAGGCAGCCCCACGAGCAGCAGCAGAGGAGCTGCCAGTAGGGGCCCTCA
GCTGGGCCCCAGCCGCGCAACGTGGCTGACACCAATCTGGAAGGCGTCGACTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_153235

Insert Size: 2055 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_153235.3</u>
RefSeq Size:	4685 bp
RefSeq ORF:	2055 bp
Locus ID:	167838
UniProt ID:	<u>Q8N3L3</u>
Cytogenetics:	6q24.1
MW:	76.5 kDa
Gene Summary:	Promotes motor nerve regeneration (By similarity). May be involved in intracellular vesicle traffic.[UniProtKB/Swiss-Prot Function]